Introduction to Surface Area Analysis



Jeffrey Bodycomb, Ph.D. Miro Pejcinovic, Amy Hou HORIBA Scientific www.horiba.com/us/particle





Why determine surface area?
Surface area and size
Determining surface area

Explore the future

HORIBA





Surface area directly correlates with desired properties.

Reactivity
Dissolution
Catalysis

Separation



Explore the future

HORIBA

© 2014 HORIBA, Ltd. All rights reserved

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Applications

- Catalysts
- Adsorbents
- Pigments
- Pharmaceutical products
- Chromatographic carriers
- Sintered materials
- Building materials
- Ceramics

Filters







When a solid is involved in a chemical reaction, either as a reagent or a catalyst, the surface area is the only accessible area for the reaction.

Dramatic example: dust explosions!



Dissolution



Dissolution of materials (e.g., API's) depends on surface area.

Noyes-Whitney Equation

$$\frac{dm}{dt} = A \frac{D}{d} \left(C_s - C_b \right)$$

m = mass dissolved material

t = time

A = <u>Surface area of interface</u>

D = Diffusion coefficient

d = Boundary layer thickness

 C_s = Concentration of substance on surface

 C_{b} = Concentration of substance in solvent





Why determine surface area?Surface area and size

Determining surface area





When a particle of a given volume is broken into two parts as shown in the figure, total volume does not change. The total surface area, however, does change. It INCREASES by the amount of the two newly-exposed edges.



This simple illustration demonstrates relationship between surface area and particle size. As ratio of surface area to volume increases, surface phenomena come into play. For this reason, measurement of surface area becomes much more important for small particles.









Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Surface Area and Size





Surface area mean size can be arithmetically calculated from a measured volume based particle size distribution.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Explore the future

© 2014 HORIBA, Ltd. All rights reserved

HORIBA

Surface roughness

- Surface area measurement using nitrogen adsorption detects total surface area of the particles.
- Conventional particle sizing detects only the overall structure of the particles.
- All of these particles will produce same result using conventional particle size measurement techniques.
- Nitrogen adsorption methods will yield progressively higher results for particles with greater irregularities.







- We can use a size analyzer, but eventually prefer to measure surface area directly.
- Specific surface area for a given size changes according to shape
- Another way to think about determining particle shape or surface roughness





Why determine surface area?
Surface area and size
Determining surface area





To find surface area, "simply" determine amount of N_2 adsorbed to the surface in a single layer. Use the projected surface area of each molecule times the number of molecules to find surface area.



Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Need a single layer, but how do we find a single layer? Need to look at adsorbed N_2 as a function of pressure.



Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Surface area





Explore the future

© 2014 HORIBA, Ltd. All rights reserved.

Multi-point method

- Used for unknown materials to determine behavior of the material.
- Measure multiple points to two parameters (a material constant and the surface area).





Single-point method



 Once that determination has been accomplished, a SINGLE POINT MEASUREMENT is quite satisfactory for many (90– 95%) practical cases.



Relative Pressure

HORIBA

Explore the future

Monolayer volume





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Monolayer volume







Measuring for BET equation

$$\frac{P/P_0}{W(1-P/P_0)} = \frac{C-1}{W_m C} (P/P_0) + \frac{1}{W_m C}$$

 $P = N_2 \text{ partial pressure } = \text{controlled by instrument (MFC)}$ $P_0 = N_2 \text{ equilibrium vapor pressure } = \text{known}$ $W = \text{mass adsorbed } N_2 = \text{measured from desorbed volume}$ of nitrogen $W_m = \text{mass of } N_2 \text{ monolayer } = \text{found from equation}$ C = BET constant = found from equation

© 2014 HORIBA, Ltd. All rights reserved



- The single-point method B.E.T. method is applicable to many, but not all samples.
- For unknown materials, the multi-point method is best. However, for most materials, the single point method applies. Single point is preferred because it is much simpler and faster to use.

Measure at a single N₂ pressure.

References:

- Brunauer, S., Emmett, P., and Teller, E., J. Amer. Chem. Soc. 60, 309 (1938)
- Gregg, S.J., and Sing, K.S.W., Adsorption, Surface Area & Porosity, Academic Press, (1967)



BET equation

$$\frac{P/P_0}{W(1-P/P_0)} = \frac{C-1}{W_m C} \left(\frac{P/P_0}{P_0}\right) + \frac{1}{W_m C}$$
$$C \to \infty \qquad \frac{1}{W_m C} \to 0 \quad (C-1) \to C$$

$$W\left(1 - \frac{P}{P_0}\right) = W_m$$

Measure W at one value of P/P_0 and obtain W_m .

Explore the future





- We obtain W_m from linear plot (multipoint) or single point.
- Use W_m (quantity of N_2) and cross sectional area of N_2 to determine surface area:

$$S_t = \frac{W_m N_A A_x}{M_w}$$

 $\begin{array}{l} S_t = \mbox{Total surface area} \\ W_m &= \mbox{Mass of N}_2 \mbox{ monolayer} \\ M_w &= \mbox{Molecular weight of N}_2 &= 28 \mbox{ g/mol} \\ N_A &= \mbox{Avagadro's Number} &= 6.02 \mbox{ x } 10^{23} \mbox{ molecules/mol} \\ A_x &= 16.2 \mbox{ x } 10^{-20} \mbox{ m}^2/\mbox{molecule} \end{array}$

Single point measurement errors



If we know c, we can predict error in single point results.



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Getting to single point when life is bad HORIBA

Single-point method offers the advantage of simplicity and speed, with acceptable accuracy.

$$V_m = V_a \left(1 - P/P^o \right)$$
 i.e. $V_m = 1/slope$

A relative pressure of 0.3 gives good general agreement with the multipoint method.

Correction of single point "error" at $P/P_0 = 0.3$ by multiplying the single point BET value by C/C-2 decreases the difference.

Sample No.	Multi-point BET (m²/g)	Uncorrected single-point (m²/g)	Uncorrected difference (%)	Corrected single – point (m²/g)	Corrected difference (%)
1	4.923	4.241	-13.9	4.948	0.51
2	4.286	3.664	-14.5	4.275	-0.26
3	8.056	6.867	-14.8	8.011	-0.56
4	5.957	5.194	-12.8	6.060	+1.73

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Not all of sample brought to lab is analyzed

Bulk or process stream	Gross sample	Laboratory sample	Test sample	Measurement sample
10 ⁿ kg	≻kg	<kg< td=""><td>g</td><td>mg</td></kg<>	g	mg

Must sub-divide sample

How to introduce representative sample into instrument?



Sampling from Drums



Powder Thief





www.samplingsystems.com

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



When a powder is stored in a container, it can be mixed by rolling and tumbling the container. The container should <u>not be more than</u> <u>half to two-thirds full</u>. It is important to perform this action before "grabbing" a sample with a spatula.



Technique: Coning & Quartering



Pile of powder is divided into 4 sections.

Two diagonal sections are discarded, and two are retained and mixed together.

Mixture is again divided into 4 sections, and two diagonal sections are mixed.

Process is repeated until remaining sample is correct amount for analysis.

Can be carried out with very small sample amount or very large samples.





Technique: Chute Riffling



Chute splitting allows sample to vibrate down a chute to vanes which separate the mass into two portions. Each portion moves further where they each are divided into two parts, now giving four parts. This may be continued until usually 8 or 16 portions are obtained.





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Sampling Technique Error Levels



Standard Deviation (σ) in % Sugar-Sand Mixture

SCOOP SAMPLING	6.31
TABLE SAMPLING	2 .11
CHUTE RIFFLER	1.10
SPINNING RIFFLER	0.27



Density of sand and sugar respectively 2.65 and 1.64 g/ml

Reference: Allen, T. and Khan, A.A. (1934), Chem Eng, 238, CE 108-112

Method	Relative Standard Deviation (%)
Cone & Quartering	6.81
Scoop Sampling	5.14
Table Sampling	2.09
Chute Riffling	1.01
Spin Riffling	0.125





Frequent source of error.

Use a balance with four places.

Carefully determine and record mass of empty sample cell. You will need this value after measurement.



Add sample



Record mass of cell+sample



Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Explore the future

© 2014 HORIBA, Ltd. All rights reserved

Degassing

Material in Sample Cell must be OUTGASSED

Sample cell containing sample must be heated to over 100 C° for at least 20 to 30 minutes to drive off impurities that adhere to surfaces.











- Need to heat sample under inert gas to drive off moisture and volatile components.
- Very sample dependent.
- High temperatures can affect some materials.
- Low temperatures require a longer time.
- Good degassing will make the results independent of sample history (e.g., humidity during manufacturing).
- Typically 1-2 hours at 200-300 C
- Inorganic oxides 120 to 150 C (NIST practice guide)
- Microporous carbons and zeolites: 300 C

Measurement



- Sample cell is partially filled with sample material, immersed into liquid nitrogen to cool sample so nitrogen gas will condense on particle surfaces.
 - Adsorption Phase: Nitrogen gas condenses on particle surfaces, resulting in a LOSS of gas.
 - Desorption Phase: Nitrogen gas evaporates from particle surfaces, resulting in a GAIN of gas.



W, mass adsorbed



To measure adsorbed mass we can measure loss of N_2 flowing across surface (adsorption process) or "extra" N_2 that appears as N_2 flows across surface (desorption process).



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Comparison









Make measurement at different N2 concentrations. Note baseline steps.



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



- Weigh after degassing (very important)
- Did you keep your tare value?
- Mass determination is a major source of uncertainty. For 0.1 gram of sample an weighing error of 0.001 gram (1 milligram) will give a 1% error in results.

Standards



NIST 2207

Certified mean value of single point surface area is 174.2

	2207	
Test 1	174.07	
Test 2	179.15	
Test 3	171.93	
Average	172.2	
S.D.	1.77	
COV	1.0%	

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Real sample



- Needle-like material
- By microscopy:
 - •Length: 50 microns
 - Diameter: 2 microns
 - Density: 1.4 g/cc



- •Estimate specific surface area of 1.5 m^2/g
- By flowing gas BET: 3.5 m²/g
 - Diameter is close to microscope resolution limit (0.5 microns)
 - Surface roughness

Explore the future

Repeatability



Two different ceramic powders, run repeats.

	P1	P2	
Test 1	90.53	3.69	
Test 2	88.21	3.59	
Test 3	88.18	3.58	
Average	88.97	3.62	
S.D.	1.34	0.0608	
COV	1.5%	1.6%	



Another sample



Sample	Mass, g	Mass loss on degassing, %	Specific Surface Area, m²/g
Sample A, Split 1	0.4372	8.0	3.52
Sample A, Split 2	0.5018	8.9	3.41

This sample had a bigger than usual loss of mass on degassing. So degassing is particularly important.







Flowing gas BET advantages

- Fast (5 to 10 minutes)
- Repeatable
- Accurate
- Direct measurement of surface area



The SA-9600 surface area analyzer





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

SA-9600 single station





•2 Outgassing Stations

•1 Analysis Station

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



SA-9603 multi-station





© 2014 HORIBA, Ltd. All rights reserved.



