### **Laser Diffraction Theory**





Jeffrey Bodycomb, Ph.D. HORIBA Scientific www.horiba.com/us/particle

© 2012 HORIBA, Ltd. All rights reserved

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



#### When a Light beam Strikes a Particle



Small particles require knowledge of optical properties:

- Real Refractive Index (bending of light)
- Imaginary Refractive Index (absorption of light within particle)
- Refractive index values less significant for large particles
- Light must be collected over large range of angles



#### **Diffraction Pattern**



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

#### HORIBA

#### **Diffraction Patterns**



Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

#### HORIBA

# Using Models to Interpret Scattering

Scattering data typically cannot be inverted to find particle shape.

We use optical models to interpret data and understand our experiments.





## **The Calculations**

- There is no need to know all of the details.
- The LA-950 software handles all of the calculations with minimal intervention.

### This talk is to improve your understanding of the results.





# **Laser Diffraction Models**

- Large particles -> Fraunhofer
  - More straightforward math
  - •Large, opaque particles
  - Use this to develop intuition
- All particle sizes -> Mie
  - Messy calculations
  - All particle sizes





Expressed in just in y-direction

$$E = E_0 \sin(ky - \omega t)$$

$$H = H_0 \sin(ky - \omega t)$$

Oscillating electric field Oscillating magnetic field (orthogonal to electric field)

Cheat by using vector notation (for 3 dimensions) and real part of complex numbers (because it makes the trig easier).

$$\mathbf{E} = \mathbf{E}_0 \exp(i\mathbf{k} \bullet \mathbf{x} - i\omega t)$$
$$\mathbf{H} = \mathbf{H}_0 \exp(i\mathbf{k} \bullet \mathbf{x} - i\omega t)$$



Complements of Lookang @ weelookang.blogspot.com



HORIBA

## **Light: Interference**

Look at just the electric field.

 $E = E_0 \sin(kx - \omega t)$  Oscillating electric field

 $E = E_0 \sin(kx - \omega t + \phi)$ 

Second electric field with phase shift



We will use interference to develop an intuitive understanding later.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

# Microscopic View of Scattering



Today we discuss the scattered photons that do not change wavelength. Also, we treat everything as a continuum.

About 1 scattered photon in a million has a wavelength shift that depends on the molecule. These extra photons are the basis for Raman spectroscopy. Horiba makes Raman instruments too!

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



A good approximation to develop intuition



Explore the future Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



#### Electric field at target from each position s: $\exp(ik(r_0 + s\sin(\theta)) - \omega t)$

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



We want to treat a hole because it is easy. We need a trick.



Babinet's principle: Scattering from an aperture is the same as scattering from its complement

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Need to figure area of each strip of aperture to integrate. For a slit it is easy (and we usually ignore length):





#### HORIBA

HORIBA

## **Fraunhofer Diffraction**



Now integrate over s from –R to R

$$\int_{-R}^{R} 2\sqrt{R^2 - s^2} \exp(iks\sin\theta) ds$$

dimensionless size parameter  $\alpha$  =  $\pi D$  /  $\lambda$ 

$$E = \frac{2E_0}{r} R^2 \frac{J_1(\alpha \sin \theta)}{\alpha \sin \theta}$$



Explore the future



### Fraunhofer Approximation

$$(S_1)^2 = (S_2)^2 = \alpha^4 \left[ \frac{J_1(\alpha \sin \Theta)}{\alpha \sin \Theta} \right]^2$$
$$I(\Theta) = \frac{I_0}{k^2 a^2} \alpha^4 \left[ \frac{J_1(\alpha \sin \Theta)}{\alpha \sin \Theta} \right]^2$$

dimensionless size parameter  $\alpha = \pi D/\lambda$ ;

 $J_1$  is the Bessel function of the first kind of order unity.

Assumptions:

a) all particles are much larger than the light wavelength (only scattering at the contour of the particle is considered; this also means that the same scattering pattern is obtained as for thin two-dimensional circular disks)

b) only scattering in the near-forward direction is considered (Q is small).

Limitation: (diameter at least about 40 times the wavelength of the light, or  $\alpha >>1$ )\* If  $\lambda$ =650nm (.65 µm), then 40 x .65 = 26 µm If the particle size is larger than about **26** µm, then the Fraunhofer approximation gives good results. Rephrased, results are insensitive to refractive index choice.

# Fraunhofer: Effect of Particle Size



Explore the future

# Diffraction Pattern: Large vs. Small Particles



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

- How many of you work with particles with sizes over 1 mm?
- How many of you work with particles with sizes over 25 microns and less than 1 mm?
- How many of you work with particles with sizes over 1 micron and less than 25 microns?
- How many of you work with particles with sizes less than 1 micron?



#### **Mie Scattering**



Explore the future Automotive Test Systems © 2012 HORIBA, Ltd. All rights reserved.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



# **Maxwell's Equations**

 Solve these and the universe's secrets are yours.
Notation from Bohren and Huffman <u>Absorption and Scattering of Light</u> <u>by Small Particles</u>

$$\nabla \bullet \mathbf{E} = \mathbf{0}$$

 $\nabla \bullet \mathbf{H} = \mathbf{0}$ 

Net charge is zero

No magnetic monopoles

$$\nabla \times \mathbf{E} = i \,\omega \mu \cdot \mathbf{H}$$

$$\nabla \times \mathbf{H} = -i\omega\varepsilon \cdot \mathbf{E}$$

Time varying magnetic field gives electric field and vice versa. No current flows through system.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



## **Boundary Conditions**

- Properties of particle different from surrounding medium
- Tangential components of E and H are continuous at boundary





# **Overview of Solution**

- Sum of incoming plane wave and scattered waves.
- We also get information on extinction and so on.





## **Mie Scattering**

$$I_{s}(m, x, \theta) = \frac{I_{0}}{2k^{2}r^{2}} \left( \left| S_{2} \right|^{2} + \left| S_{1} \right|^{2} \right)$$

Use an existing computer program for the calculations!

$$S_{1}(m, x, \theta) = \sum_{n=1}^{\infty} \frac{2n+1}{n(n+1)} \{a_{n}\pi_{n} + b_{n}\tau_{n}\}$$
$$S_{2}(m, x, \theta) = \sum_{n=1}^{\infty} \frac{2n+1}{n(n+1)} \{a_{n}\tau_{n} + b_{n}\pi_{n}\}$$

$$\tau_n = \frac{d}{d\theta} \Big( P_n^1(\cos\theta) \Big)$$

 $\pi_n = \frac{P_n^1(\cos\theta)}{\sin\theta}$ 

$$a_n = \frac{m\psi_n(mx)\psi_n'(x) - \psi_n(x)\psi_n'(mx)}{m\psi_n(mx)\xi_n'(x) - \xi_n(x)\psi_n'(mx)}$$

$$b_n = \frac{\psi_n(mx)\psi_n'(x) - m\psi_n(x)\psi_n'(mx)}{\psi_n(mx)\xi_n'(x) - m\xi_n(x)\psi_n'(mx)}$$

 $\xi, \psi$ : Ricatti-Bessel functions  $P_n^1:1^{st}$  order Legendre Functions

The equations are messy, but require just three inputs which are shown below. The nature of the inputs is important.

$$x = \pi D / \lambda$$

Decreasing wavelength is the same as increasing size. So, if you want to measure small particles, decrease wavelength so they "appear" bigger. That is, <u>get a blue light</u> <u>source for small particles</u>.



We need to know relative refractive index. As this goes to 1 there is no scattering.

Scattering Angle



## **Effect of Size**



As diameter increases, intensity (per particle) increases and location of first peak shifts to smaller angle.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

# Effect of RI: imaginary term



As imaginary term (absorption) increases location of first peak shifts to smaller angle.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



## **Effect of RI: Real Term**



#### It depends....

# Mixing Particles? Just Add



The result is the weighted sum of the scattering from each particle. Note how the first peak from the 2 micron particle is suppressed since it matches the valley in the 1 micron particle.

HORIBA

Explore the future

# Effect of Distribution Width



#### As distribution becomes wider, the peaks become less distinct.

### Iterations



#### HORIBA

HORIBA

![](_page_32_Picture_0.jpeg)

## **Pop Quiz**

As particle size increases:

- Peaks shift to smaller angles
- Peaks shift to larger angles
- •Nobody said there would be a quiz!

![](_page_32_Picture_8.jpeg)

![](_page_33_Picture_0.jpeg)

## **Pop Quiz**

As particle size increases:

- Peaks shift to smaller angles
- Peaks shift to larger angles
- •Nobody said there would be a quiz!

1/2 credit

Explore the future

![](_page_33_Picture_9.jpeg)

![](_page_34_Picture_0.jpeg)

#### **Comparison of Models**

![](_page_34_Figure_2.jpeg)

Light energy scattering patterns for an arbitrary detector configuration against particle size ( $\mu$ m) and scattering angle (°) for equal volumes of particles (Fraunhofer theory)

Light energy scattering patterns for an arbitrary detector configuration against particle size (µm) and scattering angle (°) for equal volumes of particles (Mie theory, latex particles RI 1,60 - 0,0i, in water RI 1,33)

#### Fraunhofer (left) vs. Mie (right)

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

# Comparison, Large Particles

![](_page_35_Figure_1.jpeg)

For large particles, match is good out to through several peaks.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

# Comparison, Small Particles

![](_page_36_Figure_1.jpeg)

For small particles, match is poor.

![](_page_37_Picture_0.jpeg)

#### **Practical Application: Glass Beads**

![](_page_37_Figure_2.jpeg)

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

![](_page_38_Picture_0.jpeg)

#### **Practical Application: CMP Slurry**

![](_page_38_Figure_2.jpeg)

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

# What about background (blank)?

Even pure liquids scatter. Also, there is often some constant background signal due to stray light. And, the detectors will always show a background signal. These effects are small, but they should be corrected.

The background is subtracted before applying any of analysis mentioned before.

As the background becomes smaller, the analysis becomes easier.

![](_page_39_Picture_6.jpeg)

#### HORIBA

#### LA-950 Optics

![](_page_40_Figure_2.jpeg)

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

# Short Wavelengths for Smallest Particles

![](_page_41_Figure_1.jpeg)

By using blue light source, we double the scattering effect of the particle. This leads to more sensitivity. This plot also tells you that you need to have the background stable to within 1% of the scattered signal to measure small particles accurately.

HORIBA

© 2012 HORIBA, Ltd. All rights reserved

Explore the future

![](_page_42_Picture_0.jpeg)

### Why 2 Wavelengths?

![](_page_42_Figure_2.jpeg)

Data from very small particles.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

![](_page_43_Picture_0.jpeg)

### Viewing Raw Data

![](_page_43_Figure_2.jpeg)

Explore the future

![](_page_44_Picture_0.jpeg)

# **During Data Collection**

#### Choose "Both Graphs"

![](_page_44_Figure_3.jpeg)

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

![](_page_45_Picture_0.jpeg)

### Conclusions

- It helps to know some of the theory to best use a laser diffraction particle size analyzer
  - Small particles wide angles
  - Large particles low angles
- Look at Intensity curves
- Big peaks in your background (blank) mean particles or bubbles
- Use Mie theory at all times (default whenever choosing an RI kernel other than Fraunhofer)

![](_page_45_Picture_8.jpeg)

![](_page_46_Picture_0.jpeg)

### Q&A

Ask a question at <a href="mailto:labinfo@horiba.com">labinfo@horiba.com</a>

Keep reading the monthly <u>HORIBA Particle</u> e-mail newsletter!

Visit the <u>Download Center</u> to find the video and slides from this webinar.

Jeff Bodycomb, Ph.D.

- P: 866-562-4698
- E: jeff.bodycomb@horiba.com

![](_page_46_Picture_8.jpeg)

![](_page_46_Picture_11.jpeg)