




Setting Attainable and Practical Particle Size Specifications


Mark Bumiller
HORIBA Scientific
mark.bumiller@horiba.com



Why Set Specifications?

- Product release (quality)
- Communicate grade to buyers
- Internal requirements
- Required in regulated industries (food & drugs)
- Will focus on drugs, but same concepts applicable to other industries

This is not New!



The grading of blackpowder: Earlier we mentioned the different sizes of gunpowder grains, and about how smaller grains will burn more quickly than the larger ones. The term “**grade**,” when applied to gunpowder refers to the grain size, and *not* to its quality. There are two separate categories of gunpowder grades; “C” and “F” **grade**. “C” **grade** is for cannons and large capacity explosive devices. A single “C” being the largest grain size with smaller sizes graded down as, “CC,” “CCC,” etc. **Powder** that is meant for small arms purposes uses the letter “F” to denote the grain size with a single “F” (or 1-FG for “1FG” size, “F” **grade** grain). The more “F’s” you see within the **powder grade** designation, the *smaller* will be the grain size. “FFFG,” for example, is very fine, almost a dust and was commonly used as a priming **powder** for the pans of flintlocks.

Examples (Sieves)

Retained on 8 Mesh	0.5	Maximum
Retained on 16 Mesh	2	Maximum
Through 200 Mesh	30	Maximum

Mesh	Micron
No. 20	841
No. 25	707
No. 30	595
No. 35	500
No. 40	420
No. 45	354
No. 50	297
No. 60	250
No. 70	210
No. 80	177
No. 100	149
No. 120	125
No. 140	105
No. 170	88
No. 200	74
No. 230	63
No. 270	53
No. 325	44
No. 400	37

The following convention is used to characterize particle size by mesh designation: a "+" before the sieve mesh indicates the particles are retained by the sieve; a "-" before the sieve mesh indicates the particles pass through the sieve; typically 90% or more of the particles will lie within the indicated range.

-4 +40 mesh, then 90% or more of the material will pass through a 4-mesh sieve **and** be retained by a 40-mesh sieve

-40 mesh, then 90% or more of the material will pass through a 40-mesh

Example (TiO₂)

Item No.	R-1	RM-1
TiO ₂ Content %	98.0	94.0
Whiteness		Better
Brightness	96.5	
pH Value	6.5-8.5	6.5-8.5
Mesh Residue 45 μ m max %	0.05	0.05
Tint-Reducing Power min	1700	95
Oil Absorption g/100g max	21	23
Volatile AT 105°C % max	0.5	0.80
Specific Receptivity Ω cm min	90	80
Water dispersibility, % min		90
Rutile Content min %	97.0	98.0
Application	Paints, Plastics, Rubber, Paper, Grocery	Coating, Plastic, Paper, Rubber, Printing inc.

Does this mean 0.06% fails?? They are the same numbers!



TYPICAL PROPERTIES	
TiO ₂ content	94%
Inorganic coating	Alumina
Organic treatment	Present
Crystal size	0.25μm
Specific gravity	4.05 g/cm ³
Loss at 105°C ⁽ⁱ⁾	0.6% max
Bulk density (tamped) ⁽ⁱⁱ⁾	1.3 g/cm ³
Oil absorption ⁽ⁱⁱⁱ⁾	17 cm ³ /100g pigment
Water demand ^(iv)	24 cm ³ /100g pigment
Durability	Durable
ISO 591 classification	R2
ATSM D476 designation	II, III, IV

Easier to understand, but:
How was it measured?
Dispersion?
Ultrasound?



Specification Definition*

- “a list of tests, references to analytical procedures, and appropriate acceptance criteria, which are numerical limits, ranges, or other criteria for the tests described”
- “chosen to confirm the quality of the drug substance and drug product rather than to establish full characterization, and should focus on those characteristics found to be useful in ensuring the safety and efficacy of the drug substance and drug product”

*ICH HARMONISED TRIPARTITE GUIDELINE, SPECIFICATIONS: TEST PROCEDURES AND ACCEPTANCE CRITERIA FOR NEW DRUG SUBSTANCES AND NEW DRUG PRODUCTS: CHEMICAL SUBSTANCES, Q6A, Current *Step 4* version, 6 October 1999

Pharmaceutical: FDA Guidance

ICH HARMONISED TRIPARTITE GUIDELINE

SPECIFICATIONS: TEST PROCEDURES AND ACCEPTANCE CRITERIA FOR NEW DRUG SUBSTANCES AND NEW DRUG PRODUCTS: CHEMICAL SUBSTANCES Q6A

3. GUIDELINES

3.1 Specifications: Definition and Justification

3.1.1 *Definition of Specifications*

A specification is defined as a list of tests, references to analytical procedures, and appropriate acceptance criteria which are numerical limits, ranges, or other criteria for the tests described. It establishes the set of criteria to which a new drug substance or new drug product should conform to be considered acceptable for its intended use.

FDA Guidance

3.3 Specific Tests / Criteria

3.3.1 New Drug Substances

b) *Particle size*: For some new drug substances intended for use in solid or suspension drug products, particle size can have a significant effect on dissolution rates, bioavailability, and / or stability. In such instances, testing for particle size distribution should be carried out using an appropriate procedure, and acceptance criteria should be provided.

Decision tree #3 provides additional guidance on when particle size testing should be considered.


Is the drug a solid oral dosage form or suspension?

↓ Yes


Is particle size critical to:
 Dissolution, solubility, bioavailability?
 Processability?
 Stability?
 Content uniformity?
 Maintaining appearance?

↓ Yes

Set Acceptance Criterion



FDA: Nanoparticles





Attachment A: Nanotechnology Product Evaluating Questions

1) This review contains new information added to the table below: _____ Yes _____ No Review date: _____
2) Are any nanoscale materials included in this application? (If yes, please proceed to the next questions.) Yes _____; No _____; Maybe (please specify) _____
3 a) What nanomaterial is included in the product? (Examples of this are listed as search terms in Attachment B.) _____
3 b) What is the source of the nanomaterial? _____
4) Is the nanomaterial a reformulation of a previously approved product? Yes _____ No _____
5) What is the nanomaterial functionality? Carrier _____; Excipient _____; Packaging _____; API _____; Other _____
6) Is the nanomaterial soluble (e.g., nanocrystal) or insoluble (e.g., gold nanoparticle) in an aqueous environment? Soluble _____; Insoluble _____
7) Was particle size or size range of the nanomaterial included in the application? Yes _____ (Complete 8); No _____ (Go to 9)
8) What is the reported particle size? Mean particle size _____; Size distribution _____; Other _____
9) Please indicate the reason(s) why the particle size or size range was not provided: _____ _____
10) What other properties of the nanomaterial were reported in the application (see Attachment E)? _____
11) List all methods used to characterize the nanomaterial. _____

Particle Size and Processability

- Milling/size reduction
- Mixing/blending
- Separation
- Filtration
- Granulation
- Homogenization
- Crystallization



Narrow particle size distributions minimize segregation problems during mixing – more homogeneous distribution of components in final product

Particle Size and Dissolution


Equation 1.

$$\frac{dX_s}{dt} = -\frac{3DX_0^{1/3}X_s^{2/3}}{\rho hr_0} \left(C_s - \frac{X_d}{V} \right)$$

X_s is the mass of solid drug (mg),
 t is time (minutes),
 D is the drug diffusivity (cm²/min),
 X_0 is the initial drug mass (mg),
 ρ is the drug density (mg/mL),
 h is the diffusion layer thickness (cm),
 r_0 is the initial particle radius (cm),
 C_s is the drug solubility (mg/mL),
 X_d is the mass of dissolved drug (mg),
 V is the volume of dissolution media (mL).



FIGURE 4

Ondansetron Dissolution as a Function of Particle Size Fractions at pH 6.8. Data are From Model Predictions (Solid Lines) and Data Collected in Dissolution Experiments (Data Points).



David R. Friend, PhD; Gregory E. Parry, PhD; T. Francis, PhD; Gary Kupperblatt, PhD; Suggy S. Chrai, PhD; and Gerald Slack, Mathematical Modeling of a Novel Controlled-Release Dosage Form Drug Delivery Technology,

Content Uniformity



Idealized concept:

- Powder particles compacted into tablet
- All particles are active ingredients
- Small particles = specified dose
- Large particle = over dose

Content Uniformity*




Fig. 1. Particle size distributions of Bantam-milled drug (○) and Jet-milled drug (●). Simulated lines were drawn using geometric mean particle sizes of 18.5 and 6.1 μm and geometric standard deviation of 1.7 and 1.6 for the Bantam and Jet-milled drug, respectively.

$$\text{mass}_{ni} = \frac{\text{mass}_i}{\sum_{i=1}^{100} \text{mass}_i} \times \text{dose}$$

The volume v_i and mass m_i of a single particle of radius was calculated as follows:

$$v_i = \frac{4}{3}\pi r_i^3 \tag{7}$$

$$m_i = v_i \rho \tag{8}$$


where ρ is the drug density.

Percent of intent	Number of unit doses within the given percent of intent range			
	Experimental		Simulated	
	Bantam	Jet	Bantam	Jet
70-74				
74-78				
78-82				
82-86				
86-90	1			
90-94	5			
94-98	13	5	483 249	
98-102	20	59	313 436	1000 000
102-106	16		145 337	
106-110	7		45 537	
110-114	1		9744	
114-118			2191	
118-122	1		447	
122-126			58	
126-130	1		1	

*Zhang, Y, Johnson, K, Effect of drug particle size on content uniformity of low-dose solid dosage forms, International Journal of Pharmaceutics 154 (1997) 179-183

Tablets

- Size of active ingredient effects content uniformity
- Size influences tablet hardness
- Size and shape effects packing
- Size and shape effect powder flow



Suspensions


- Dissolution and absorption
- Content uniformity
- Ability to stay in suspension
- Feel in mouth



Only API's? No, Also Excipients

- Particle size and physical characteristics critical in selection and performance

HPC Grade		SSL	SL	L	
Material	Lactose (g)	700	700	700	
	Corn starch (g)	300	300	300	
	8% HPC aqueous solution (g)	375	375	375	
Property of granule	Particle size distribution (%)	1400µ on	0.2	0.2	-
		500µ	0.2	0.2	0.4
		355µ	0.4	0.6	1.7
		250µ	1.3	2.1	8.4
		180µ	4.0	6.0	14.4
		150µ	8.0	9.6	15.6
		106µ	22.1	22.1	24.7
		75µ	30.3	26.9	19.6
		75µ under	33.6	32.3	15.0
	Bulky density (kg/cm ³)	Loose	90	95	130
packed		0.5	0.47	0.46	
Property of tablet	Hardness (kg)	14	14	12	
	Friction loss (%)	0.2	0.2	0.2	
	Disintegration time (min)	6	8	9	



Techniques and Specifications



Explore the future


Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA


© 2010 HORIBA, Ltd. All rights reserved.

DLS: Brownian Motion

- Particle is randomly diffusing from Brownian motion
 - Larger particles will diffuse more slowly
 - Smaller particles will diffuse more quickly
- Scatter light off this diffusing particle
- Measure the frequency shift of the signal



Two Approaches: Correlator or Power Spectrum



$$R_H = \frac{kT}{6\pi\eta D}$$

ISO 22412: PSA by DLS



- Correlation method or frequency analysis
- Measure at least 3 times
- Record average particle size \overline{X}_{DLS} & PI
- Repeatability better than 5%
- Check for sedimentation
- Dilution study, use most dilute if changes
- Verification: within 2% of 100nm standard, repeatability < 2%, PI < 0.2
 - All specs too tight for reality
 - Within 5% OK

DLS Specifications

- Only z average & PI defined in ISO
- Can convert intensity to volume distribution
 - Requires RI of sample
- Can dramatically alter reported values
- Typically used for industrial products where D10, D50 & D90 familiar

Laser Diffraction

Particle size 0.01 – 3000 μm



- Converts scattered light to particle size distribution
- Quick, repeatable
- Most common technique


Distribution Parameters

Symmetric Distribution
 Mean = Median = Mode


Mode - Peak of the distribution

Median - 50% Point
 (50% above - 50% below)


Mean – Weighted Average - Center of Gravity of Distribution




Other Data Points: D10, D50, D90



Other Data Points: D10, D50, D90



Distribution Parameters



Volume Mean Diameter

- D[4,3] which is often referred to as the Volume Mean Diameter [VMD]

$$D [4, 3] = \frac{\sum D_i^4 n_i}{\sum D_i^3 n_i}$$

Setting a D [4,3] specification will emphasize the presence of large particles

Mean Size

The frequency distribution is found using the arithmetical mean diameter, as shown in the formula below.


$$\text{Mean Diameter} = \frac{\sum \{q(J) \times X(J)\}}{\sum \{q(J)\}}$$

J : Particle Diameter Division Number

q(J) : Frequency Distribution Value (%)

X(J) : Jth Particle Diameter Range's Representative Diameter (μm).

D 4,3 Volume Mean




15% difference


Classic example: CMP slurries

Bimodal Distribution

Which numbers to use for specifications?
 D50 still an option, but some prefer finer details



Bimodal Distribution Result Details



Laser Diffraction Standards

- Standards good source for guidance to methods & specifications
 - ISO13320
 - Pharmaceutical: EP 2.9.31 & USP<429>
- All based on ISO standard
- Test Reproducibility
 - Don't believe anything unless it's reproducible
- Verify your system
 - On a regular basis using polydisperse standards

Reproducibility

- Prepare & measure sample 3 times
- Record D10, D50, D90
- Calculate average D10, D50, D90 and COV
 - ISO: COV < 3% at median x_{50}
 COV < 5% at x_{10} & x_{90}
 - EP/USP: COV < 10% at median x_{50}
 COV < 15% at x_{10} & x_{90}
 - Can double COV values when $D50 < 10 \mu\text{m}$
- Part of specification? Perhaps internal

Calculation Automation

From LA-950 Software

Select Summary Items

Item List

- Test or Assay Number
- Remarks 1
- Remarks 2
- Remarks 3
- Remarks 4
- Remarks 5
- Remarks 6
- Remarks 7
- Remarks 8
- Remarks 9
- Remarks 10

Summary Items

- Sample Name
- Material
- Source
- Lot Number
- D(v,0.1)
- D(v,0.5)
- D(v,0.9)

Add >> Delete

Clear Up Down

Font: MS Sans Serif Font Open

Orientation: Portrait Landscape

Show Summary Averages Show Summary Std. Dev.

Show Coefficient of variation(Relative Std. Dev.)

Validation
Specification: USP 429

	D(v,0.1) Range (± %)	D(v,0.5) Range (± %)	D(v,0.9) Range (± %)
D(v,0.5) >= 70µm	15	10	15
D(v,0.5) < 10µm	30	20	30

Save As Cancel OK

Summary Report

Export Summary Print Summary Edit Layout Best Fit Columns Hide Selected Exit


Sample Name	Material	Source	Lot	D(v,0.1)	D(v,0.5)	D(v,0.9)
Sample 4	PinnoThin TG Powde	Herbalife		0.052	0.052	0.052
Sample 4	PinnoThin TG Powde	Herbalife		0.052	0.052	0.052
Sample 4	PinnoThin TG Powde	Herbalife		0.052	0.052	0.052
Sample 4	PinnoThin TG Powde	Herbalife		0.045	0.045	0.045
Sample 4	PinnoThin TG Powde	Herbalife		0.045	0.045	0.045
Sample 4	PinnoThin TG Powde	Herbalife		0.045	0.045	0.045
Sample 4	PinnoThin TG Powde	Herbalife		0.045	0.045	0.045
Sample 4	PinnoThin TG Powde	Herbalife		0.040	0.040	0.040
Sample 4	PinnoThin TG Powde	Herbalife		0.039	0.039	0.039
Sample 4	PinnoThin TG Powde	Herbalife		0.040	0.040	0.040
Sample 4	PinnoThin TG Powde	Herbalife		0.048	0.048	0.048
Sample 4	PinnoThin TG Powde	Herbalife		0.048	0.048	0.048
Sample 4	PinnoThin TG Powde	Herbalife		0.048	0.048	0.048
Sample 4	PinnoThin TG Powde	Herbalife		0.048	0.048	0.048
Sample 4	PinnoThin TG Powde	Herbalife		0.045	0.045	0.045
Average				0.046	0.046	0.046
Std. Dev.				0.005	0.005	0.005
CV (%)				9.805	9.805	9.805
USP 429 (30.0, 20.0, 30.0)				PASSED	PASSED	PASSED

System Verification

- Use polydisperse standard
 - Whitehouse, NIST
- 3 independent measurements, calculate mean
- Accuracy:
 - X50 < 3% certified range of values
 - X10 & X90 < 5% certified range of values
- Repeatability
 - COV X50 < 3%
 - COV X10 & X90 < 5%

Calculation Automation

From LA-950 Software



Distribution Graph | Data Table | Result Data

Mean Size : 0.18408(μm)
Variance : 1.8988E-3(μm²)
Median Size : 0.17730(μm)
Mode Size : 0.1649(μm)
Std.Dev. : 0.0436(μm)
Chi Square : 4.162519
R Parameter : 3.7379E-1
Diameter on Cumulative % : (2)10.00 (%) - 0.1345(μm)
 : (9)90.00 (%) - 0.2450(μm)
Cumulative % on Diameter : (1)850.0 (μm) - 100.000(%)
 : (2)600.0 (μm) - 100.000(%)
 : (3)425.0 (μm) - 100.000(%)
 : (4)300.0 (μm) - 100.000(%)
 : (5)212.0 (μm) - 100.000(%)
 : (6)150.0 (μm) - 100.000(%)
 : (7)106.0 (μm) - 100.000(%)
 : (8)75.00 (μm) - 100.000(%)
 : (9)53.00 (μm) - 100.000(%)
 : (10)38.00 (μm) - 100.000(%)
Verification : 1.OK 4.3% [D(v,0.5) 0.170 (μm)(± 6.000%)]
 : 2.OK 3.5% [D(v,0.1) 0.130 (μm)(± 10.00%)]
 : 3.OK 6.5% [D(v,0.9) 0.230 (μm)(± 10.00%)]

Data Name	Graph Type	Transmittance(R)	Median Size	R Parameter
andy1		88.3(%)	0.17730(μm)	0.373795
200801181026014		81.1(%)	9.35329(μm)	0.069234
andy1		88.3(%)	0.17730(μm)	0.373795

Pharmaceutical Specifications

- What Information Should be Included?

Analytical procedures (sampling, dispersion, system suitability, etc.)


- Method Validation (precision, ruggedness, dispersion stability, robustness, etc.)
- Acceptance Criteria (upper and lower limits)
- Representative plots of particle size distribution measurements should be included well as the method validation report.

Justification

- Justification should be presented for each procedure and each acceptance criterion included
- Should refer to relevant development data, pharmacopoeial standards, test data for drug substances and drug products used in toxicology and clinical studies, and results from accelerated and long term stability studies, as appropriate.
- Additionally, a reasonable range of expected analytical and manufacturing variability should be considered.

Example Specifications

- | | |
|---|--|
| <ul style="list-style-type: none"> ■ Dv50 = 25 μm ± 15% <ul style="list-style-type: none"> ● 21.2 to 28.8 μm ■ Dv10 = 8 μm ± 25% <ul style="list-style-type: none"> ● 6 to 10 μm ■ Dv90 = 75 μm ± 25% <ul style="list-style-type: none"> ● 56.2 to 93.7 μm ■ Maybe D[4,3] | <p>Dv50 NMT 25 μm</p> <p>Dv10 NMT 8 μm</p> <p>Dv90 NMT 75 μm</p> |
|---|--|


 Looks busy


 Looks cleaner
 Sufficient?

Both better than: Average = 25

Example Specification from FDA Website

Tests	Acceptance criteria	Analytical procedure	Test results for Lot#15531
Appearance	A white, crystalline powder.	Visual	Complies
Identification A: IR B: UV	A. IR: Corresponds to RS B. UV: Absorptivities at xxx nm, do not differ by more than 3.0% from the reference standard.	USP<197M> USP<197U>	Complies Complies
Heavy metals	NMT 20 ppm	USP<231>	LT 20 ppm
Assay	98.0-102.0%	USP method	99.5%
Residual solvents	Methanol: NMT 3000 ppm Methylene Chloride: NMT 600 ppm Toluene NMT 890 ppm	USP <467>	300 ppm 150 ppm 80 ppm
Related Substances	Specified Impurities* RC 1: NMT 0.15% RC 2: NMT 0.25% RC 3: NMT 0.25 % Any unspecified impurity: NMT 0.10% (each) Total impurities: NMT 0.75%	method #41	LT 0.05% LT 0.05% 0.10% LT 0.05% 0.30 %
Polymorphic Form (XRD)	Ratio of peak at 2θ= xx to peak at 2θ=yy: LT 5%	method #47	LT 1%
Particle size (Laser Diffraction)	D90: NMT 30 μm D50: NMT 15 μm D10: NMT 5 μm	method #48	20 μm 10 μm 2.5 μm



*RC 1: Finest identified

Pharmaceutical Reference*

*John, E, How to Set Specifications for the Particle Size Distribution of a Drug Substance?, American Pharmaceutical Review, April 2009, 72-77

	x10	x50	x90
mean, n = 6	0.8µm	1.8µm	3.9µm
STDrel, n = 6	13%	9%	13%
min	0.7µm	1.7µm	3.5µm
max	0.9µm	1.9µm	4.8µm
DEVrel (min, max)	29%	12%	37%

Table 1. Characteristic values of the particle size distributions (PSD) of six batches of Drug Substance A in terms of mean values, relative standard deviations, min. max. values and relative deviations between them.

Includes manufacturing and analytical variability over 3 years

FDA Reference*

Regulatory Issues on Particle Size Specifications



Zhigang Sun, Ph.D.

Office of Generic Drugs
Center for Drug Evaluation and Research
Food and Drug Administration

Zhigang.Sun@FDA.HHS.GOV



Opinions expressed in this presentation are those of the speaker and do not necessarily reflect the views or policies of the FDA

45th Annual Pharmaceutical Technologies Arden Conference


February 4, 2010

1

*http://www.aapspharmaceutica.com/meetings/workshops/Arden/presentations/Regulatory_Issues_on_Particle_Size-_Sun.pdf

FDA Reference*

Acceptance Criteria for Laser Diffraction



Upper and lower limits of D10, D50, and D90 are established based upon prior knowledge or design of experiment (DOE) studies.

Inadequate LD Criteria


Spec. (Laser Diffraction)

D10: NMT 35 μm
 D50: NMT 75 μm
 D90: NMT 120 μm



Spec. (Laser Diffraction)


90% particles: NMT 100 μm
Testing results
 D90: 3.7 – 6.2 μm



One-sided limit is not adequate without scientific justification.


Specification on X or Y Axis

Errors on X & Y axis not the same



Specification with Error

Must tighten internal spec by lab error %
 Then product always within performance specification



http://www.spcpress.com/pdf/Manufacturing_Specification.pdf, By David Wheeler

Specification with Error

- Must tighten internal spec by lab error
- Therefore minimize lab error makes life easier
- How to minimize error?
 - Get sampling right
 - Structured method development
 - Eye on the goal: reproducibility

Conclusions

- Specifications based on product performance
- Tighten internal specification to include measurement error
- Report results in format created by instrument
 - Zave & PI intensity from DLS
 - Volume results from laser diffraction
- Standards provide help to set specifications
- Avoid Dv100 with laser diffraction
- Use D 4,3 if performance sensitive to small amount of large particles

<http://www.horiba.com/scientific/products/particle-characterization/>

To Learn More

<http://www.horiba.com/scientific/products/particle-characterization/>

HORIBA Scientific

United States

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | **Scientific** | All Segment Product Browser

Products | Markets & Industries | Service & Support | News & Events | About Us | Employment

Particle Characterization

Home > Scientific > Products > Particle Characterization

Particle Characterization

HORIBA designs, manufactures, and supplies state of the art particle characterization instruments.

Every instrument across the five business segments must meet stringent requirements before the HORIBA name is attached. The Particle Characterization group of analyzers has incorporated this principle into each new design since entering the business in 1979. Relentless innovation united with high performance to attain the ultimate goal: a new standard in usability.

Particle Characterization Products

HORIBA offers instruments for particle size, particle shape, zeta potential, and surface area analysis. Measurable particle size range is from 1 nanometer to 30 millimeters, at concentrations ranging from 1 ppm to 50 vol% with shape determination available starting at 1 micrometer. A range of analytical techniques are employed including laser diffraction (Mie Theory), dynamic light scattering, acoustic and electroacoustic spectroscopy, and dynamic and static image analysis. (measuring both particle size and shape information).

Request Information

Quick Request

Newsletter

Download Center

User login

Thank-you