



# The Importance of Refractive Index When using Laser Diffraction

Mark Bumiller  
[mark.bumiller@horiba.com](mailto:mark.bumiller@horiba.com)




Explore the future


**HORIBA**

© 2011 HORIBA, Ltd. All rights reserved.


# Fraunhofer Approximation



# Mie Theory






# Mie vs. Fraunhofer





**Figure A.3 -- Comparison of scattering patterns of non-absorbing particles according to Fraunhofer and Mie calculations ( $N_p = 1,59 - 0,0$ ;  $n_{\text{water}} = 1,33$ ; wavelength = 633 nm)**



# Practical Application of Theory: Mie vs. Fraunhofer for Glass Beads


Data Name	Graph Type	Refractive Index (R)
Standard Glass Beads Mie		STD-GLASSBEADS[STD-GLASSBEADS( 1.510 - 0.000i),
Standard Glass Beads Fraunhofer		Fraunhofer Kernel[Fraunhofer Kernel( 0.000 - 0.000i]




Graph Type	D(v,0.1)	D(v,0.5)	D(v,0.9)
	8.98783(μm)	13.47741(μm)	18.8536
	2.58072(μm)	3.62044(μm)	22.3174

# Practical Application of Theory: Mie vs. Fraunhofer for CMP Slurry


Data Name	Graph Type	Refractive Index (R)
CMP Slurry Mie	 —	2.20-0.0i[2.20-0.0i( 2.200 - 0.000i),Water( 1.333)]
CMP Slurry Fraunhofer	 —	Fraunhofer Kernel[Fraunhofer Kernel( 0.000 - 0.000i)]



# Influence of Particle Size on Angular Light Intensity Scattering Patterns



# Influence of Imaginary Parts of RI (Absorbancies)





# Refractive Index

- Refractive Index is defined by two components – Real and Imaginary

$$RI = n + ik$$

Where:

$N$  = real component, velocity of light in a vacuum/ velocity of light in the particle

$$= c/v_p$$

$c$  = speed of light in vacuum

$V_p$  = speed of light in particle (liquid, air)

$k$  = the extinction coefficient of the material

$$i = \sqrt{-1}$$

# RI: Real Component

- Real component from published tables or measured using Snell's Law

$$n \sin \theta = n' \sin \theta'$$


where:

$n$  = RI of air


$\theta$  = angle of incidence

$n'$  = RI of particle

$\theta'$  = angle of refraction



# n is Frequency Dependent



Remember the  
LA-950 has  
two light sources  
Red = 650 nm  
Blue = 405 nm

In this example  
1.53 vs. 1.51  
Not very important

Could matter for  
sub-micron pigments

Sellmeier equation

$$n^2(\lambda) = 1 + \frac{B_1\lambda^2}{\lambda^2 - C_1} + \frac{B_2\lambda^2}{\lambda^2 - C_2} + \frac{B_3\lambda^2}{\lambda^2 - C_3},$$

# RI: Imaginary Component

- Imaginary component ( $k$ ) in Mie equation is the extinction coefficient of the particle. Reduction of transmission of optical radiation caused by absorption and scattering of light.


$$k = ( \lambda / 4 \pi ) \alpha$$

Where:

$k$  = the extinction coefficient


$\alpha$  = the absorption coefficient

$\lambda$  = the wavelength of light used



- Absorption coefficient ( $\alpha$ ) = reciprocal distance light penetrates surface & be attenuated to  $1/e$  of original intensity, about 37%. Opaque particles have high coefficient.

# RI: Imaginary Component



- $1 / d = \alpha$
- $k = (\lambda / 4 \pi) \alpha$
- Example: for a value of  $k = 0.1$
- $k = (\lambda / 4 \pi) [\alpha] = 0.05 \times [1 / d] = 0.1 = k$
- $d = 0.05 / 0.1 = 0.5 \text{ micron}$

# RI: Imaginary Component

- For transparent particles use 0 for the imaginary component
- For slightly opaque materials use 0.01 or 0.1
- For opaque materials use 1.0 or higher
- Values can exceed 1.0 (see below\*)

## Complex Index of Refraction Look-up Utility

### Instructions

Choose a material from the drop down menu, enter the incident wavelength, and click on "Calculate." The material's index of refraction and extinction coefficient will appear below.

### Index of Refraction Explanation

The Index of Refraction of a material is the ratio of phase velocity of an electromagnetic wave in free space to the phase velocity in the material. The index of refraction of two materials can be used to predict how light will pass from one medium to the other. The Extinction Coefficient is the imaginary part of the index of refraction.

These optical constants are calculated via linear interpolation on the optical constant values reported in the book [Handbook of Optical Constants of Solids](#) by Edward D. Palik. The graphs and numerical data can be found on the [Tabulated Optical Constants](#) page.

Calculate

Material

Copper

Material  
Index:

1.6623006134

Material  
Extinction  
Coefficient

2.0176993865

Incident Wavelength:

650 (in nanometers)

\*<http://www.ee.byu.edu/photonics/opticalconstants.phtml>

# Refractive Index Effect

Most pronounced when:

- Particles are spherical
- Particles are transparent
- RI of particle is close to RI of fluid
- Particle size is close to wavelength of light source

Least pronounced when:

- Particles are not spherical
- Particles are opaque
- RI of particle is larger than RI of the fluid
- Particle size is much larger than wavelength of the light source


# Older HORIBA Software: RI Kernels

If  $n = 1.33$  (water) and  $n' = 1.60$  (particle)

Relative refractive index (RRI) :  $1.60/1.33 = 1.203$

Transparent particles:  $k = 0$

Kernel function: 120-000rri



Particles somewhat opaque: then  $k > 0$

KERNEL function: 120-020rri

If powder, relative index:

$n = 1.0$  (air) and  $n' = 1.60$  (particle), RRI:  $1.60/1.0 = 1.60$

Kernel function: 160-000rri or 160-020rri




# Practical Approach

- Measure sample once
- Change RI
- Does result change significantly?
- No: stop worrying
- Yes: keep working
- For LA-950 users: determine real component, vary imaginary component to minimize residual R value.



Oops, gave the plot away

# Abbe Refractometer

- Dissolve sample at different concentrations
- Plot conc. vs. RI
- Extrapolate to infinite concentration




# Becke Lines




Bright line is called the **Becke line** and will always occur closest to the substance with a higher refractive index

# Becke Line Test

As you move *away* from the thin section (raising the objective or lowering the stage), the Becke Line appears to move *into the material with greater refractive index*.



A particle that has greater refractive index than its surroundings will refract light inward like a crude lens.



A particle that has lower refractive index than its surroundings will refract light outward like a crude diverging lens.

# Becke Line Test

## Cargille Labs

Services to the sciences since 1924

### Refractive Index (Matching) Liquids

Also see our **Specialized Optical Coupling Liquids**

Cargille Refractive Index Liquids have become standard items in many laboratories as their applications have expanded from routine mineralogical identifications and quality control. New and broader uses in many more fields such as chemicals, engineering, medical, forensic, optics, instrumentation are continuously discovered. Many special requirements for specific applications have created a need for more technical data, new formulations, extended ranges, smaller increments and higher degrees of precision.

### Price Schedule

All Refractive Index Values are standardized at 5893 angstroms and 25°C.  
For Custom Made Liquids please fill in our **Optical Liquids Worksheet**  
All of the Refractive Index Liquid Sets are sold in 7.4cc bottles (1/4 fl oz)

**COMBINED SETS:  $n_D$  1.400-1.700**

The range of many minerals, most chemicals and virtually all biological materials are covered by the three Certified Series AA, A and B. These three Series are available as complete sets:

Cat. #		Price
18001	RF-1.....Full Set; Intervals 0.002; 151 liquids	\$1,865.00
18002	RF-1/2.....Half Set; Intervals 0.004; 76 liquids	1,025.00
18005	RF-1/5.....Fifth Set; Intervals 0.01; 31 liquids	415.00

Individual Series, fractions of Series or selected liquids in the Standard Group may also be purchased separately.

McCRONE  
ASSOCIATES

Services  
Techniques  
Case Studies  
Publications  
Courses  
Open Positions  
Contact Us



## WELCOME

[SEARCH](#) | [FEEDBACK](#) | [HOME](#)

McCrone Associates' team of scientists use state-of-the-art instrumentation to solve difficult and unique particle identification, materials characterization and analysis problems. Our extensive array of microscopy tools along with the knowledge, experience, technical skills and creative enthusiasm of our staff, provide an unequalled combination of analytical and problem solving capabilities.



Explore the future

**HORIBA**

© 2011 HORIBA, Ltd. All rights reserved.

# Luxpop.com

## Index of Refraction


- Return the **refractive index** of a substance at a given wavelength,  $\lambda$  (nm). Further information in addition to the index of refraction also may be given.
- Luxpop returns the absolute refractive index (i.e. with resp. to vacuum), unless stated otherwise. [Click here for more index of refraction terminology.](#)
- To facilitate search, select the input box and type in the first letter of the desired substance. [Contact](#) Luxpop to request more materials.
- See also our [long list of Index of Refraction Values \(A-Z\)...](#) for other materials.

Substance:   $\lambda$ :  nm temperature (certain substances only)  deg C

At a wavelength of 650 nm (1.908 eV), the index of refraction of ZnO is  
 $n = 1.98$ ,  $k = 0.000$ .

At a wavelength of 405 nm (3.061 eV), the index of refraction of ZnO is  
 $n = 2.23$ ,  $k = 0.010$ .

Note RI is dependent on wavelength of light.  
Can adjust RI for red & blue light, but only need to for small, absorbing particles.



# RI of Polymers

[http://www.texloc.com/closet/cl\\_refractiveindex.html](http://www.texloc.com/closet/cl_refractiveindex.html)

[Previous Page](#)

[Home Page](#)

[Search](#)

## TexLoc Refractive Index of Polymers

### NOTES:

Products supplied by TexLoc are designated in **bold green text**.

Water has a refractive index of 1.33.

PVDF	<b>Poly(vinylidene fluoride)</b>	<b>1.4200</b>
ECTFE	<b>Ethylene Chlorotrifluoroethylene</b>	<b>1.4470</b>
	Poly(trifluoroethyl methacrylate)	1.4370
	Poly(methyl octadecyl siloxane)	1.4430
	Poly(methyl hexyl siloxane)	1.4430
	Poly(methyl octyl siloxane)	1.4450
	Poly(isobutyl methacrylate)	1.4470
	Poly(vinyl isobutyl ether)	1.4507
	Poly(methyl hexadecyl siloxane)	1.4510
PEO	Poly(ethylene oxide)	1.4539
	Poly(vinyl ethyl ether)	1.4540
	Poly(methyl tetradecyl siloxane)	1.4550
	Poly(ethylene glycol mono-methyl ether)	1.4555
	Poly(vinyl n-butyl ether)	1.4563
PPOX	Poly(propylene oxide)	1.4570
	Poly(3-butoxypropylene oxide)	1.4580
	Poly(3-hexoxypropylene oxide)	1.4590
	Poly(ethylene glycol)	1.4590
	Poly(vinyl n-pentyl ether)	1.4590



# Google Search

refractive index magnesium stearate - Google Search - Windows Internet Explorer

http://www.google.com/search?q=refractive+index+magnesium+stearate&rlz=p.com.microsoft.\*:IE-SearchBox&ie=UTF-8&oe=UTF-8&sourceid=ie7&rlz=117GFRC\_en

File Edit View Favorites Tools Help

Links Customize Links

Google index magnesium stearate Search ABC Check AutoFill refractive index magnesium stearate Sign In

McAfee SiteAdvisor


refractive index magnesium stearate - Google Search

Web Images Maps News Video Gmail more

Google refractive index magnesium stearate Search Advanced Search Preferences

Web Results 1 - 10 of about 7,550 for refractive index magnesium stearate. (0.23 seconds)

24 results stored on your computer - Hide - About

 Diffraction Performance W... - Amount: 132 mg Magnesium Stearate Median (D50) March 31 Webinar.ppt - Sample Amount: 132 mg Magnesium Stearate Median (D50) 9.33

**MAGNESIUM STEARATE (Octadecanoic Acid, Magnesium Salt)** ✓  
Magnesium Distearate; Dibasic Magnesium Stearate; Magnesiumdistearat (German); Diestearato de ... REFRACTIVE INDEX. FLASH POINT. STABILITY. APPLICATIONS ...  
chemicaland21.com/industrialchem/organic/MAGNESIUM%20STEARATE.htm - 40k -  
Cached - Similar pages

**Resin having a large refractive index, lenses comprising the resin ...** ✓  
The present invention relates to a resin having a high refractive index that ..... calcium laurate, magnesium stearate, magnesium oleate, magnesium laurate, ...  
www.freepatentsonline.com/5084545.html - Similar pages  
by T Nagata - 1992 - Cited by 1 - Related articles - All 2 versions

**Plastic lenses having a high-refractive index, process for the ...** ✓  
Exemplary suitable metallic salts include zinc stearate, zinc oleate, ..... The refractive index, Abbe's number, mold release characteristics and appearance ...  
www.freepatentsonline.com/5594088.html - Similar pages  
by T Nagata - 1997 - Cited by 1 - Related articles - All 4 versions  
More results from www.freepatentsonline.com »

**STOCHEM - Specialty Chemical Distribution** ✓  
Size = 0.3 microns, pH = 7, Oil absorption = 43, refractive index = 1.56, ... Highly pulverized powder, Magnesium Stearate Modified, GE brightness = 87%, ...  
www.stochem.com/searchdb.asp?searchStr2=&supplier=511&page=2 - 17k -  
Cached - Similar pages

Sponsored Links

**Magnesium stearate** ✓  
Get the Answers You're Looking For.  
Magnesium stearate  
www.RightHealth.com/Nutrition





# Mixtures

$$n_{mix} = \sum_{i=1}^z n_i V_i$$

Where  $n_{mix}$  is the refractive index of the mixture,  $n_i$  is the refractive index of the  $i$ th component,  $V_i$  is the volume fraction of that component, and  $z$  is the total number of components. The real


$$n_{mix} = 1 + \rho_{mix} \sum_{i=1}^z m_i k_i$$

The term  $m_i$  is the molar fraction of the component in the mixture and  $k_i$  is its Gladstone-Dale constant. These



Larsen, E.S., Berman, H., *The Microscopic Determination of the Non-Opaque Minerals, Second Edition*, United States Department of the Interior, Geological Survey Bulletin 848, 1934, US Government Printing Office, Washington, DC.

# Mixtures



**Table 1**

	Median Size ( $\mu\text{m}$ )	
	Volume Number	
SEM	0.75	0.27
2.05-3.0i	0.64	0.24
2.35-3.0i	0.65	0.24

# Anisotropy

Also called birefringence

Many crystals have different RI as inspected through different planes.

Use average, consider flow alignment



I would rather eat French food

# Uniaxial vs. Biaxial Materials

Uniaxial materials, at 590 nm<sup>[3]</sup>

Material	$n_o$	$n_e$	$\Delta n$
beryl $\text{Be}_3\text{Al}_2(\text{SiO}_3)_6$	1.602	1.557	-0.045
calcite $\text{CaCO}_3$	1.658	1.486	-0.172
calomel $\text{Hg}_2\text{Cl}_2$	1.973	2.656	+0.683
ice $\text{H}_2\text{O}$	1.309	1.313	+0.004
lithium niobate $\text{LiNbO}_3$	2.272	2.187	-0.085
magnesium fluoride $\text{MgF}_2$	1.380	1.385	+0.006
quartz $\text{SiO}_2$	1.544	1.553	+0.009
ruby $\text{Al}_2\text{O}_3$	1.770	1.762	-0.008
rutile $\text{TiO}_2$	2.616	2.903	+0.287
peridot $(\text{Mg, Fe})_2\text{SiO}_4$	1.690	1.654	-0.036
sapphire $\text{Al}_2\text{O}_3$	1.768	1.760	-0.008
sodium nitrate $\text{NaNO}_3$	1.587	1.336	-0.251
tourmaline (complex silicate)	1.669	1.638	-0.031
zircon, high $\text{ZrSiO}_4$	1.960	2.015	+0.055
zircon, low $\text{ZrSiO}_4$	1.920	1.967	+0.047

Biaxial materials, at 590 nm<sup>[3]</sup>

Material	$n_\alpha$	$n_\beta$	$n_\gamma$
borax	1.447	1.469	1.472
epsom salt $\text{MgSO}_4 \cdot 7(\text{H}_2\text{O})$	1.433	1.455	1.461
mica, biotite	1.595	1.640	1.640
mica, muscovite	1.563	1.596	1.601
olivine $(\text{Mg, Fe})_2\text{SiO}_4$	1.640	1.660	1.680
perovskite $\text{CaTiO}_3$	2.300	2.340	2.380
topaz	1.618	1.620	1.627
ulexite	1.490	1.510	1.520

Biaxial birefringence, also known as trirefringence, describes an anisotropic material that has more than one axis of anisotropy


$$\Delta n = n_e - n_o$$

where  $n_o$  and  $n_e$  are the refractive indices for polarizations perpendicular (ordinary) and parallel (extraordinary) to the axis of anisotropy respectively.

3 Elert, Glenn. "Refraction". *The Physics Hypertextbook*. <http://hypertextbook.com/physics/waves/refraction/>.

# Scary Example: Is the Small Peak Real?


**Starch sample**  
**Low imaginary**



$$\text{INDEX} = 1.22 - 0.1i$$

# Scary Example: Is the Small Peak Real?


Starch sample  
High imaginary




INDEX = 1.23 - 4.13i

# Scary Example: Is the Small Peak Real?


Starch sample  
Zero imaginary



# Effect of RI: Cement



Fixed absorbance, vary real



Fixed real, vary absorbance



# Chi Square and R Parameter

$$\chi^2 = \sum \left\{ \frac{1}{\sigma_i^2} [y_i - y(x_i)]^2 \right\}$$

$$R = \frac{1}{N} \sum_{i=1}^N \left\{ \frac{1}{y(x_i)} |y_i - y(x_i)| \right\}$$





$y_i$  The measured scattered light at each channel (i) of the detector.

$y(x_i)$  The calculated scattered light at each channel (i) of the detector based on the chosen refractive index kernel and reported particle size distribution.

$\sigma_i$  The standard deviation of the scattered light intensity at each channel (i) of the detector. A larger  $\sigma_i$  indicates lower reliability of the signal on a given detector.

N The number of detectors used for the calculation

# Using R Value for i




Explore the future

© 2011 HORIBA, Ltd. All rights reserved.

HORIBA

# Changing RI



Explore the future

© 2011 HORIBA, Ltd. All rights reserved.

HORIBA

# Changing RI

**Measurement Dialog**

Load

Sample information

Sample Name: LUDOX

Material: Colloidal silica

Source: LUDOX

Lot Number: 01-01183

Test or Assay. Number: F0706U09-IT

Calculation Data Setting

☐ ActiveData

☒ Select Data in Memory

Select Data

☐ 200711201424025

☐ 200711201437026

☐ 200711201443027

☐ 200711261953029

☐ 200711271359036

☐ 200711271404037

☐ 200711271405038

☐ 200711271410039

☐ Red Iron Oxide 170 Max Wet - 01

Refractive Index

File Name: 1.45 1.0

Comment:

Form of Distribution

☐ Manual ☒ Auto

Condition Iteration Number 15

Distribution base

☒ Volume ☐ Area

☐ Length ☐ Numbers

Advanced

ReCalc

Cancel

**Select Kernel**

Folder: C:\Program Files\HORIBA\LA-950E\LAACQUISITION\Kernel

Select Folder

File Name	File Comment	Sample Name	Sample Comment	Sam...	Sam...	Sam...	Sam...	Dispersion Name	Dispersio...	Disp...	Disp...
1.33 1.0 in 1.385		1.33 1.0 in 1.385		1.3300	1.0000	1.3300	1.0000	Heptane	Heptane	1.3850	1.38
1.45 1.0 in 1.33		1.45 1.0		1.4500	1.0000	1.4500	1.0000	Water	Water	1.3330	1.33
1.45 1.0		1.45 1.0		1.4500	1.0000	1.4500	1.0000	Water	Water	1.3330	1.33
1.51 1.0 in 1.33		1.51 1.0		1.5100	1.0000	1.5100	1.0000	Water	Water	1.3330	1.33
1.55 1.0 in 1.33		1.55 1.0 in 1.33		1.5500	1.0000	1.5500	1.0000	Water	Water	1.3330	1.33
1.57 0		1.57 0		1.5700	0.0000	1.5700	0.0000	Water	Water	1.3330	1.33
1.59 0.1 in 1.378		1.59 0.1 in 1.378		1.5900	0.1000	1.5900	0.1000	Isopropanol	Isopropanol	1.3780	1.37
1.6 0.1 in 1.33		1.6 0.1 in 1.33		1.6000	0.1000	1.6000	0.1000	Water	Water	1.3330	1.33
1.60-0 in water		RI=1.60		1.6000	0.0000	1.6000	0.0000	Water	Water	1.3330	1.33
1.70-0.1 IPA		1.70-0.1i		1.7000	0.1000	1.7000	0.1000	Isopropanol	Isopropanol	1.3780	1.37
Alumina	water	Alumina	Alumina	1.6600	0.0000	1.6600	0.0000	Water	Water	1.3330	1.33
Aluminum	water	Aluminum	Aluminum	1.6000	5.4000	1.6000	5.4000	Water	Water	1.3330	1.33
Amber	water	Amber	Amber	1.5400	0.0000	1.5400	0.0000	Water	Water	1.3330	1.33
Antimony	water	Antimony	Antimony	3.2000	5.0000	3.2000	5.0000	Water	Water	1.3330	1.33
Asphalt	water	Asphalt	Asphalt	1.6300	0.0000	1.6300	0.0000	Water	Water	1.3330	1.33
Barium carbonate	water	Barium carbonate	Barium carbonate	1.6000	0.0000	1.6000	0.0000	Water	Water	1.3330	1.33
Barium fluochloride	water	Barium fluochloride	Barium fluochloride	1.6400	0.0000	1.6400	0.0000	Water	Water	1.3330	1.33
Barium fluoride	water	Barium fluoride	Barium fluoride	1.4700	0.0000	1.4700	0.0000	Water	Water	1.3330	1.33
Barium phosphate	water	Barium phosphate	Barium phosphate	1.6200	0.0000	1.6200	0.0000	Water	Water	1.3330	1.33
Barium sulfate	water	Barium sulfate	Barium sulfate	1.6200	0.0000	1.6200	0.0000	Water	Water	1.3330	1.33
Barium sulfide	water	Barium sulfide	Barium sulfide	2.1600	0.0000	2.1600	0.0000	Water	Water	1.3330	1.33
Barium yellow	water	Barium yellow	Barium yellow	1.6300	0.0000	1.6300	0.0000	Water	Water	1.3330	1.33
Cadmium sulfide	water	Cadmium sulfide	Cadmium sulfide	2.4200	0.0000	2.4200	0.0000	Water	Water	1.3330	1.33
Calcium alminate	water	Calcium alminate	Calcium alminate	1.7100	0.0000	1.7100	0.0000	Water	Water	1.3330	1.33
Calcium borate	water	Calcium borate	Calcium borate	1.6000	0.0000	1.6000	0.0000	Water	Water	1.3330	1.33
Calcium carbonate	water	Calcium carbonate	Calcium carbonate	1.5800	0.0000	1.5800	0.0000	Water	Water	1.3330	1.33
Canadian balsam	water	Canadian balsam	Canadian balsam	1.5200	0.0000	1.5200	0.0000	Water	Water	1.3330	1.33
Carbon	water	Carbon	Carbon	1.9200	0.0000	1.9200	0.0000	Water	Water	1.3330	1.33
Celunene	water	Celunene	Celunene	1.8400	0.0000	1.8400	0.0000	Water	Water	1.3330	1.33
Chrome green	water	Chrome green	Chrome green	2.4000	0.0000	2.4000	0.0000	Water	Water	1.3330	1.33
Chromium oxide	water	Chromium oxide	Chromium oxide	2.5000	0.0000	2.5000	0.0000	Water	Water	1.3330	1.33
Cobalt blue	water	Cobalt blue	Cobalt blue	1.7400	0.0000	1.7400	0.0000	Water	Water	1.3330	1.33
Cobalt green	water	Cobalt green	Cobalt green	1.9700	0.0000	1.9700	0.0000	Water	Water	1.3330	1.33

OK

Cancel

Explore the future

© 2011 HORIBA, Ltd. All rights reserved.

**HORIBA**

# Changing RI

Measurement Dialog

Load

Sample information

Sample Name  
LUDOX

Material  
Colloidal silica

Source  
LUDOX

Lot Number  
01-01183

Test or Assay. Number  
F0706U09-IT

Refractive Index

File Name: 1.45 1.0 Select

Comment: Create

Form of Distribution

☐ Manual ☒ Auto

Condition Iteration Number 15

Distribution base

☒ Volume ☐ Area

☐ Length ☐ Numbers

Advanced

Calculation Data Setting

☐ ActiveData


☒ Select Data in Memory

Select Data


- ☐ 200711201424025
- ☐ 200711201437026
- ☐ 200711201443027
- ☐ 200711261953029
- ☐ 200711271359036
- ☐ 200711271404037
- ☐ 200711271405038
- ☐ 200711271410039
- ☐ Red Iron Oxide 170 Max'Wet - 01

ReCalc

Cancel




# Using R Value for i



# Using R Value for i

Real component = 1.57 via Becke line

Vary imaginary component, minimize Chi square & R parameter



# Automation by Method Expert

## ■ Analytical conditions

Circulation Pump Speed

Circulation Pump Speed

Particle Concentration

Particle Concentration

Ultrasonic Treatment

Ultrasonic Treatment

Measurement Duration

Measurement Duration

## ■ Calculation conditions

Real Refractive Index Wizard

Real Refractive Index Wizard

Imaginary Refractive Index Wizard

Imaginary Refractive Index Wizard

View Method Expert webinar  
TE004 in [Download Center](#)



# Automated RI Computation

## ■ Real part study

- Need to fix imaginary part
- Set up to 5 real parts
- Software will compute all RI and display R parameter variation with RI selection

The screenshot displays a software interface for automated RI computation, organized into five steps:

- Step 1: Select measurement data for test**  
Options: ☒ Select Active Memory Data, ☐ Select DataFile [ ] [Select File]
- Step 2: Choose RI for liquid dispersant**  
Input: 1.333 [Open List]
- Step 3: Input RI imaginary component for test**  
Input: 0
- Step 4: Input RI real component for test**  
Test Value 1: 1.5, Test Value 2: 1.6, Test Value 3: 1.7, Test Value 4: 1.8, Test Value 5: 1.9
- Step 5: Push "Execute ..." button.**  
Text: This wizard is temporarily closed, and the test sequence is executed.  
Button: Execute Test Sequence >>


# Automated RI Computation

## Calculation Optimization


### Real Refractive Index Wizard -Result-



Graph Summary | Distributions | D90,D50,D10 | R Parameter



# Automated RI Computation



# Automated RI Computation

## ■ Imaginary part study


- Need to fix real part
- Set up to 5 imaginary parts
- Software will compute all RI and display R parameter variation with RI selection

The screenshot displays a multi-step wizard for automated RI computation. The steps are as follows:

- Step 1: Select measurement data for test**  
Options: ☒ Select Active Memory Data, ☐ Select DataFile [ ] [Select File]
- Step 2: Choose RI for liquid dispersant**  
Input: 1.33 [Open List]
- Step 3: Input RI real component for test**  
Input: 2.75
- Step 4: Input RI imaginary component for test**  
Test Value 1: 0.01, Test Value 2: 0.1, Test Value 3: 0.5, Test Value 4: 0.7, Test Value 5: 1
- Step 5: Push "Execute ..." button.**  
Text: This wizard is temporarily closed, and the test sequence is executed.  
Button: Execute Test Sequence >>

# Automated RI Computation

## ■ Imaginary study




# Summary

- Measure sample, recalculate w/different RI – see how important it is
- Use one of the described approaches to determine the real component
- Recalculate using different imaginary component
- Choose result that minimizes R parameter, but also check if result makes sense
- You wish you had Method Expert by your side

# It's All on the Website

- Register to visit Download Center
- Application Notes
- Technical Notes
- Past webinars
- Bibliography
- And soon... today's webinar
- mark.bumiller@horiba.com
- ian.treviranus@horiba.com

# Thank-you



Explore the future

© 2011 HORIBA, Ltd. All rights reserved.

**HORIBA**