




# Interpreting and Understanding Dynamic Light Scattering Data



Jeffrey Bodycomb, Ph.D.  
HORIBA Scientific  
[www.horiba.com/us/particle](http://www.horiba.com/us/particle)

# Outline

---

“Kitchen sink” talk. I cover a range of topics with an emphasis on conclusions rather than derivations.


- Introduction
- Looking at the ACF
- Effect of temperature and what it means to you.
- The Z-average
- Effect of concentration
- Hydrodynamic size

# What is Dynamic Light Scattering?

- Dynamic light scattering refers to measurement and interpretation of light scattering data on a microsecond time scale.
- Dynamic light scattering can be used to determine
  - Particle/molecular size
  - Size distribution
  - Relaxations in complex fluids

# Particle Diameter ( $\mu\text{m}$ )

**HORIBA**  
Scientific



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

**HORIBA**

© 2012 HORIBA, Ltd. All rights reserved.

# Use the Right Tool

- It is a struggle to use a micrometer with this steel block.



- It's easy with a steel rule (or calipers).




- The same applies for particles

# Particle Sizes over 1 micron

- Note that the upper limit of DLS is at 8 microns. This depends on particle density and other factors.
- Rule of thumb: If your particles are routinely bigger than 1 micron, consider laser diffraction.
- Particles that are too big often appear as 10~20 microns in DLS results. You will not be trapped by big particles that seem small.

# Brownian Motion

Particles in suspension undergo **Brownian motion** due to solvent molecule bombardment in random thermal motion.




## ■ Brownian Motion

- Random
- Related to Size
- Related to viscosity
- Related to temperature


# The SZ-100

- Single compact unit that performs size, zeta potential, and molecular weight measurements.



# Correlation Function


- Random fluctuations are interpreted in terms of the autocorrelation function (ACF).



$$C(\tau) = 1 + \beta \exp(-2\Gamma \tau)$$



# Smooth Autocorrelation Function

- These look good.
- As size increases, decay moves to longer times.
- Not enough data to decide if concentration is too high.



# Effect of Dust and Contamination

- These are examples of questionable data.
- Either the particle of interest is too large or there are too many large particle impurities.
- Filter samples or use software noise cut function.



# Gamma to Size

$$\Gamma = D_m q^2$$

$$q = \frac{4\pi n}{\lambda} \sin\left(\frac{\theta}{2}\right)$$

$$D_h = \frac{k_B T}{3\pi\eta(T)D_t}$$

$\Gamma$  decay constant  
 $D_t$  diffusion coefficient  
 $q$  scattering vector  
 $n$  refractive index  
 $\lambda$  wavelength  
 $\theta$  scattering angle  
 $D_h$  hydrodynamic diameter  
 $\eta$  viscosity  
 $k_B$  Boltzman's constant

Note effect of temperature!


# Effect of Temperature (and trends)

on(Sec)	Z-Average(nm)	PI	Peak1 Size (
30	94.3	0.132	
30	102.5	0.054	
30	110.9	0.043	
30	109.7	0.014	
30	108.7	0.056	
30	105.2	0.060	
0	6.9	0.044	

- Look at Z-average size. Data is OK.
- Is there a trend?
- Probably sample is not to temperature and viscosity value used in calculation is incorrect.
- In this case, I set up measurement conditions to force this event.

# Viscosity

- Get your viscosity correct.
- Choose the right liquid.
- Use viscosity at temperature of measurement.



# Why Z-average?

- Numerically stable
  - Result is not overly sensitive to noise in the data.
  - Important for QC work
- Described in detail in ISO standards
  - ISO-22412:2008 as  $\overline{d_{DLS}}$
  - ISO-13321:2004 as  $x_{PCS}$
- It is a useful measure of size since as average size increases, so does Z-average.

# Huh? What is the Z-average?

Determined by a mathematical method known as  
cumulants.

# The equations


$$\frac{1}{D_z} = \frac{\sum f_i D_i^6 P(\theta) D_i^{-1}}{\sum f_i D_i^6 P(\theta)}$$

Assume small angle compared to size so  $P(\theta)=1$

$$D_z = \frac{\sum f_i D_i^6}{\sum f_i D_i^5}$$


# Particle Size Distributions

- Particle size distributions can be plotted in several ways.
- Most often you see volume (mass) and number distributions




# A new distribution: Intensity

- Scattering goes by  $\sim d^6$
- The exponent works for small particles. We do the full calculations.




# Look at a linear scale

- These are lognormal distributions, so asymmetric.



# Need to use Diffusion Coefficient

- These are lognormal distributions, so asymmetric.



# Z-average

- As average size increases, so does Z-average.
- Tends to weight larger particles more than smaller (due to the physics of the measurement).

# Reproducibility

- PSL standards: you can get results better than 1%
- Don't expect this all the time.
- Expect 3~5%
  
- This is for Z-average.
  - Other average values (e.g., volume weighted mean) tend to vary more.
  - PI varies more.


# Comparing Techniques

- Always a good idea to check your results.
- Don't expect an exact match.
- Differences of 10~20% between laser diffraction (LD) and DLS are to be expected.

	D50 (vol. basis), nm LD	D50 (vol. basis), nm DLS	% diff in D50 (DLS/LD-1)	Z-avg. Diameter, nm DLS	% diff in size (DLS Z-avg./LD-1)
100 nm PSL	101	102.1	1.1	103.2	2.2
1 micron PSL	1059	1039.5	1.8	1112.7	5.1
E-1	129.8	146.6	12.9	118.3	-8.9
E-2	149.8	170.5	13.8	138.7	-7.4
E-3	110.0	100.2	-8.9	112.7	2.5
E-4	49.4	45.5	-7.9	32.4	-34.4
Ludox + 0.01 M KCl	36	21.2	-41.1	31.8	-11.7
Coffee Creamer wet	354	215.8	-39.1	336.9	-4.8




# Effect of Concentration

- Best is to make a study of measured size vs. concentration
- Note range of concentrations for which data is independent of concentration.
- Example below is “fake” data.



# What is Hydrodynamic Size?

- DLS gives the diameter of a sphere that moves (diffuses) the same way as your sample.



# Why DLS?

- Non-invasive measurement
- Fast results
- Requires only small quantities of sample
- Good for detecting trace amounts of aggregate
- Good technique for macro-molecular sizing



The SZ-100 from HORIBA



Questions?

[labinfo@horiba.com](mailto:labinfo@horiba.com)

[www.horiba.com/us/particle](http://www.horiba.com/us/particle)

[Jeff.Bodycomb@horiba.com](mailto:Jeff.Bodycomb@horiba.com)



**Thank you**

ありがとうございました

**Cảm ơn**

ขอบคุณครับ

谢谢

اشكر

**Gracias**

**Grazie**

Σας ευχαριστούμε

धन्यवाद

நன்றி

**Tacka dig**

**Danke**

**Merci**

**Obrigado**

감사합니다

Большое спасибо

おかしく おかし

*Omoshiro Okashiku*