



How to Select the Best Refractive Index



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA





Laser Diffraction Calculations

Importance of Refractive Index

Choosing Refractive Index

Comparing Methods for Choosing Refractive Index

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA





Laser Diffraction Calculations

Importance of Refractive Index

Choosing Refractive Index

Comparing Methods for Choosing Refractive Index

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

LA-950 Optics





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



When a Light beam Strikes a Particles



- Small particles require knowledge of optical properties:
 - Real Refractive Index (bending of light)
 - Imaginary Refractive Index (absorption of light within particle)
 - Refractive index values less significant for large particles
- Light must be collected over large range of angles

Diffraction Pattern





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Using Models to Interpret



Scattering

Scattering data typically cannot be inverted to find particle shape.

We use optical models to interpret data and understand our experiments.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



- There is no need to know all of the details.
- The LA-950 software handles all of the calculations with minimal intervention.





- Large particles -> Fraunhofer
 - More straightforward math
 - •Large, opaque particles
 - Use this to develop intuition
- All particle sizes -> Mie
 - Messy calculations
 - •All particle sizes





$$(S_1)^2 = (S_2)^2 = \alpha^4 \left[\frac{J_1(\alpha \sin \Theta)}{\alpha \sin \Theta} \right]^2$$
$$I(\Theta) = \frac{I_0}{k^2 a^2} \alpha^4 \left[\frac{J_1(\alpha \sin \Theta)}{\alpha \sin \Theta} \right]^2$$

dimensionless size parameter $\alpha = \pi D/\lambda$;

 J_1 is the Bessel function of the first kind of order unity.

Assumptions:

a) all particles are much larger than the light wavelength (only scattering at the contour of the particle is considered; this also means that the same scattering pattern is obtained as for thin two-dimensional circular disks)

b) only scattering in the near-forward direction is considered (Q is small).

Limitation: (diameter at least about 40 times the wavelength of the light, or $\alpha >>1$)* If λ =650nm (.65 µm), then 40 x .65 = 26 µm If the particle size is larger than about **26** µm, then the Fraunhofer approximation gives good results. Rephrased, results are insensitive to refractive index choice.

© 2013 HORIBA, Ltd. All rights reserved

Mie Scattering



$$I_{s}(m, x, \theta) = \frac{I_{0}}{2k^{2}r^{2}} \left(\left| S_{2} \right|^{2} + \left| S_{1} \right|^{2} \right)$$

Use an existing computer program for the calculations!

$$S_{1}(m, x, \theta) = \sum_{n=1}^{\infty} \frac{2n+1}{n(n+1)} \{a_{n}\pi_{n} + b_{n}\tau_{n}\}$$
$$S_{2}(m, x, \theta) = \sum_{n=1}^{\infty} \frac{2n+1}{n(n+1)} \{a_{n}\tau_{n} + b_{n}\pi_{n}\}$$

$$\tau_n = \frac{d}{d\theta} \Big(P_n^1(\cos\theta) \Big)$$

 $\pi_n = \frac{P_n^1(\cos\theta)}{\sin\theta}$

$$a_n = \frac{m\psi_n(mx)\psi_n'(x) - \psi_n(x)\psi_n'(mx)}{m\psi_n(mx)\xi_n'(x) - \xi_n(x)\psi_n'(mx)}$$

$$b_n = \frac{\psi_n(mx)\psi_n'(x) - m\psi_n(x)\psi_n'(mx)}{\psi_n(mx)\xi_n'(x) - m\xi_n(x)\psi_n'(mx)}$$

 ξ, ψ : Ricatti-Bessel functions $P_n^1:1^{st}$ order Legendre Functions



The equations are messy, but require just three inputs which are shown below. The nature of the inputs is important.

$$x = \pi D / \lambda$$

Decreasing wavelength is the same as increasing size. So, if you want to measure small particles, decrease wavelength so they "appear" bigger. That is, <u>get a blue light</u> <u>source for small particles</u>.



We need to know relative refractive index. As this goes to 1 there is no scattering.

Scattering Angle

Effect of Size





As diameter increases, intensity (per particle) increases and location of first peak shifts to smaller angle.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Mixing Particles? Just Add



HORIBA



The result is the weighted sum of the scattering from each particle. Note how the first peak from the 2 micron particle is suppressed since it matches the valley in the 1 micron particle.

Explore the future





Laser Diffraction Calculations

Importance of Refractive Index

Choosing Refractive Index

Comparing Methods for Choosing Refractive Index

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

What do we mean by RI?



- Optical properties of particle different from surrounding medium
- Note that intensity and wavelength of light changes in particle (typical dispersants do not show significant absorption)
- Wavelength changes are described by real component
- Intensity changes are described by imaginary component



Effect of RI: imaginary term





As imaginary term (absorption) increases location of first peak shifts to smaller angle.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Effect of RI: Real Term



It depends....

HORIBA

ORIBA

Scientific





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific







Refractive Index Effect



Most pronounced when:

- Particles are spherical
- Particles are transparent
- RI of particle is close to RI of fluid
- Particle size is close to wavelength of light source

Least pronounced when:

- Particles are not spherical
- Particles are opaque
- RI of particle is larger than RI of the fluid
- Particle size is much larger than wavelength of the light source





Laser Diffraction Calculations

Importance of Refractive Index

Choosing Refractive Index

Comparing Methods for Choosing Refractive Index

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Abbe Refractometer



- Dissolve sample at different concentrations
- Plot conc. vs. RI
- Extrapolate to infinite concentration









HORIBA





Bright line is called the *Becke line* and will always occur closest to the substance with a higher refractive index

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

© 2013 HORIBA, Ltd. All rights reserved.

Explore the future

Becke Line Test



As you move *away* from the thin section (raising the objective or lowering the stage), the Becke Line appears to move *into the material with greater refractive index.*



A particle that has greater refractive index than its surroundings will refract light inward like a crude lens.

A particle that has lower refractive index than its surroundings will refract light outward like a crude diverging lens.



Becke Line Test



Cargille Labs

Services to the sciences since 1924

Refractive Index (Matching) Liquids

Also see our Specialized Optical Coupling Liquids

Cargille Refractive Index Liquids have become standard items in many laboratories as their applications have expanded from routine mineralogical identifications and quality control. New and broader uses in many more fields such as chemicals, engineering, medical, forensic, optics, instrumentation are continuously discovered. Many special requirements for specific applications have created a need for more technical data, new formulations, extended ranges, smaller increments and higher degrees of precision.

Price Schedule

All Refractive Index Values are standardized at 5893 angstroms and 25°C. For Custom Made Liquids please fill in our Optical Liquids Worksheet All of the Refractive Index Liquid Sets are sold in 7.4cc bottles (1/4 fl oz)

COMBINED SETS: n_D 1.400-1.700

The range of many minerals, most chemicals and virtually all biological materials are covered by the three Certified Series AA, A and B. These three Series are available as complete sets:

Cat. #		Price
18001	RF-1Full Set; Intervals 0.002; 151 liquids	\$1,865.00
18002	RF-1/2Half Set; Intervals 0.004; 76 liquids	1,025.00
18005	RF-1/5Fifth Set; Intervals 0.01; 31 liquids	415.00
Individual	Series, fractions of Series or selected liquids in the	Standard

Group may also be purchased separately.





McCrone Associates' team of scientists use state-of-the-art instrumentation to solve difficult and unique particle identification, materials characterization and analysis problems. Our extensive array of microscopy tools along with the knowledge, experience, technical skills and creative enthusiasm of our staff, provide an unequaled combination of analytical and problem solving capabilities.





Services

Courses

Techniques

Case Studies

Publications

Contact Us

Open Positions

Luxpop.com



Index of Refraction

- Return the refractive index of a substance at a given wavelength, λ (nm). Further information in addition to the index of refraction also may be given.
- Luxpop returns the absolute refractive index (i.e. with resp. to vacuum), unless stated otherwise. Click here for more index of refraction terminology.
- To facilitate search, select the input box and type in the first letter of the desired substance. Contact Luxpop to request more materials.
- See also our long list of Index of Refraction Values (A-Z)... for other materials.



At a wavelength of 650 nm (1.908 eV), the index of refraction of ZnO is n = 1.98, k = 0.000.

At a wavelength of 405 nm (3.061 eV), the index of refraction of ZnO is n = 2.23, k = 0.010.

Note RI is dependent on wavelength of light. Can adjust RI for red & blue light, but only need to for small, absorbing particles.



Google Search



🖉 refractive index magnesium stearate - Google Search - Windows Internet Explorer	
COC + Khtp://www.google.com/search?q=refractive+index+magnesium+stearate&rls=p,com.microsoft:*:IE-SearchBox&ie=UTF-8&	pe=UTF-8&sourceid=ie7&rlz=117GFRC_en
File Edit View Favorites Tools Help Links Customize Links Customize Links Customize Links Customize Links Google Index magnesium stearate Search * Image: Steadwiser Search * Image: Steadwiser Tools Complex SiteAdvisor * * Image: Steadwiser * * Search * *	• 🌛 🔩 refractive 🔩 index 🔩 magnesium 🔩 stearate 🛛 🔌 • 🔵 Sign In •
😪 🏟 🎦 refractive index magnesium stearate - Google Search	🚹 👻 🔝 👻 🖶 Page 💌 🎯 Tools 🔹
Web Images Maps News Video Gmail more Google refractive index magnesium stearate Search Advanced Search Preferences Preferences Preferences	Sign in 💆
Web	Results 1 - 10 of about 7,550 for refractive index magnesium stearate. (0.23 seconds)
24 results stored on your computer - Hide - About	Sponsored Links Magnesium stearate Get the Answers You're Looking For. Magnesium stearate www.RightHealth.com/Nutrition
<u>Plasucienses naving a nigh-refractive index, process for the</u> Exemplary suitable metallic salts include zinc stearate, zinc oleate, The refractive index, Abbe's number, mold release characteristics and appearance www.freepatentsonline.com/5594088.html - <u>Similar pages</u> by T Nagata - 1997 - <u>Cited by 1 - Related articles - All 4 versions</u> More results from www.freepatentsonline.com » <u>STOCHEM - Specialty Chemical Distribution</u> Size = 0.3 microns, pH = 7, Oil absorption = 43, refractive index = 1.56, Highly pulverized powder, Magnesium Stearate Modified, GE brightness = 87%, www.stochem.com/searchdb.asp?searchStr2=&supplier=511&page=2 - 17k - Cached - Similar pages	

Using R Value for i







$$R = \frac{1}{N} \sum_{i=1}^{N} \frac{|y_i - y(x_i)|}{y(x_i)}$$

y_i is the measured scattering at detector I

 $y(x_i)$ is the scattering data at detector i calculated from the size distribution

N is the number of detector channels used in the calculation.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Changing RI



Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Changing RI



Measurement Dialog	×												
Calculation Data Setting C ActiveData C Select Data in Memory	Sample information Sample Name LUDOX Material Colloidal silica Source LUDOX Lot Number [01-01183 Test or Assay, Number F0705U09+IT Refractive Index File Name: 1.451.0 Comment:	Select Kernel Folder: C. VProgram Files/H File Name 1.33 1.0 in 1.385 1.45 1.0 in 1.33 1.45 1.0 in 1.33 1.55 1.0 in 1.33 1.55 0.0 in 1.33 1.50 0.1 in 1.33 1.50 0.1 in 1.33 1.50 0.0 in 1.37 1.60 0.1 in 1.33 1.80 0.1 in 1.34 1.80 0.1 in 1.34 1.80 0.1 in 1.34 Aumine Aluminum Amber Antimory Asphal Barium Rubchloride	DRIBALLA-950E LLAA(File Comment water water water water water water water water water water water water	Sample Name 1.33 1.0 in 1.385 1.45 1.0 1.51 1.0 1.51 1.0 1.55 1.0 in 1.33 1.57 0 1.50 0.1 in 1.378 1.6.0.1 in 1.378 1.6.0.1 in 1.378 1.70.0.1 i Alumina Alumina Alumina Alumina Alumina Bairum Catbonate Bairum Ruschloride	IISTION\Kernel Sample Name Sample Comment 1.33 1.0 in 1.395 1.45 1.0 1.45 1.0 1.51 1.0 1.55 1.0 in 1.33 1.57 0 1.55 1.0 in 1.33 1.57 0 1.59 0.1 in 1.33 1.61 1.0 1.61 1.0 1.51 0 1.57 0 1.59 0.1 in 1.33 1.61 0.1 1.61 1.1 Alumina Alumina Alumina Alumina Alumina Alumina Anber Amber Antimorry Antimory Asphalt Asphalt Barium carbonate Barium carbonate	Sam 1.3300 1.4500 1.5500 1.5500 1.5700 1.5900 1.6000	Sam 1.0000 1.0000 1.0000 1.0000 0.0000 0.1000 0.1000 0.1000 0.0000 0.0000 5.0000 0.0000 0.0000 0.0000	Selec Sam 1.3300 1.4500 1.5100 1.5500 1.5500 1.5500 1.6000 1.6000 1.6000 1.6000 1.6000 1.6400 1.6400	Sam. Dispersion Name 1.0000 Heptane 1.0000 Water 1.0000 Water 1.0000 Water 1.0000 Water 0.0000 Water	Dispersion Name Heptane Water Water Water Sopropanol Water Water Water Water Water Water Water Water Water Water Water	Dispersio Heptane Water Water Water Isopropanol Water Sater Water Water Water Water Water Water Water Water Water Water	Disp. Disp. Disp. 1.3850 1.38 1.3330 1.33 1.3330 1.33 1.3330 1.33 1.33	
Calculation Data Setting ActiveData Select Data in Memory Select Data 200711201424025 200711201437026 200711201437026 2007112714593303 200711271404037 200711271404037 200711271404037 200711271404038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 200711271410038 20071127141038 20071127141038 20071127141038 20071127141038 20071127141038 2007112714104 2007112714104 200711271404 200711271404 200711271404 200711271404 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 2007112014 200711271404 200711271404 2007112714 200711271404 200711271404 2007112714 2007112	Coloidal silica Source LUDOX Lot Number [01-01183 Test or Assay, Number F0706U03+1T Refractive Index File Name: [1.451.0 Comment: Condition Iteration Number 15 Distribution base © Volume © Area © Length © Numbers Advanced	1.45 1.0 in 1.33 1.45 1.0 in 1.33 1.55 1.0 in 1.33 1.55 0.0 in 1.33 1.55 0.0 in 1.33 1.50 0 in 1.50	water water	1.451.0 1.451.0 1.551.0 in 1.33 1.570 in 1.33 1.570 in 1.33 RI-1.60 1.700.1 in 1.378 Aluminum Aluminum Antimory Asphalt Barium cationate Barium Muchioride Barium Muchioride Barium Sulfate Barium sulfate Calcium brate Calcium Bariate Calcium Calcium Bariate Calciu	Alumina Aluminum Amber Antimory Asphalt Barium rathonate Barium fluochloide Barium phosphate Barium sulfide Barium sulfide Barium sulfide Barium sulfide Barium sulfide Calcium animate Calcium animate Calciu	1.4500 1.4500 1.5500 1.5500 1.5500 1.5500 1.5500 1.5600 1.5400 1.5400 1.5400 1.5400 1.5400 1.5400 1.5400 1.5400 1.5400 1.5200 1.		1.4500 1.5100 1.5700 1.5700 1.5700 1.6000 1.80000 1.80000 1.80000 1.80000 1.80000 1.80000 1.80000 1.80000000000	1.0000 1.0000 0.00000 0.00000 0.00000 0.000000	Water Water Water Water Isopropanol Water	Water Water Water Water Isopropanol Water Isopropanol Water	1.3330 1.33 1.3300 1.33 1.	
	alc Cancel												





HORIBA

Using R Value for i



Real component = 1.57 via Becke line Vary imaginary component, minimize Chi square & R parameter



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Automation by Method Expert



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Real part study

- Need to fix imaginary part
- Set up to 5 real parts
- Software will compute all RI and display R parameter variation with RI selection

Step 1: Select measure Select Active M	ment data for test lemory Data C Select [DataFile Select File			
Step 2: Choose RI for lie	quid dispersant	<u>Step 3</u> : Input RI imaginary component for test			
1.333	Open Lis	t 0			
Step 4: Input RI real con	nponent for test	Step 5: Push "Execute" button			
Test Value 1: 1.5	Test Value 4: 1.8	This wizard is temporarily closed,			
Test Value 2: 1.6	Test Value 5: 1.9	and the test sequence is executed.			
Test Value 3: 17	-	Execute Test Sequence >>			

Automated RI Computation



Calculation Optimization Real Refractive Index Wizard -Result-Graph Summary Distributions D90,D50,D10 R Parameter Distribution Graph 16 Real_RI-1.80- 0.00i Real_RI-1.70-0.00i Real RI-1.60-0.00i q(96) 10-Real RI-1.50-0.00i 0-. C F T T 0.01 100 1000 3000 0.1 Diameter(micro_m) D90.D50.D10 R Parameter 20 **D90** 0.134 Dia(micro_m) Parametei D50 D10 0.120-Q, 3 0.114 ********** (. 1.80- 0.001 1.50- 0.001 1.80- 0.00i 1.50- 0.001 Refractive Index Refractive Index

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Automated RI Computation





Scientific

Summary



- Measure sample, recalculate w/different RI see how important it is
- Use one of the described approaches to determine the real component
- Recalculate using different imaginary component
- Choose result that minimizes R parameter, but also check if result makes sense
- You wish you had Method Expert by your side





Laser Diffraction Calculations

Importance of Refractive Index

Choosing Refractive Index

Comparing Methods for Choosing Refractive Index

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Study on TiO₂



Look up real refractive index, use R parameter to find imaginary.



Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Effect of RI on R parameter





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Effect of RI on measured BA Scientific **D50** 10⁰ 0.38 Calculated Value of D50 0.5 Imaginary refractive index, dimensionless 0.00 V V V V 0.38 10⁻¹ - 0.4 10⁻² -0.44 - 0.3 0.2 10⁻³ 0.1 0.46 0.42 0.4 10⁻⁴ 2.8 3.0 2.4 2.6 3.2 3.4 3.6 Real refractive index, dimensionless

Effect of RI on measured



D10



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Effect of RI on measured







Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Results



Model	D10, microns	Difference from correct value	D50, microns	Difference from correct value	D90, microns	Difference from correct value
2.7 - 0.001i (correct value)	0.30		0.42		0.80	
2.6 - 0.001i	0.30	0%	0.45	7.1%	0.78	2.5%
Fraunhofer	0.41	37%	0.57	36%	0.81	1.2%
Minimize R parameter (3.25 - 0.1i)	0.26	13%	0.37	12%	0.72	10%

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





- Use the Mie model when evaluating laser diffraction data.
- Search for literature for real and imaginary refractive index values or measure your sample yourself.
- If literature values are not available, use the data to estimate values, it is better than guessing or using the Fraunhofer model.
- Once you choose a refractive index value, stick with it.



Ask a question at labinfo@horiba.com

Keep reading the monthly <u>HORIBA Particle</u> e-mail newsletter!

Visit the <u>Download Center</u> to find the video and slides from this webinar.

Jeff Bodycomb, Ph.D.

- P: 866-562-4698
- E: jeff.bodycomb@horiba.com









