

Trusting Your Zeta Potential Results



lan Treviranus ian.treviranus@horiba.com www.horiba.com/particle

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Intro to zeta potential

Factors which build trust/distrust

What to look for before measuring

Evaluating measurements





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2013 HORIBA, Ltd. All rights reserved

mar i medicar i Gerniconu

- Instrument is fully functional
- Suitable sample for technology
 - Counters are classic example
- Stay current with free webinars
 - Particle size essentials, method development, sampling & dispersion advice, DLS and zeta potential technology reviews







General concerns



- What is zeta potential? (webinar TE013)
- Dilution, DILUTION, DILUTION!
 - Zeta potential result interpretation
- 🎻 🗖 Assess data quality
 - How can z.p. be used in conjunction with particle size (larger topic, IEP, etc.)





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Featured technologies



LA-950



Dynamic Light Scattering & Zeta Potential

CAMSIZER & CAMSIZER XT

Dynamic Image Analysis

PSA300

Static Image Analysis

SA-9600

Flowing Gas BET Surface Area



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

SZ-100: Dynamic Light Scattering

- Particle size: 0.3 nm 8 µm
- Zeta potential: -200 +200 mV
- Molecular weight: 1x10³ 2x10⁷ Da
- Patented ultra long-life graphite electrodes
- Lowest total cost of ownership
- Optional autotitrator

nano partica SZ-100 series



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Intro to zeta potential

Factors which build trust/distrust

What to look for before measuring

Evaluating measurements





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Zeta potential is the charge on a particle at the shear plane.



How do Surfaces Acquire Charge?

Ionization of surface groups





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

How do Surfaces Acquire Charge?

Specific adsorption of ions, e.g. ionic surfactants





© 2013 HORIBA, Ltd. All rights reserved

Why Zeta Potential?



- Good way of evaluating electrostatic stabilization of suspensions
- Can use to predict interactions





Emulsion IEP Study: Stability



© 2013 HORIBA, Ltd. All rights reserved.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Intro to zeta potential

Factors which build trust/distrust

What to look for before measuring

Evaluating measurements





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Factors creating trust



- Excellent repeatability, reproducibility
- Stable results over time
- Long lifetime parts
- Robust method (avoid finicky results)
- A bit of training to know what results to expect when conditions change



Precision = Confidence



- Real life situation: historic, reproducible data is trusted data... though it may not be accurate data
- Is the data actionable?





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Factors creating <u>dis</u>trust



- Poor precision
 - Resort to "cherry pick" data
 - Is it a regular practice to discard 10%+ of z.p. results?
- Trending results
 - Indicates changing conditions
- Replacing cell too often
 - •Ten measurements?
- Lack of technical support, customer care



HORIBA

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Question



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Intro to zeta potential

Factors which build trust/distrust

What to look for before measuring

Evaluating measurements





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Practical Tips



- Visual inspection
 - Floaters?
- Minimum concentration
 - Concentrate with a centrifuge
- Maximum concentration
 - Dilute using supernatant (or similarly close)
 - Take care to match:
 - Salt concentration
 - pH
 - Surfactant concentration



Dilution



- Avoid dilution if possible

 simplest approach

 Don't dilute with DI water
 - •No ions, changes surface chemistry & ZP
- Best: equilibrium dilution with same liquid as sample, but with no particles
 - Use supernatant after sedimentation or centrifugation



Otherwise, dilute with 10 mM KCL solution







Question



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Eyeballing it



Samples up to little bit cloudy are OK.

Pick samples that are clear or just slightly cloudy.





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Intro to zeta potential

Factors which build trust/distrust

What to look for before measuring

Evaluating measurements





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



- Have existing method, does it produce "good" data?
- Stress the method
 - Measure at multiple pH values
 - Measure at multiple concentrations
- Always! collect multiple samplings
 - Reproducibility > challenge than repeatability
 - Evaluate with COV of chosen metric(s)







Question



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Good zeta potential data?



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Look at frequency data





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Look at frequency data





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

New ISO Standards



www.iso.org

ISO	International Organization for Standardization	International Standards for Business, Go	Search >>			
Home	Products Stan	dards development News and media	About ISO	For ISO Members	FAQs	Fr Ru 🔛 ISO Store



🖉 Standards under development





- 8.2.1 reference Materials
- Be sufficiently homogeneous and stable for the stated time and temperature range
- The accepted electrophoretic mobility value was obtained by several operators and rigorously proven.
- The material should be well documented in terms of sampling procedure, dilution, if required, and measurement protocol.



ISO13099: Verification



8.2.2 Repeatability

- Prepare sample following provided procedure
- Measure <u>same portion</u> three times
- Pass if mean value CV <10%</p>
 - •Assuming $2 \times 10-8 \text{ m}^2/\text{V} \cdot \text{s}$
- Note: expect most customers to use zeta potential values



Calculate COV



Liposome sample

6 repeats, 30 sec measurement, 5 sec delay

Zeta Potential	-					
Legend	No.	Date		Measurement Type	Sample Name	Zeta Potential (Mean)(mV)
	1	Friday, December 14, 2012	8:51:04 AM 2	Zeta Potential	FSS Liposome	-43.4
	2	Friday, December 14, 2012	8:52:20 AM 2	Zeta Potential	FSS Liposome	-44.2
	3	Friday, December 14, 2012	8:53:35 AM 2	Zeta Potential	FSS Liposome	-46.3
	4	Friday, December 14, 2012	8:54:51 AM 2	Zeta Potential	FSS Liposome	-46.1
	5	Friday, December 14, 2012	8:56:06 AM 2	Zeta Potential	FSS Liposome	-44.4
	6	Friday, December 14, 2012	8:57:39 AM 2	Zeta Potential	FSS Liposome	-41.7
Average	<u>(</u>				COV = 3	-44.3
S. D.					COV - J	1.7
30 30 (1) 25 (1) 25 (1) 25 (1) 25 (1) 10 (1)	-100	50 100 150 200 Zeta Deteriol (m/)				TRUST



- 8.2.3 Intermediate Precision
- Same procedure as 8.2.2 but using <u>different</u> <u>portions</u> of sample
- Pass if Repeatability; CV <15%</p>
 - •Assuming $2 \times 10-8 \text{ m}^2/\text{V} \cdot \text{s}$
- Doubt this is done often
- In pharmaceutical industry intermediate precision implies multiple systems, multiple operators, different days



Explore the future

Calculate COV



Liposome sample, two samplings

12 repeats, 30 sec measurement, 5 sec delay

Zeta Potential	-				
Legend	No.	Date	Measurement Type	Sample Name	Zeta Potential (Mean)(mV)
	1	Friday, December 14, 2012 8:51:04 AM	Zeta Potential	FSS Liposome	-43.4
	2	Friday, December 14, 2012 8:52:20 AM	Zeta Potential	FSS Liposome	-44.2
	3	Friday, December 14, 2012 8:53:35 AM	Zeta Potential	FSS Liposome	-46.3
	4	Friday, December 14, 2012 8:54:51 AM	Zeta Potential	FSS Liposome	-46.1
	5	Friday, December 14, 2012 8:56:06 AM	Zeta Potential	FSS Liposome	-44.4
<u> </u>	6	Friday, December 14, 2012 8:57:39 AM	Zeta Potential	FSS Liposome	-41.7
	7	Friday, December 14, 2012 9:45:24 AM	Zeta Potential	FSS Liposome	-49.1
	8	Friday, December 14, 2012 9:46:41 AM	Zeta Potential	FSS Liposome	-50.2
	9	Friday, December 14, 2012 9:47:58 AM	Zeta Potential	FSS Liposome	-51.3
	10	Friday, December 14, 2012 9:49:15 AM	Zeta Potential	FSS Liposome	-48.7
	11	Friday, December 14, 2012 9:50:33 AM	Zeta Potential	FSS Liposome	-51.5
	12	Friday, December 14, 2012 9:52:08 AM	Zeta Potential	FSS Liposome	-48.6
Average					-47.1
S. D.					0.0% 3.2



Zeta Potential Cells





Gold coated electrodes (ruined)

Carbon coated electrodes



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA





Rinse cell between 1st and 2nd series.

Data taken from poster submission to Techconnect Nanotech 2013



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Soft Nanoparticle



HORIBA



Data taken from Technical Note 167 available at www.horiba.com/particle

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



- 8.2.4 **Trueness** (accuracy)
- Prepare, measure same portion 3 times
- Pass if mean value within 10% of published electrophoretic mobility value, assuming > 2 × 10-8 m2/V•s.
- Note: when calculating pass/fail criteria it is OK to include the uncertainty of the sample
- Note: most customers use the zeta potential value, not the mobility



Instrument validation



Ideally, run traceable standard reference material(s)



NIST SRM 1980 is candidate, but not terribly popular

Can run internal reference and/or transfer "standard"



Explore the future

NIST SRM1980 on SZ-100



Pass criteria: Mean within 10% of reference value COV<10%

Reference value 2.53 µm·cm/V·s ± 0.12 µm·cm/V·s

Upper limit: (2.53 + 0.12)*1.1 = 2.92

Lower limit: (2.53 - 0.12) * 0.9 = 2.17

Pass



Sources of Error



- Contamination from previous sample
- Poor sample preparation
- Inappropriate sample
- Inappropriate liquid medium
- Poor temperature stabilization
- Condensation on the illuminated surfaces
- Too large a potential applied







- Particles, fingerprints or scratches on the optical surfaces
- Incorrect entry of parameters by the operator
- Air bubbles
- Cell coating damage
- Inappropriate theory for calculating zeta-potential from the measured electrophoretic mobility



Explore the future

Summary



- Zeta potential is less well understood
- Look at the sample before measurement
- (If possible) Know about the sample's chemistry before measurement
- Try to maintain sample integrity
- Look beyond the final z.p. result
- When all else fails... minimize the sample amount ③





Intro to zeta potential

Factors which build trust/distrust

What to look for before measuring

Evaluating measurements





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA





www.horiba.com/particle





Ian Treviranus Product Line Manager ian.treviranus@horiba.com

- Talk to us, ask questions
 labinfo@horiba.com
- Receive news of updates
- View application & technical notes (170+), webinars (70+), white papers.