Why the University of Washington chose the HORIBA LA-950 Laser Scattering Particle Size Distribution Analyzer

Presented by Michelle Katz and Tiffany Tang

COLLEGE OF ENGINEERING





- > About us
- > Why we chose the HORIBA LA-950
- > Intro to Laser Scattering Particle Size Distribution Analysis (LSPA)
- > How we use the HORIBA LA-950
 - Training users
 - Maintenance and Operation
- > Complimentary Instruments
- > Example analysis using tin oxide
- > Concluding Remarks

OUR ROLE AT UW

- Graduate students in Materials Science & Engineering (MSE)
- > Lab Assistants in the MSE User Facility
 - Instruments and equipment available for research and teaching
 - We maintain and train students/faculty to use these resources
- > Previously taught undergraduate classes involving the HORIBA LA-950



Tiffany



Michelle

OUR CRITERIA FOR CHOOSING AN INSTRUMENT

- > Analysis software that is free to students
- > Technical support after installation
- > Representative of current state of the art
- > Cost effective

THE HORIBA LA-950



SCIENCE OF LSPA & THE HORIBA LA-950

- > Particles in water or solvent are hit with light
- > Uses static light scattering (a.k.a. laser diffraction)
- > Calculates particle size using angle and intensity of scattered light
 - Large particles: higher intensity and smaller angle
 - This calculation is dependent on refractive index
- > HORIBA LA-950 particle size range: nm to mm

IMPORTANCE OF REFRACTIVE INDEX

- > Refractive index of many materials is already in software
 - If not, can be found in literature
- > Size calculation uses Mie Scattering Theory
 - Requires user to input refractive index
- > Refractive index is irrelevant for particles larger than ~30 microns
 - The intensity of refracted light is much smaller than reflected and diffracted light

HOW WE USE THE HORIBA LA-950

> Part of undergraduate instructional curriculum

- Students measure size of polymer spheres they learn to synthesize
- Students measure size of ceramic materials they process in class
- Students get hands-on experience from day 1

> Used by research groups across UW

- Example: used by Geology Department to measure rock samples

Students can use for independent projects at no cost

THE HORIBA LA-950 IS FASTER THAN...

- > Sieving
- > Gravity Settling
- > Photo-sedimentation
- > Electrical sensing zone techniques

The LA-950 is intuitive to use and has better repeatability.

TRAINING USERS

> Mandatory 2-part training session & written test

- 1st session: we show them how to use instrument & software
- 2nd session: they test their sample and demonstrate competence
- Short written test to prove thorough understanding of SOP

> We test samples for students occasionally

- Usually for one-time guests or short-term projects

OPERATION GUIDELINES

- > Analysis software available for free download
- > Reservation system: CORAL
 - Manages training, scheduling and user fees (for funded-research & industry)
- > SOP & step-by-step instructions provided for reference
- > Users may not troubleshoot, calibrate, or adjust the instrument in any way
 - Must contact us or our manager for assistance
- > Non-hazardous materials only

MANAGING MULTIPLE SAMPLES

- > Rinse cycle to clean sample holder
- > Change suspension fluid if necessary
- > Update the refractive index
- > Take a new background measurement

COMMON CAUSES OF BAD DATA

- > Surfactants can cause micro-bubbles
- > Some materials are very prone to agglomeration
- > Agitation can break particles into smaller pieces

PREVENTING AGGLOMERATION

- > Surfactants increase the surface area of the solvent, but can also cause bubbles
- > Agitation, Circulation, and Ultrasonic settings on LA-950 keep samples moving
- > Also a *de-bubble* function

LIMITATIONS TO LSPA

> Some materials do not work very well with LSPA

- Unstable nanoparticles
- Magnetic particles

> There is a minimum amount of material needed

- Especially a concern for small-batch synthesis
- Specific amount of sample dependent on material properties; can be tricky to predict in advance

TECHNIQUES TO COMPLIPENT LSPA

> X-ray Diffraction (XRD)

- Measures crystallite size (not particle size)
- Provides insight into crystal structure

> Optical Microscopy

- Provides insight into shape of large particles

> Scanning electron microscopy (SEM)

- Provides insight into shape of small particles
- Can confirm elemental composition

XRD AND SEM



 \leftarrow XRD

 $SEM \rightarrow$



EXAMPLE: TIN (II) OXIDE WITH LSPA

2013.03.19 14:56:45

University of Washington-MSE



Sample Name Data Name Transmittance(R) Transmittance(B) Circulation Speed Agitation Speed Ultra Sonic	: Tin Oxide : Average of 15): 95.3(%) : 96.5(%) d: 15 : 5 : 00:30 (7)	Median Size Mean Size Std.Dev. Geo.Mean Size Geo.Std.Dev. Span Diameter on Cumulative %	6	28.60883(µm) 31.67864(µm) 14.3838(µm) 29.1207(µm) 1.4945(µm) 1.0832 (2)10.00 (%)- 17.8674(µr (9)90.00 (%)- 48.8557(µr	Form of Distribution Refractive Index (R) Refractive Index (B) Distribution Base Span	: Auto : SnO[Stannous : SnO[Stannous : Volume : YES (10.00-	Oxide(1.840 - Oxide(1.840 - 90.00)	0.000i),SnOi 0.000i),SnOi
---	--	---	---	--	---	--	---	------------------------------





EXAMPLE: TIN (IV) OXIDE WITH LSPA





EXAMPLE: TIN (IV) OXIDE WITH LSPA

2013 03 20 17:49:30

University of Washington-MSE

Laser Scattering Particle Size Distribution Analyzer LA-950



TIN (II) OXIDE SEEN WITH SEM





TIN (IV) OXIDE SEEN WITH SEM





CONFIRMATION WITH XRD

- > XRD can determine crystal structure
- > Scherrer's equation estimates crystallite size
 - Uses full-width half-max of spectra peaks
- > Crystallite size is different than particle size

CONFIRMATION WITH XRD



FINAL REMARKS

- > HORIBA LA-950 is user-friendly for all experience levels
- > It is a multi-purpose tool to enhance education and research
 - Much faster and more reliable than many other methods
- > Pairs well with other characterization methods
 - XRD , OM, SEM

HORIBA LA-950 vs. LA-960



Sept. 2007 – Dec. 2014 0.01 - 3000µm for wet 0.1-3000µm for dry



July, 2013 0.01 - 3000µm for wet 0.1-5000µm for dry Improved user interface

CELL OPTIONS

Standard cell: up to 3000µm



Cell option: up to 5000µm

Specific Dry cell for size higher than 3mm

Free fall nozzle inside, for larger opening

Labinfo@horiba.com