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January, 2017

© 2017 HORIBA Advanced Techno Co., Ltd.

Conformable standards
This equipment conforms to the following standards:
RoHS: EN50581
9. Monitoring and control instruments

Authorised representative in EU
HORIBA UK Limited
2 Dalston Gardens, Stanmore, Middx HA7 1BQ, UK

The sodium ion-selective electrode uses a special glass that responds selectively to the sodium ion for its ion sensitive membrane. It may be used to measure the sodium ion concentration in water, food, boiler water, etc. Be sure to read the operating instructions carefully before usage.

1 Contents of the Package

Electrode main unit1
Instruction manual1

2 Specifications

Concentration range	2.3mg/L ~ 23,000mg/L Ag ⁺ (10 ⁻⁴ ~ 1mol/L-10mol/L)
Temperature range	0 ~ 60°C
pH range ^{(*)1}	greater than pH4.5 max.
Slope(mV/activity decade) ^{(*)2}	59mV
Response time T ₉₀ ^{(*)3}	Less than 30 seconds
Material used for the body	Glass
Supported reference electrodes	Double junction 1512A-06T 2565A-10T
Reference electrode	Inner tube 3.33mol/L KCl (248.3g/L KCl)
Type of internal filling solution	Outer tube 0.1mol/L KCl (7.5g/L KCl)

(*)1 This is the permitted pH range for the sample liquid ion concentration value given inside the brackets.
(The pH range varies slightly depending on the ion concentration).
(*)2 This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C).
(*)3 The response time when the concentration is changed from 10⁻⁴ mol/L to 10⁻⁵ mol/L.

3 The Electrode


Figure 1.

4 Electrode Characteristics

1. Calibration characteristics

Figure 2 shows the relationship between the concentration of the sodium ion and the output electrode potential.


Figure 2.

2. Influence of pH

The sodium ion electrode is affected by the pH in acidic regions. Figure 3 shows the relationship between the pH and the difference in output electrode potential for the sodium ion electrode.


Figure 3.

3. Influence of coexisting ions

The sodium ion electrode is influenced by various other ions. The influence of the coexisting ions can be shown in the following equation. K is the selectivity constant and is used as a standard for quantitatively indicating the influence of the coexisting ions.

$$E = E_0 + \frac{2.303RT}{F} \log \left(a_{Na^+} + \frac{1}{K} a_x \right) \quad \text{--- ②}$$

Where a_{Na^+} is the sodium ion concentration, a_x is the concentration of interfering ions. The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the test sample increases.

The selectivity constant K for various ions is listed in Table-1.

Table-1. Selectivity constants

Ca ²⁺	500
NH ₄ ⁺	20
K ⁺ , Li ⁺	10

4. Influence of temperature

The calibration curve of the electrode varies with temperature. You must maintain the standard solution and the sample solution at the same temperature during measurement.

5 Preparation Before Measurement

1. Preparation of the standard solution

Place 58.40 g of commercially available guaranteed reagent grade sodium chloride and 1 mL of 1 mol/L sodium hydroxide in a 1 L flask. Add enough pure water to make 1 L. This solution is your 1 mol/L Na⁺ (23,000 mg/L Na⁺) standard solution. Dilute this to produce 10⁻⁴ mol/L Na⁺ (2.3 mg/L Na⁺). Use pure water with a relative resistance of less than 0.5 μS/cm.

Use a plastic container for less than 10⁻³ mol/L Na⁺ standard solutions.

2. The reference electrode

Use the double-junction reference electrode (2565A-10C) in combination with the sodium ion electrode.

3. Before measurement

If the glass electrode has been left for a long time in a dry condition, soak it for a day and a night (at least more than 12 hours) in 0.01 to 1 mol/L sodium chloride solution before usage.

6 Measuring

1. Equipment

Assemble the electrode and ion meter (potentiometer) as shown in figure 4.


Figure 4.

- Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.
- Because of heat generated from the stirrer, the temperature of the solution is slowly rising, causing a drift in the electrode potential. If this is the case, lay a sheet of adiabatic material such as styrofoam between the stirrer and beaker in order to prevent temperature changes.
- The stirring may generate noise in the electrode potential. If this happens, and you are using an ion meter with a smoothing function, enable the smoothing function.
- Please use the recorder if necessary.

2. Operation

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- Direct potentiometric Method
- Known Addition Method
- Known Subtraction Method
- Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) Direct Potentiometric Method.

(a) Direct potentiometric Method

With the Direct Potentiometric Method, first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. Then measure the electrode potential in the sample solution, and arrive at the quantification value using the calibration curve.

3. Precautions

- Do not leave the sodium ion electrode soaking in pure water for long periods. It may reduce the response speed of the electrode.
- A special glass is used as the ion response element of the ion sensitive membrane in the sodium ion electrode. Do not use in high temperature high alkaline conditions. It may damage the glass ion sensitive membrane.
- The glass ion sensitive membrane is made of extremely thin glass and is very vulnerable to strong impact. Take appropriate care when handling the electrode.
- The electrode connector requires a high degree of insulation. Take care not to touch this part

7 Maintenance

- When storing the sodium ion electrode for long periods, pour 0.01 to 1 mol/L sodium chloride solution in the protection cap and attach it to the glass ion sensitive membrane section before storage.
- If the end of the electrode is very dirty and cannot be cleaned with pure water, soak some gauze in an appropriate organic solvent such as acetone or alcohol and use it to clean the electrode.
- If there is a drastic reduction in the response speed or reproducibility of the electrode, perform the operation described in item (2) of the maintenance procedures, then wash the end of the electrode with pure water, and soak it for a day in 1 mol/L sodium chloride solution. Check the response.