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The copper ion-selective electrode may be used to measure the copper ion concentration in any aqueous solution including rivers, lakes and marshes, or used in plant management in the plating, metals or cosmetics industries. In order to ensure correct usage of the ion electrode, be sure to read the operating instructions carefully before usage.

1 Contents of the Package

Copper Ion-selective Electrode 1
 Water proof abrasive paper (#1200) 1
 Instruction Manual 1

2 Specifications

Concentration range	0.06mg/L ~ 6,350mg/L Cu ²⁺ (10 ⁻⁶ ~ 10 ⁻¹ mol/L Cu ²⁺)	
pH range ^{※1}	pH2 ~ 6 (6.35mg/L Cu ²⁺)	
Temperature range	0 ~ 50°C	
Slope(mV/activity decade) ^{※2}	29mV	
Response time T90 ^{※3}	Less than 10 seconds	
Material used for the body	Epoxy resin	
Supported reference electrodes 8006-06T 8006-10C	Double junction referense electrode 2535A-06T 2565A-10T	
Reference electrode Type of internal filling solution	Inner tube	3.33 mol/L KCl : Cat. No. 300 (248.3g/L KCl)
	Outer tube	1mol/L KNO3 (100g/L KNO3)

(※1)

This is the permitted pH range for the sample liquid ion concentration value given inside the brackets. (The pH range varies slightly depending on the ion concentration).

(※2)

This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C.)

(※3)

The response time when the concentration is changed from 10⁻⁴ mol/L to 10⁻² mol/L.

3 The Electrode


Fig.1

4 Electrode Characteristics

1. Detection characteristic

Figure 2 shows the relationship between the copper ion concentration and the electrode potential.


Fig.2

2. The influence of pH

In alkaline regions, the copper ion becomes a hydroxide that is difficult to dissolve. The electrode must be used in a neutral to acidic solution. Although the pH level that gives rise to hydroxides varies with the copper concentration and the temperature, it is close to pH6.

Figure 3 shows the relationship between the pH and the difference in output electrode potential for the copper ion electrode.


Fig.3

3. Influence of coexisting ions

The copper ion electrode is influenced by various other ions. The influence of the coexisting ions can be shown in the following equation. K is the selectivity constant and is used as a standard for quantitatively indicating the influence of the coexisting ions.

$$E = E_0 - \frac{2.303RT}{2F} \log \left(a_{Cu^{2+}} + \frac{1}{K} a_x^{\frac{2}{n}} \right)$$

Where $a_{Cu^{2+}}$ is the cyanide ion concentration a_x is the concentration of interfering ions, n is the charge number of interfering ions. The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the test sample increases.

The selectivity constant K for various ions is listed in Table-1.

Table-1. Selectivity constants

Interfering Ion	Selectivity constant (Limit capacity)
Ni ²⁺ , Na ⁺	1000
Fe ²⁺	0.1

The copper ion forms a difficult to dissolve precipitate with the halogen ions (Cl⁻, Br⁻, I⁻). Note that if these anions exist in the test sample, there will be a reduction in the copper ion concentration reading.

When the concentration of coexisting ions exceeds the selectivity constant, try using a masking agent.

(Masking agents that could be considered are a chelating reagent, a precipitating agent, or complex compound.)

4. Influence of temperature

The calibration curve of the electrode varies with temperature. You must maintain the standard solution and the sample solution at the same temperature during measurement.

5 Preparations Before Measurement

1. Preparation of the standard solution

Accurately weigh out 3.928g of commercially available guaranteed reagent grade copper sulfate ($\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$) and dissolve it in pure (deionized) water. Add enough pure water to make 1 L. This solution is your 1000mg/L Cu^{2+} standard solution. Dilute 100mL, 10mL, and 1mL of this solution with pure water to produce 100mg/L, 10mg/L and 1 mg/L Cu^{2+} standard solutions.

Add 0.1mol/L of potassium nitrate (10g/L KNO_3) as a supporting electrolyte to each of these solutions. In addition, add the same concentration of supporting electrolyte as the standard solution to the test sample to obtain an even higher level of accuracy in your measurements.

2. Preparation of the reference electrode

A reference electrode is necessary for measurement. A double junction electrode with a sleeve type of liquid junction is used. Fill the outer tube with 1mol/L potassium nitrate as the internal solution. An accuracy level of $100 \pm 2\text{g/L}$ is adequate.

3. Electrode conditioning

Prepare the electrode before taking measurements by soaking it for approximately 10minutes in the 100mg/L Cu^{2+} standard solution.

6 Measuring

1. Equipment

Assemble the electrode and ion meter (potentiometer) as shown in figure 4.


Fig.4

- Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.
- Because of heat generated from the stirrer, the temperature of the solution is slowly rising, causing a drift in the electrode potential. If this is the case, lay a sheet of adiabatic material such as styrofoam between the stirrer and the beaker in order to prevent temperature changes.
- The stirring may generate noise in the electrode potential. If this happens, and you are using an ion meter with a smoothing function, enable the smoothing function.
- Use a recorder if required.

2. Preparation of the test sample

Add enough commercially available guaranteed reagent grade potassium nitrate to the test sample to make 0.1mol/L (10 g/L KNO_3).

3. Operation

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- (a) Direct potentiometric Method
- (b) Known Addition Method
- (c) Known Subtraction Method
- (d) Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) Direct Potentiometric Method.

With the Direct Potentiometric Method, you first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. You then measure the electrode potential in the sample solution, and arrive at the quantification value using the calibration curve.

7 Precautions

1. Do not rub or hit the ion sensitive membrane of the electrode. It may damage the membrane and degrade its special characteristics.
2. If direct sunlight falls on the ion sensitive membrane, the light may cause changes in the electrode potential. Use a brown beaker when using the equipment under such conditions.
3. Before soaking the electrode in a liquid with a different concentration, be sure to wash it first with deionized water.
4. The electrode connector requires a high degree of insulation. Take care not to touch this part with dirty hands or dampen it.

8 Troubleshooting

1. Check whether the problem is in the ion meter (potentiometer).
2. If the response from the electrode is slow, or its sensitivity has degraded, take the following steps to revive the ion sensitive membrane.
 - Drop a few drops of water onto the water proof abrasive paper included in the package. Lightly polish the ion sensitive membrane by holding the electrode vertically over the paper and rubbing it in circles. Afterwards, soak the electrode for approximately 10 minutes in 100mg/L Ca^{2+} standard solution before taking measurements.

9 Storage

- When storing the electrode for long periods, wash the ion sensitive membrane on the electrode with deionized water, and wipe it with some gauze. Fit the rubber cap included with the package over the response membrane and store the electrode in a cool dark place with a dry atmosphere.

Main Office/Factory: 2-Banchi, Kisshouinmiya-no-higashi-machi Minami-ku, Kyoto 601
Tel: (075) 313-8121
Tokyo Sales Office: 1-7-8 Higashi-kanda Chiyoda-ku, Tokyo 101
Tel: (03) 3861-8231