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The fluoride ion-selective electrode may be used to measure the fluoride ion concentration in any aqueous solution as well as used as a sensor for an automatic fluorine gas analyzer in the solution extracted from the atmosphere or furnaces of iron foundries, aluminum refineries, ceramics and cement factories.

The electrode uses a tip type ion electrode for easy replacement of the ion sensitive membrane. (Only 8010-10C electrode)

In order to ensure correct usage of the ion electrode, be sure to read the operating instructions carefully before usage.

1 Contents of the Package

Fluoride Ion-selective Electrode 1
Instruction Manual 1

2 Specifications

Concentration range	0.02mg/L ~ 19,000mg/L F ⁻ (10 ⁻⁶ ~ 1mol/L F ⁻)
pH range* ¹	pH4 ~ 10 (20mg/L F ⁻)
Temperature range	0 ~ 50°C
Slope(mV/activity decade)* ²	- 59mV
Response time T90* ³	Less than 5 seconds
Tip type electrode model (apply to 8010-10C)	#7661
Material used for the body	Poly vinyl chloride resin
Supported reference electrodes 8010-06T 8010-10C	2080A-06T or 2535A-06T 2060A-10T or 2565A-10T
Reference electrode Type of filling solution	3.33 mol/L KCl : Cat.No.#300 (248.3 g/L KCl)

(※1)

This is the permitted pH range for the sample liquid ion concentration value given inside the brackets. (The pH range varies slightly depending on the ion concentration).

(※2)

This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C.)

(※3)

The response time when the concentration is changed from 10⁻⁴ mol/L to 10⁻² mol/L.

3 The Electrode


Fig.1

4 Replacing the electrode

• Replacement method

Attachment and removal of the tip electrode from the main unit is via a single screw. Grip the rollette part of the tip electrode with your fingers and rotate it.

5 Electrode Characteristics

1. Detection characteristic

Figure 2 shows the relationship between the chloride ion concentration and the electrode potential.


Fig. 2

2. The influence of pH

The fluoride ion electrode is affected by the pH in strongly acidic or highly alkaline regions. The influence of pH varies with the fluoride ion concentration.

Figure 3 shows the relationship between the pH and the difference in output electrode potential for the fluoride ion electrode.


Fig.3

3. Influence of coexisting ions

The fluoride ion electrode is not affected by any anions other than the OH⁻ anion. The influence of the coexisting OH⁻ anion can be shown in the following equation. K is the selectivity constant and is used as a standard for quantitatively indicating the influence of the coexisting ions.

$$E = E_0 - \frac{2.303RT}{F} \log \left(a_{F^-} + \frac{1}{K} a_{OH^-} \right)$$

Where a_{F^-} is the fluoride ion concentration a_{OH^-} is the concentration of interfering ions. The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the test sample increases.

The selectivity constant K for various OH⁻ anion is equal to 10.

4. Influence of temperature

The calibration curve of the electrode varies with temperature. You must maintain the standard solution and the sample solution at the same temperature during measurement.

6 Preparations

1. Preparations Before Measurement

Accurately weigh out 3.058 g of commercially available guaranteed reagent grade potassium fluoride and dissolve it in pure (deionized) water. Add enough pure water to make 1L. This solution is your 1000mg/L F^- standard solution. Dilute 100 mL, 10mL, and 1mL of this solution with pure water to produce 100mg/L, 10mg/L and 1mg/L F^- standard solutions.

Add 0.1mol/L of potassium nitrate (10 g/L KNO_3) as a supporting electrolyte to each of these solutions. In addition, add the same concentration of supporting electrolyte as the standard solution to the test sample to obtain an even higher level of accuracy in your measurements.

Transfer the standard solutions to sealed plastic containers for storage.

2. Electrode conditioning

Prepare the electrode before taking measurements by soaking it for approximately 10 minutes in the pure (deionized) water.

- Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.
- Because of heat generated from the stirrer, the temperature of the solution is slowly rising, causing a drift in the electrode potential. If this is the case, lay a sheet of adiabatic material such as styrofoam between the stirrer and the beaker in order to prevent temperature changes.
- The stirring may generate noise in the electrode potential. If this happens, and you are using an ion meter with a smoothing function, enable the smoothing function.
- Use a recorder if required.

2. Preparation of the test sample

Add enough commercially available guaranteed reagent grade potassium nitrate to the test sample to make 0.1 mol/L (10 g/L KNO_3).

If the pH of the test sample is not within the correct range, or metal ions (Fe^{3+} , Al^{3+} , etc.) coexist in the test sample, add some buffer solution (Total Ionic strength adjustment buffer) to both the standard solution and the test sample. Use a buffer solution given in JIS K0102-1993, or a commercially available TISAB solution.

When using a buffer solution, there is no need to add the supporting electrolyte (potassium nitrate) described above.

3. Operation

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- (a) Direct potentiometric Method
- (b) Known Addition Method
- (c) Known Subtraction Method
- (d) Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) Direct Potentiometric Method.

With the Direct Potentiometric Method, you first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. You then measure the electrode potential in the sample solution, and arrive at the quantification value using the calibration curve.

7 Measuring

1. Equipment

Assemble the electrode and ion meter (potentiometer) as shown in figure 4.


Fig.4

8 Precautions

1. Do not rub or hit the ion sensitive membrane of the electrode. It may damage the membrane and degrade its special characteristics.
2. Before soaking the electrode in a liquid with a different concentration, be sure to wash it first with deionized water.
3. The electrode connector requires a high degree of insulation. Take care not to touch this part with dirty hands or dampen it.

9 Troubleshooting

1. Check whether the problem is in the ion meter (potentiometer).
2. If there is a drift in the electrode potential, or you are unable to get an electrode reading, take the following steps and try again.
 - Oil may have become attached to the membrane. Clean the resin by lightly wiping the ion sensitive membrane with gauze or tissue paper that has been soaked in ethyl alcohol.
 - Air bubbles may have become attached to the ion sensitive membrane inside the electrode. Grip the cap at the top of the tip firmly with your fingers and shake the electrode as you would shake a thermometer to release the bubbles.

10 Storage

- When storing the electrode for long periods, wash the ion sensitive membrane on the electrode with deionized water, and wipe it with some gauze. Fit the rubber cap included with the package over the response membrane and store the electrode in a cool dark place with a dry atmosphere.

The logo for HORIBA, consisting of the word "HORIBA" in a bold, white, sans-serif font centered on a black rectangular background.

Main Office/Factory: 2-Banchi, Kisshounmiya-
no-higashi-machi Minami-ku, Kyoto 601
Tel: (075) 313-8121
Tokyo Sales Office: 1-7-8 Higashi-kanda
Chiyoda-ku, Tokyo 101
Tel: (03) 3861-8231