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The silver ion-selective electrode may be used to measure the silver ion concentration in any aqueous solution as well as used to measure the silver ion concentration in the pharmaceutical and photographic industries. It may also be used as the indicator electrode for electrical potential difference titration.

In order to ensure correct usage of the ion electrode, be sure to read the operating instructions carefully before usage.

**1 Contents of the Package**

Silver Ion-selective Electrode .....1  
Water proof Abrasive paper (#1200) .....1  
Instruction Manual .....1

**2 Specifications**

Concentration range	0.01mg/L~110,000mg/L Ag <sup>+</sup> (10 <sup>-7</sup> ~1mol/L Ag <sup>+</sup> )	
pH range*1	pH2~10 (1mg/L Ag <sup>+</sup> )	
Temperature range	0~50°C	
Slope(mV/activity decade)**2	59mV	
Response time T90**3	Less than 10 seconds	
Material used for the body	Epoxy resin	
Supported reference electrodes	Double junction referense electrode	
8011-06T	2535A-06T	
8011-10C	2565A-10T	
Reference electrode Type of internal filling solution	Inner tube	3.33 mol/L KCl (Cat. No. 300) (248.3g/L KCl)
	Outer tube	1mol/L KNO3 (100g/L KNO3)

(※1)

This is the permitted pH range for the sample liquid ion concentration value given inside the brackets. (The pH range varies slightly depending on the ion concentration).

(※2)

This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C.)

(※3)

The response time when the concentration is changed from 10<sup>-4</sup> mol/L to 10<sup>-2</sup> mol/L.

### 3 The Electrode


Fig.1

### 4 Electrode Characteristics

#### 1. Detection characteristic

Figure 2 shows the relationship between the chloride ion concentration and the electrode potential.


Fig.2

#### 2. The influence of pH

The silver ion electrode is affected by the pH in strongly acidic or highly alkaline regions. Use the electrode within the 2 to 12 pH range.

Figure 3 shows the relationship between the pH and the difference in output electrode potential for the silver ion electrode.

#### 3. Influence of coexisting ions


Fig.3

The silver ion electrode is influenced by various other ions. The influence of the coexisting ions can be shown in the following equation. K is the selectivity constant and is used as a standard for quantitatively indicating the influence of the coexisting ions.

$$E = E_0 - \frac{2.303RT}{F} \log \left( a_{Ag^+} + \frac{1}{K} a_x \right)$$

Where  $a_{Ag^+}$  is the cyanide ion concentration,  $a_x$  is the concentration of interfering ions. The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the test sample increases.

The selectivity constant K for various ions is listed in Table-1.

Table-1. Selectivity constants

Interfering Ion	Selectivity constant (Limit capacity)
$Cu^{2+}, Cd^{2+}, Pb^{2+}, Zn^{2+}, Mg^{2+}, Ca^{2+}, Na^+, K^+$	1000
$Hg^{2+}$	Not possible

Note : "Not possible" means that these ions cause a chemical change in the ion sensitive membrane.

#### 4. Influence of temperature

The calibration curve of the electrode varies with temperature. You must maintain the standard solution and the sample solution at the same temperature during measurement.

### 5 Preparations Before Measurement

#### 1. Preparation of the standard solution

Accurately weigh out 1.575g of commercially available guaranteed reagent grade silver nitrate and dissolve it in pure (deionized) water. Add enough pure water to make 1L. This solution is your 1000 mg/L  $\text{Ag}^+$  standard solution. Dilute 100mL, 10mL, and 1mL of this solution with pure water to produce 100mg/L, 10mg/L and 1mg/L  $\text{Ag}^+$  standard solutions.

Add 0.1mol/L of potassium nitrate (10g/L  $\text{KNO}_3$ ) as a supporting electrolyte to each of these solutions. In addition, add the same concentration of supporting electrolyte as the standard solution to the test sample to obtain an even higher level of accuracy in your measurements.

#### 2. Preparation of the reference electrode

A reference electrode is necessary for measurement. A double junction electrode with a sleeve type of liquid junction is used. Fill the outer tube with 1mol/L potassium nitrate as the internal solution. An accuracy level of  $100 \pm 2\text{g/L}$  is adequate.

#### 3. Electrode conditioning

Prepare the electrode before taking measurements by soaking it for approximately 10 minutes in the 100 mg/L  $\text{Ag}^+$  standard solution.

### 6 Measuring

#### 1. Equipment

Assemble the electrode and ion meter (potentiometer) as shown in figure 4.


Fig.4

- Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.
- Because of heat generated from the stirrer, the temperature of the solution is slowly rising, causing a drift in the electrode potential. If this is the case, lay a sheet of adiabatic material such as styrofoam between the stirrer and the beaker in order to prevent temperature changes.
- The stirring may generate noise in the electrode potential. If this happens, and you are using an ion meter with a smoothing function, enable the smoothing function.
- Use a recorder if required.

#### 2. Preparation of the test sample

Add enough commercially available guaranteed reagent grade potassium nitrate to the test sample to make 0.1mol/L (10 g/L  $\text{KNO}_3$ ).

#### 3. Operation

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- (a) Direct potentiometric Method
- (b) Known Addition Method
- (c) Known Subtraction Method
- (d) Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) Direct Potentiometric Method.

With the Direct Potentiometric Method, you first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. You then measure the electrode potential in the sample solution, and arrive at the quantification value using the calibration curve.

### 7 Precautions

1. Do not rub or hit the ion sensitive membrane of the electrode. It may damage the membrane and degrade its special characteristics.
2. Before soaking the electrode in a liquid with a different concentration, be sure to wash it first with deionized water.

3. The electrode connector requires a high degree of insulation. Take care not to touch this part with dirty hands or dampen it.

## **8** Troubleshooting

1. Check whether the problem is in the ion meter (potentiometer).
2. If the response from the electrode is slow, or its sensitivity has degraded, take the following steps to revive the ion sensitive membrane.
  - Drop a few drops of water onto the water proof abrasive paper included in the package. Lightly polish the ion sensitive membrane by holding the electrode vertically over the paper and rubbing it in circles. Afterwards, soak the electrode for approximately 10 minutes in 10 mg/L  $\text{Cl}^-$  standard solution before taking measurements.

## **9** Storage

- When storing the electrode for long periods, wash the ion sensitive membrane on the electrode with deionized water, and wipe it with some gauze. Fit the rubber cap included with the package over the response membrane and store the electrode in a cool dark place with a dry atmosphere.

The logo for HORIBA, consisting of the word "HORIBA" in a bold, white, sans-serif font centered within a solid black square.

Main Office/Factory: 2-Banchi, Kisshounmiya-  
no-higashi-machi Minami-ku, Kyoto 601-8510  
Tel: (075) 313-8121  
Tokyo Sales Office: 1-7-8 Higashi-kanda  
Chiyoda-ku, Tokyo 101-0031  
Tel: (03) 3861-8231