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The logo for HORIBA, consisting of the word "HORIBA" in a bold, white, sans-serif font centered on a black rectangular background.

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The potassium ion-selective electrode may be used to measure the potassium ion concentration in water, food, and body tissue test samples, etc. This electrode uses a special plastic membrane, so be sure to read the operating instructions carefully before usage.

1 Contents of the Package

Electrode main unit	1
package - Electrode tips	2
Instruction manual	1

2 Specifications

Concentration range	0.04mg/L ~ 39,000mg/L K ⁺ (10 ⁻⁶ ~ 1mol/L Ag ⁺)	
pH range ^{※1}	pH5 ~ 11 (3.9mg/L K ⁺)	
Temperature range	0 ~ 50°C	
Slope(mV/activity decade) ^{※2}	59mV	
Response time T ₉₀ ^{※3}	Less than 15 seconds	
Material used for the body	Poly acetal resin	
Supported reference electrodes 8202-06T 8202-10C	Double junction reference electrode 2535A-06T 2565A-10T	
Reference electrode Type of internal filling solution	Inner tube	3.33mol/L KCL (248.3g/L KCL)
	Outer tube	3.33mol/L NaCl (195g/L NaCl)

- (※1) This is the permitted pH range for the sample liquid ion concentration value given inside the brackets.
(The pH range varies slightly depending on the ion concentration).
- (※2) This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C.)
- (※3) The response time when the concentration is changed from 10⁻¹ mol/L to 10⁻² mol/L.

POTASSIUM ION-SELECTIVE ELECTRODE TYPE 8202 INSTRUCTION MANUAL CODE : I040175200

3 The Electrode

1. Electrode main unit


Figure 1.

2. Structure of the Tip Type Ion Electrode


Figure 2.

4 Replacing the Electrode

Replacement method

Attachment and removal of the tip electrode from the main unit is via a single screw. Grip the rollette part of the tip electrode with your fingers and rotate it.

5 Electrode Characteristics

1. Calibration characteristics

Figure 3 shows the relationship between the concentration of the potassium ion and the output electrode potential.


Figure 3.

2. Influence of pH

The potassium ion electrode is affected by the pH in strongly acidic or alkaline regions. Use the electrode within the 5 to 11 pH range. In particular, the electrode potential is greatly influenced at pH 11 or higher.

3. Influence of coexisting ions

The potassium ion electrode is influenced by various other ions. The influence of the coexisting ions can be shown in the following equation. K is the selectivity constant and is used as a standard for quantitatively indicating the influence of the coexisting ions.

$$E = E_0 + \frac{2.303RT}{F} \log \left(a_{K^+} + \frac{1}{K} a_x \right) \quad \text{--- ①}$$

Where a_{K^+} is the potassium ion concentration, a_x is the concentration of interfering ions. The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the test sample increases.

The selectivity constant K for various ions is listed in Table-1.

Table-1. Selectivity constants

Interfering Ion	Selectivity constant (Limit capacity)
Li^+	greater than 1000
Na^+	greater than 1000
Rb^+	0.4
Cs^+	3
NH_4^+	70
Mg^{2+}	greater than 1000
Ca^{2+}	greater than 1000
Sr^{2+}	greater than 1000
Ba^{2+}	greater than 1000

(in 10^{-4} mol/L K^+)

Table 1.

4. Influence of Temperature

The calibration curve of the electrode varies with temperature. You must maintain the standard solution and the sample solution at the same temperature during measurement.

6 Preparations

1. Preparations of the standard solution

Accurately weigh out 1.907 g of commercially available guaranteed reagent grade potassium chloride and dissolve it in pure (deionized) water. Add enough pure water to make 1 L. This solution is your 1000 mg/L K⁺ standard solution. Dilute 100 mL, 10 mL, and 1 mL of this solution with pure water to produce 100 mg/L, 10 mg/L and 1 mg/L K⁺ standard solutions.

Add 0.1 M of sodium chloride as a supporting electrolyte to each of these solutions.

2. The reference electrode

A reference electrode is necessary for measurement.

A double junction electrode

with a sleeve type of liquid junction is used. Use 1 mol/L potassium chloride solution in the outer tube.

3. Electrode conditioning

Soak the electrode for 2 to 3 hours in 1 mol/L potassium nitrate solution before taking measurements.

7 Measuring

1. Equipment

Assemble the electrode and ion meter (potentiometer) as shown in figure 4.


Figure 4.

- Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.
- Because of heat generated from the stirrer, the temperature of the solution is slowly rising, causing a drift in the electrode potential. If this is the case, lay a sheet of adiabatic material such as styrofoam between the stirrer and beaker in order to prevent temperature changes.
- The stirring may generate noise in the electrode potential. If this happens, and you are using an ion meter with a smoothing function, enable the smoothing function.
- Use a recorder if required.

2. Operation

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- (a) Direct potentiometric Method
- (b) Known Addition Method
- (c) Known Subtraction Method
- (d) Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) Direct Potentiometric Method.

(a) Direct potentiometric Method

With the Direct Potentiometric Method, you first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. You then measure the electrode potential in the sample solution, and arrive at the quantification value using the calibration curve.

3. Precautions

- Precautionary notes
 - The ion sensitive membrane section is made from a special plastic. Take care not to scratch it or touch it with your fingers. In addition, do not wipe the ion sensitive membrane with an organic solvent.
- Before soaking the electrode in a liquid with a different concentration, be sure to wash it first with deionized water.
- If the test sample is strongly acidic (less than pH 5) or highly alkaline (greater than pH 11), add sodium hydroxide or chloride acid to change the solution to within a 5 to 11 pH range.

4. Troubleshooting

- Check whether the problem is in the potentiometer.
- If there is a drift in the indicator (in the standard solution)
- When there is no response (in the standard solution)

Make a thorough check to ensure that there is no error in the other circuitry before considering the ion sensitive membrane. Next, check whether there are any bubbles on the outside of the ion sensitive membrane. Shake the ion sensitive membrane tip. If bubbles in the internal filling solution are degrading the contact between the membrane and the liquid, the shaking will release the bubbles.

- If there is a lot of noise
Make sure that more than 5 cm of the electrode is immersed in the sample solution. If noise occurs in a concentration of 1000 mg/L, try grounding the liquid to remove induction.

8 After measurement

If you will be reusing the electrode the next day, store it in 10-1 mol/L potassium chloride solution. For long-term storage, wash the electrode well with pure water and wipe it well with gauze. Store it dry in a cool dark place.

9 Other precautionary notes

1. The electrode is composed of organic materials inside plastic and therefore has a limited life span. The life span is over a year when stored in open air.
2. Although sodium chloride was used for the supporting electrolyte in the standard solution, magnesium chloride, barium chloride, or calcium chloride may be used instead.

HORIBA

Main Office/Factory:
2-Banchi, Kisshounmiya-no-higashi-machi
Minami-ku, Kyoto 601
Tel: (075) 313-8121