

# ELECTRODES

рН	ORP	lon	Conductivity
Resistivity	Total Dissolved Solids	Dissolved Oxygen	Salinity







HORIBA

# LAQUA Electrode Technology

# Born from the fusion of our technical expertise and state-of-the-art manufacturing

As a leading pH electrode manufacturer, HORIBA uses the latest technology for all your measurement needs.

Since the development of Japan's first glass electrode for pH meter, HORIBA has focused on continually improving our electrode technology, especially in materials and manufacturing. HORIBA is committed to continually explore and employ groundbreaking solutions in manufacturing next-generation electrodes so that we always provide you with the newest and best electrodes.

pH Ele	octrode	•						3-in-1 ELECT	RODES		
-					PLASTIC			STANDARD ToupH	LONG	MICRO ToupH	SLEEVE ToupH
Select	tion Gl	lide	9651-10D / 9652-10D	9625-10D	9630-10D	9631-10D	9632-10D	9615S-10D	ToupH 9680S-10D	9618S-10D	9681S-10D
	Applicable tempe		0-60/0-80	0-100	0-100	0-60	0-100	0-100	0-100	0-60	0-60
Specification	range (°C) Diameter (mm)		16	16	16	16	16	12	8	3	12
_	Length (mm)		150	150	150	155	150	198	283	185	203
pH - Samp					1			· · · · · · · · · · · · · · · · · · ·			
		rmal (over 100 /m)	۲	۲	۲	۲	۲	۲	۲	۲	۲
		v (approx.10 10 mS/m			۲						۲
	Conductivity Ver	y low (approx.			0						۲
		100 mS/m h (approx.					-				•
Aqueous	5 5	/m)	0	0	0	0	0	0	0		
Solution	Strong alkaline (p Strong acidity (pl						۲	0	0		
	HF sample	I U-2) " Except				۲		۲			
	Quick heat change		۲	۲	۲	۲	۲				
-	High viscosity (ap Containing non-a							-			۲
	solvent	406003						0	0	0	۲
	Suspension							0	0	0	۲
	Inside Surface										
	Gunado							1			
	Microtube/plate (	> 50 μL)								۲	
	-	ø4 mm								۲	
-	Micro container (:	> 2 mL) 3 mm, L:100 ~							0	۲	
Sample Containers	Tube	150 mm							۲		
L		) mL~1L	۲	۲	۲	۲	0	۲	0	0	0
	Large container (: Petri dish	•1L)	0	0	0	0	0	0	۲		
-	Droplet										
	Pure/ion-exchang (approx. 0.1 mS/n water (approx. 0.5	n)/ Distilled						0			۲
Water	Tap/drinking wate 10 mS/m)	r (approx.	0	0	۲			0			۲
-	Surface water Pharmaceutical w	ater/	-		۲			0			۲
	Enviromental wat	er/acid rain	0	0	0			0			0
	Caustic/strong ac HF sample)	id (Except				۲		۲			
Chemical	Hydrofluoric acid					۲					
reagent/ solvent	Surfactant							0			0
-	Water-based pain Dye/coloring age							0			0
	Protein-containin							0		0	
	Medicinal prepara							Ŭ		0	0
Pharmaceutical/ biological	Enzyme solution								0	۲	
sample	Tris buffer							۲		0	
-	Suspension Agar medium							0			۲
	Jam							0			0
	Meat/fish/Fruit/v Dough	egetable/									
Food	Honey										۲
	Cheese/butter Yogurt		0	0	0			0			
	Beer		0	0	0			0			۲
Beverage/	Milk/Carbonated	drink/juice/		~				0			0
seasoning	sauce/soy sauce Mayonnaise/ketc	מער						0			0
	Beauty cream/ma							0			0
Cosmetic/	Gel/soap/shampo							0			0
lotion	lotion Emulsified liquid							0			•
	EIIIUISIIIEU IIIIIII							· · · ·			

			COMBINATION ELECTRODES									
SLEEVE	NEEDLE	PLASTIC	STANDARD ToupH	MICRO ToupH	SLEEVE ToupH	LONG	LONG ToupH	FLAT	GENERAL			
6367-10D	6252-10D	9425-10C	9415-10C	9418-10C	9481-10C	6069-10C	9480-10C	6261-10C	0040-10D			
0-60	0-60	0-100	0-100	0-60	0-60	0-60	0-100	0-50	0-60			
12	12	16	12	3	12	3	8	12	16			
150	150	150	198	185	203	291	283	150	190			

۲	۲	۲	۲	۲	۲	۲	۲	۲	۲
					0				
					0				
		0	0		۲		0		
0			0		0		0		
			۲						
		۲							
0					۲				
0			0	0	0		0		0
			0	0	۲		0		0
	0								
								۲	۲

				۲					
				۲		0			
				۲		0	0		
						۲	۲		
0	0	۲	۲	0	0	0	0	0	0
		0	0				۲		
								۲	۲
								۲	۲

					0			
			0		0	0		
			0		0			
			0		0	0		
			0		۲			
			۲		0			
			۲		0			
			۲					
			۲	0	0			0
			0	0				
		0					0	
			0	0	۲			
			۲		0			
•	•							
(surface)			۲		0		0	
(surface)	0						۲	
O (surface)								
O (surface)	0						0	
(surface)	0		0		0	0	0	0
			۲		0	0		0
			۲		0			0
			۲		0			
			۲		0		0	
			۲		0			
			0		0			

Recommended
 Can be measured
 02 | LAQUA



#### Thick membrane technology

HORIBA's glass moulding technology allows the manufacture of tougher pH glass bulbs.

Double-junction electrodes

All HORIBA pH combination electrodes

are double-junction electrodes. Flexible

to use in a wide-range of applications.

# **Expertise in Manufacturing**

Sophisticated processing technology

HORIBA's in-house expertise in the manufacture of electrodes is the accumulation of more than 60 years of experience. Our sophisticated electrode processing technology provides flexibility in designing various shapes of the electrode bulb and different structural designs of the electrodes.



Unique flat electrode design as well as 3mm diameter micro-electrode with integrated temperature sensor (US Patent No. 7314541/ China Patent No. ZL0315796)

#### **Convenient slider**

Refillable electrodes are equipped with a slider to open or close the refilling port easily.

Fast response & highly accurate

ToupH glass bulb does not compromise responsiveness and sensitivity (US Patent No. 8262877). Specially designed electrodes are available for hydrofluoric acid & strong alkaline application.

Built-in clip for hooking onto electrode stand arm

Top housing of electrodes is designed with a built-in clip to hook onto HORIBA's electrode stands.

#### **ORP Electrodes**

ORP Elect	rodes												
Model			Temp. Rang	ge (°C)		Application							
9300-10D	30140	046710	Pt / Gl	lass	0 - 60	)	Waterproo	f & refillable; Flat plat	inum tip; Sเ	uitable for lab us	e.		
9301-10D	3200	922105	Pt / Pla	astic	0 - 80	)	Waterproo	f & maintenance-free	; Flat platin	tinum tip; Suitable for field use.			
on Selecti	ve Electro	des (IS	Es)										
Model	Part No.	Con	nbination I	ISE	Temp. Range (°C)		Measurement Range		Replacement Tip		Part No.		
5002S-10C	3200698386	Ammo	onia (NH <sub>3</sub> )	)	0 - 50		0.01 - 1	3,000 mg/L NH <sub>4</sub> +	NH <sub>3</sub> Me	mbrane Caps	3200705774		
6583S-10C	3200697410	Calciu	um (Ca <sup>2+</sup> )		0 - 50		0.4 - 40	),080 mg/L Ca <sup>2+</sup>	-	7683S	3200697414		
6560S-10C	3200697407	Chlori	ide (Cl-)		0 - 50		0.35 - 3	35,000 mg/L Cl	-	7660S	3200697411		
6561S-10C	3200693774	Fluori	de (F <sup>-</sup> )		0 - 50		0.02 -	19,000 mg/L F	-	7661S	3200693606		
6581S-10C	3200697408	Nitrate	e (NO <sub>3</sub> -)		0 - 50		0.62 - 6	2,000 mg/L NO <sub>3</sub> -		7681S	3200697412		
6582S-10C	3200697409		sium (K+)		0 - 50			39,000 mg/L K+	-	7682S	3200697413		
Conductivi	itv Cells								r				
Туре	Model	Part I	No.	Temp	. Range (°C)	Cell	Constant	Measurement F	Range	Ap	plication		
	3551-10D	301408	01710	0 - 60 0.1 cm <sup>-1</sup> 0.1 µS/cm - 10 mS/cm Low conductiv		ty water (e.g.,							
	3551-100	301400	51712		0 - 00	10 m <sup>-1</sup>		10 µS/m - 1 S/m		deionized, dist	illed)		
	9371-10D	320087	78882		0 - 100			0.01 μS/cm - 500 μS	S/cm	Low conductiv	ity water and ultra		
	3371-100	520007	0002		0 - 100	10 m <sup>-1</sup>		1 µS/m - 50 mS/m		pure water			
	9382-10D	301404	6700		0 - 80			1 µS/cm - 100 mS/c	m	Gonoral purpo	se use; Waterproc		
Submersible	9302-100	301404	0709		0 - 80	100 m <sup>-1</sup>		0.1 mS/m - 10 S/m			se use, waterproc		
Submersible	9383-10D	320078	0927		0 - 80	1 cm <sup>-1</sup>	1 µS/cm - 100 mS		S/cm - 100 mS/cm		General purpose use; Waterprod		
		020070	10021		0 00	100 m <sup>-1</sup>		0.1 mS/m - 10 S/m		General pulpose use, waterpro			
	3552-10D	301408	1545		0 - 100	1 cm <sup>-1</sup>		1 µS/cm - 100 mS/cm		General purpo			
	0002 100	001400	1040		0 100	100 m <sup>-1</sup>		0.1 mS/m - 10 S/m			30 030		
	3553-10D	301408	31714		0 - 60	10 cm <sup>-1</sup>		10 µS/cm - 1 S/cm		High conductiv	vity water		
	0000 102					1000 m		1 mS/m - 100 S/m		- iigii conadoai	ity mator		
	3561-10D	301408	2350		0 - 60	0.1 cm <sup>-1</sup>		0.1 µS/cm - 10 mS/c					
						10 m <sup>-1</sup>		10 µS/m - 1 S/m		deionized, dist	llied)		
	3562-10D	301408	32513		0 - 60	1 cm <sup>-1</sup>		1 µS/cm - 100 mS/c	m	General purpo	se use		
Flow						100 m <sup>-1</sup>		0.1 mS/m - 10 S/m					
	3573-10C	301408	2590		0 - 60	10 cm <sup>-1</sup> 1000 m	.1	10 µS/cm - 1 S/cm 1 mS/m - 100 S/m		High conductiv	vity water		
						1000 m	•	1 mS/m - 100 S/m 10 µS/cm - 100 mS/		0 11 1			
	3574-10C	301408	32592		0 - 60	10 cm <sup>-1</sup>	-1	1 mS/m - 10 S/m	GHI	Small volume s	ample (e.g., colur		
						1000 11		1 113/11 - 10 3/11			17		

• Material: All have platinum-platinum black / glass-body, except 9382-10D and 9383-10D (titanium-platinum black / plastic-body) and 9371-10D (stainless steel).

Dissolved	Oxygen Probes					
Туре	Model	Part No.	Temp. Range (°C)	Measurement Range	Replacement Tip	Part No.
Field	9551-20D / 9551-100D	3014047090 / 3014047091	0 - 40	0 - 19.99 mg/L DO	5401	3014072770
Field	9552-20D / 9552-50D	3200780939 / 3200780941	0 - 50	0 - 20.00 mg/L DO	5402	3200781553
Lab	9521-10D	3200891722	0 - 50	0 - 20.00 mg/L DO	7544	3200891724

# pH Combination Electrodes

HORIBA pH Combination electrodes manufactured with 1 meter cable terminating in BNC connector allow these electrodes to be used with any pH meter<sup>1</sup>. Enjoy the full spectrum of features and benefits of these electrodes on your existing pH meter<sup>1</sup>. (For applications where temperature measurement and compensation is required, please refer to the 3-in-1 pH electrodes).

<sup>1</sup> pH meters must have BNC connector

Model	pH Range	Operating Temperature Range (°C)	Liquid Junction	Application
ToupH Standard Electrode         9415-10C       General laboratory application         Overall length: 151 mm         S200611623	0-14	0-100	Ceramic	<ul> <li>The electrode offers quick stability and drift reduction.</li> <li>Constructed with responsive glass that is 10X stronger than JIS standards</li> <li>The one-touch refilling port slider allows one-hand operation</li> <li>Waterproof, Pb-free glass</li> <li>Perfect for preparing pH buffers and other aqueous test solutions.</li> </ul>
Standard Plastic Electrode 9425-10C General field application	0-14	0-100	Ceramic	<ul> <li>The electrode has plastic body, which is ideal for field measurement.</li> <li>Can be submerged up to 1m depth and 30mins (with refilling port closed)</li> <li>Waterproof, Pb-free glass</li> <li>Recommended for field use. For measurement of tap water and drinking water.</li> </ul>
ToupH Sleeve Electrode         9481-10C       High viscosity application         Overall length: 151 mm         Diameter of probe: 12 mm         Connector: BNC	0-14	0-60	Movable sleeve	<ul> <li>The electrode gives stable readings in highly viscous samples.</li> <li>The liquid junction is designed with a movable sleeve that can be cleaned easily and prevents clogging</li> <li>Waterproof, Pb-free glass</li> <li>For measurement of highly viscous samples and samples containing non-aqueous solvents (e.g., cosmetics, paints).</li> </ul>
ToupH Micro Electrode         9418-10C Precious trace amount sample         Overall length: 151.5 mm Diameter of probe: 3 mm Connector: BNC	0-14	0-60	Ceramic	<ul> <li>The electrode can measure samples as small as 50µL.</li> <li>Compatible with extremely small containers (e.g., micro tubes)</li> <li>Temperature sensor is placed next to the bulb for quick response</li> <li>Waterproof</li> <li>Suitable for low-volume samples and wide range of aqueous solutions.</li> </ul>
ToupH Long Electrode         9480-10C For large containers and long test tubes         Overall length: 251 mm         3200611628	0-14	0-100	Ceramic	<ul> <li>The long, thin body of the electrode is perfect for large containers and test tubes.</li> <li>283mm length, 8mm diameter</li> <li>Constructed with responsive glass that is 10X stronger than JIS standards</li> <li>Waterproof, Pb-free glass</li> <li>For measuring samples (e.g., microbial culture fluids) in test tubes and tall beakers.</li> </ul>
Long Electrode 6069-10C For very slender test tubes	0-14	0-60	Ceramic	<ul> <li>The long, thin body of the electrode is perfect for very slender test tubes.</li> <li>291mm length, 3mm diameter</li> <li>Waterproof</li> <li>For measuring samples in slender tubes (e.g., NMR test tube).</li> </ul>

	Model	pH Range	Operating Temperature Range (°C)	Liquid Junction	Application
Flat Electrode 6261-10C	Overall length: 150 mm Diameter of probe: 12 mm Connector: BNC	0-12	0-50	Sleeve	<ul> <li>The sensor is located on the flat surface of the tip.</li> <li>Measurement can be made from minute amount of moisture on solid sample surface</li> <li>Pure water can be applied for samples with no moisture</li> <li>Waterproof</li> <li>Perfect for measuring samples in shallow containers (e.g., petri dishes) and gelatinous materials (e.g., nutrient agar). For surface measurement of meat, paper, skin, and cloth.</li> </ul>

# 3-in-1 pH Glass Body Electrodes<sup>2</sup>

HORIBA pH Combination electrodes with an integrated thermistor offer higher accuracy as these electrodes measure temperature concurrently with pH. The pH meter is able to continuously monitor and compensate for temperature effects automatically.

<sup>2</sup>Only compatible with HORIBA pH meters

Operating Application Range (°C) The electrode offers quick stability and **ToupH Standard Electrode** drift reduction. 9615S-10D General laboratory application Constructed with responsive glass that is 10x stronger than JIS standards The one-touch refilling port slider allows 0-14 0-100 Ceramic one-hand operation Waterproof, Pb-free glass Overall length: 151 mm Diameter of probe: 12 mm Connectors: BNC & phono jack Perfect for preparing pH buffers and other 3200585428 aqueous test solutions. **ToupH Sleeve Electrode** The electrode gives stable readings in highly viscous samples. 9681S-10D High viscosity application The liquid junction is designed with a movable sleeve that can be cleaned easily Movable and prevents clogging 0-14 0-60 sleeve • Waterproof, Pb-free glass For measurement of highly viscous samples and Overall length: 151 mm 3200585463 samples containing non-aqueous solvents (e.g. Diameter of probe: 12 mm nectors: BNC & phono jack Conn cosmetics, paints). The electrode can measure samples as small as **ToupH Micro Electrode** 50µL. 9618S-10D Precious trace amount sample Compatible with extremely small containers . (e.g. micro tubes) Temperature sensor is placed next to the 0-14 0-60 Ceramic bulb for quick response Waterproof Overall length: 151.5 mm Suitable for low-volume samples and a wide 3200585447 Diameter of probe: 3 mm Connectors: BNC & phono jack range of aqueous solutions. The long, thin body of the electrode is perfect for **ToupH Long Electrode** large containers and test tubes. 9680S-10D For large containers and long test tubes 283mm length, 8mm diameter . Constructed with responsive glass that is • 10x stronger than JIS standards 0-14 0-100 Ceramic Waterproof, Pb-free glass • Overall length: 251 mm Diameter of probe: 8 mm Connectors: BNC & phono jack For measuring samples (e.g. microbial culture 3200585428 fluids) in test tubes and tall beakers. **Needle Electrode** 6252-10D For food application Needle electrode allows measurement of food LAQUA 0-12 0-60 Ceramic samples and aqueous solutions. Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack 3014080850

	Model	pH Range	Operating Temperature Range (°C)	Liquid Junction	Application
Standard Sleeve E 6367-10D	Electrode				
		0-14	0-60	Sleeve	Uses a sleeve at the liquid junction for improved stability and repeatability. For measuring pH at high accuracy.
3014079136	Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack				

# 3-in-1 pH Plastic Body Electrodes<sup>2</sup>

<sup>2</sup> Only compatible with HORIBA pH meters

				<sup>2</sup> Only compatible with HORIBA pH meters		
Model	pH Range	Operating Temperature Range (°C)	Liquid Junction	Applications		
Gel-filled pH Electrode 9651-10D For Field 3200642020 Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack	0-14	0-80	Porous sintered polyethylene	<ul> <li>The plastic body of the electrode is filled with gel electrolyte. Less maintenance is needed as refilling is not required.</li> <li>Can be submerged up to 1m depth of water for 30mins.</li> <li>Waterproof, Pb-free glass</li> <li>Recommended for field use.</li> </ul>		
Gel-filled pH Electrode 9652-10D; 9652-20D For Field	0-14	0-80	Porous sintered polyethylene	<ul> <li>The plastic body of the electrode is filled with gel electrolyte. Less maintenance is needed as refilling is not required.</li> <li>Can be submerged up to 1m depth of water for 30mins.</li> <li>Waterproof, Pb-free glass</li> <li>Recommended for field use.</li> </ul>		
Standard Plastic Electrode 9625-10D; 9625-20D; 9625-30D For Field Constant of protection of the standard of th	0-14	0-100	Ceramic	<ul> <li>The electrode has a plastic body which is ideal for field measurement.</li> <li>Can be submerged up to 1m depth of water for 30mins. (with refilling port closed)</li> <li>Waterproof, Pb-free glass</li> <li>Recommended for field use. For measurement of tap water and drinking water.</li> </ul>		
Hydrofluoric Acid Resistant Electrode 9631-10D	2-12	0-60	Ceramic	<ul> <li>The electrode can measure 1% hydrofluoric acid solution (at 25°C, immersed at 1min.) for about 1000 times.</li> <li>Rolled glass design for long-term reliable measurement and easy maintenance</li> <li>Compliant with Japan's Measurement Act Certification</li> <li>Waterproof, Pb-free glass</li> <li>Suitable for drain water measurement after etching process.</li> </ul>		
Strong Alkali Resistant Electrode 9632-10D	0-14	0-100	Ceramic	The alkali-resistant glass membrane has higher resistance and longer stability (about 5X in 0.1mol/L sodium at 60°C, pH 13) than conventional electrodes. • Waterproof, Pb-free glass Suitable for strong alkali samples such as plating solutions.		
Standard Plastic Electrode 9630-10D For tap water Overall length: 150 mm Diameter of probe: 16 mm Connectors: BNC & phono jack	0-14	0-100	Ceramic	<ul> <li>The electrode can measure samples with low conductivity or buffering capacity.</li> <li>Made of high purity multicomponent lithium series glass</li> <li>Waterproof, Pb-free glass</li> <li>Suitable for tap water measurement and quality control in water purification plant. Recommended to use with cleaning solution 230.</li> </ul>		

# 3-in-1 ORP Electrodes

HORIBA ORP electrodes have flat platinum tip that allows measurement of low-volume samples and thermistor that ensures accurate temperature reading during ORP measurement.

Model	Operating Temperature Range (°C)	Electrode Material	Liquid Junction	Applications
Standard ORP Electrode 9300-10D For Laboratory Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack	0-60	Pt / Glass	Ceramic	Waterproof & refillable with 3.33M KCI; Suitable for lab use.
Gel-filled ORP Electrode 9301-10D For Field LAOUA Diameter of probe: 16 mm 2200922105	0-80	Pt / Plastic	Porous sintered polyethylene	Waterproof & maintenance-free; Suitable for field use.

# **Dissolved Oxygen (DO) Electrodes & Tips**

HORIBA Dissolved Oxygen (DO) electrodes are galvanic probes with integrated temperature sensors. With galvanic DO probes, calibration can be performed immediately and in air. The HORIBA DO probes use unique and innovative tips which are replaceable. No need to replace membranes or refill electrolytes.

Three models are available: 9520 that can be used for lab BOD measurements, and 9551 / 9552 housed in a rugged casing available in 2m, 5m, and 10m cable configurations for field use. The 9520 model is fitted with a rotor as well as an adaptor to facilitate BOD measurements.

#### **DO Electrodes**

Ма	odel	Measurement Range	Response Time	Temperature Range (°C)	Features
<b>9521-10D</b> For laborato	Verall length: 183 mm Diameter of probe: 15.6 mm Connectors: BNC & phono jack	0-20.00mg/L DO	20 seconds (90% response time at constant temperature)	0-50	Waterproof; Comes with replaceable DO tip 7544.
<b>9551-20D; 9551-100D</b> 3014047090; 3014047091	For field Overall length: 165 mm Diameter of probe: 32 mm Connectors: BNC & phono jack	0-19.99mg/L DO	30 seconds (90% response time at constant temperature)	0-40	Waterproof; Comes with replaceable DO tip 5401.
9552-20D; 9552-50D 3200780939; 3200780941	For field Overall length: 165 mm Diameter of probe: 30 mm Connectors: BNC & phono jack	0-20.00 mg/L DO	30 seconds (90% response time at constant temperature)	0-50	Waterproof; Comes with replaceable DO tip 5402.

#### **DO Electrode Tips**

Model	Description
7541 3014074145 Overall length: 26.5 mm Diameter: 15 mm	Replacement DO electrode tip for 9520-10D
7544 3200891724 Overall length: 43 mm Diameter: 15 mm	Replacement DO electrode tip for 9521-10D
<b>5401</b>	Replacement DO electrode tip for 9551-20D and 9551-100D
5402 3200781553 Overall length: 19 mm Diameter: 21 mm	Replacement DO electrode tip for 9552-20D and 9552-50D

# Conductivity Electrode Cells

HORIBA Conductivity cells are available as Submersible type and Flow type, as well as in a variety of cell constants ranging from 0.1 to 10.0.

The HORIBA Conductivity cells are integrated with temperature sensor (except for 3573 & 3574) and the wetted material is either **Stainless Steel or Platinum / Titanium coated with Platinum black**. Rugged Titanium allows cell to be used in a wide range of applications, including highly corrosive samples such as concentrated acids and sea water. Maintenance is simple – soak in deionized/demineralized water or with the conditioning solution. The 9371-10D stainless steel conductivity cell is compatible with the glass flow cell.

#### Conductivity Cells (Submersible Type)

Model	Cell Constant	Measurement Range	Temp. Range (°C)	Cell Material	Thermistor	Minimum Sample Volume (ml)	Application				
3551-10D	0.1 cm <sup>-1</sup>	0.1 µS/cm - 10 mS/cm	0 - 60	Pt-Pt black /	Built-in	50	Low conductivity water (e.g.,				
Overall ler Diameter of p 3014081712 Connectors: BNC	ngth: 175 mm probe: 23 mm & phono jack	10 µS/m - 1 S/m	0 - 00	Glass	Duitein	50	deionized, distilled)				
3552-10D	1 cm-1	1 µS/cm - 100 mS/cm	0 - 100	Pt-Pt black /	Built-in	15	General				
Overall ler Diameter of p 3014081545 Connectors: BNC	ngth: 150 mm probe: 12 mm & phono jack	0.1 mS/m - 10 S/m	0 - 100	0 - 100 black / Glass					Duit-in	15	purpose use
3553-10D	10 cm <sup>-1</sup>	10 µS/cm - 1 S/cm	0 - 60	Pt-Pt black /	Built-in	50	High conductivity				
Overall ler Width of r 3014081714 Connectors: BNC	ngth: 175 mm orobe: 28 mm & phono jack	1 mS/m - 100 S/m	0-00	Glass	Duit-in	50	water				
9371-10D	0.1 cm <sup>-1</sup>	0.01 µS/cm - 500 µS/cm	0 - 100	Stainless	Built-in	20-30	Low conductivity				
Overall len Diameter of p 3200878882 Connectors: BNC &	gth: 180 mm robe: 16 mm & phono jack	1 µS/m - 50 mS/m	0 - 100 Steel		Steel	Steel	20-30	water and ultra- pure water			
9382-10D	1 cm-1	1 µS/cm - 100 mS/cm	0.00	Ti-Pt black		00.00	General				
Overall len Diameter of p 3014046709 Connectors: BNC 8	gth: 150 mm robe: 16 mm & phono jack	0.1 mS/m - 10 S/m	0 - 80	0 - 80 / Plastic	Built-in	20-30	purpose use; Waterproof				
9383-10D	1 cm-1	1 µS/cm - 100 mS/cm	0.90	Ti-Pt black	Duiltin	20.20	General				
Overall len Diameter of p 3200780927 Connectors: BNC 8	gth: 150 mm robe: 16 mm & phono jack	0.1 mS/m - 10 S/m	0 - 80 / Plastic	n <u>- 80   tateren  </u>	0 - 80		Built_in	20-30	purpose use; Waterproof		

#### **Conductivity Cells (Flow Type)**

	Model	Cell Constant	Measurement Range	Temp. Range (°C)	Cell Material	Thermistor	Minimum Sample Volume (ml)	Application
3561-10D		0.1 cm <sup>-1</sup>	0.1 µS/cm - 10 mS/cm	0 - 60	Pt-Pt black /	Built-in	10	Low conductivity water (e.g.,
3014082350	Overall length: 143 mm Diameter of probe: 18 mm Connectors: BNC & phono jack	10 m <sup>-1</sup>	10 µS/m - 1 S/m		Glass		10	deionized, distilled)
3562-10D		1 cm <sup>-1</sup>	1 µS/cm - 100 mS/cm	0 - 60	Pt-Pt black /	Built-in	16	General
3014082350	Overall length: 205 mm Diameter of probe: 18 mm Connectors: BNC & phono jack	100 m <sup>-1</sup>	0.1 mS/m - 10 S/m	0-00	Glass	Duit in	10	purpose use
3573-10C	uan per	10 cm <sup>-1</sup>	10 µS/cm - 1 S/cm	0 - 60	Pt-Pt black /		4	High conductivity
3014082590	Overall length: 222 mm Diameter of probe: 18 mm Connector: BNC	1000 m <sup>-1</sup>	1 mS/m - 100 S/m	0-00	Glass		4	water
3574-10C		10 cm <sup>-1</sup>	10 µS/cm - 100 mS/ cm	0 - 60	Pt-Pt black /		0.25	Small volume sample (e.g.,
3014082592	Overall length: 136 mm Diameter of probe: 66 mm Connector: BNC	1000 m <sup>-1</sup>	1 mS/m - 10 S/m	0 - 60	Glass		0.25	column chromatography)

# **Combination ISE**

lon-selective electrodes are responsive to concentration of particular ions in the test liquid and are variable-potential electrodes. They are used in conjunction with reference electrodes to measure the concentration of particular ions. HORIBA's years of experience and know-how in this field are behind the wide range of ion electrodes we offer.

When measurements are made using an ion meter, calibrating it with various standard solutions will give direct readings of the ion concentration. Note that since volume-detection level changes with temperature, measurements must be taken at a fixed temperature.

Model	Accessories Included	Temp. Range (°C)	Measurement Range	pH Range
Ammonia (NH <sub>3</sub> ) electrode 5002S-10C 3200698386 Overall length: 150 mm Diameter of probe: 15 mm Connector: BNC	<ul> <li>membrane cap, 3pcs</li> <li>1000mg/L ammonium ion standard solution, 50ml</li> <li>100mg/L ammonium ion standard solution, 50ml</li> <li>ammonia electrode filling solution, 50ml</li> <li>syringe</li> <li>dropper</li> <li>protective pipe</li> <li>manual</li> </ul>	0 - 50	0.01 - 18,000 mg/L NH <sub>4</sub> * (5 x 10 <sup>-7</sup> to 1 mol/L NH <sub>4</sub> *)	pH 12 or more
Calcium ion (Ca <sup>2+</sup> ) electrode 6583S-10C 3200697410 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC	<ul> <li>calcium electrode tip, 2pcs</li> <li>1000mg/L calcium ion standard solution, 50ml</li> <li>100mg/L calcium ion standard solution, 50ml</li> <li>calcium electrode filling solution, 50ml</li> <li>calcium ionic strength adjustor, 50ml</li> <li>syringe</li> <li>dropper</li> <li>protective pipe</li> <li>manual</li> </ul>	0 - 50	0.4 - 40,080 mg/L Ca <sup>2+</sup> (10 <sup>-5</sup> to 1 mol/L Ca <sup>2+</sup> )	4.0 mg/L (10 <sup>-4</sup> mol/L) Ca <sup>2+</sup> , pH 5 to 11
Chloride ion (CI <sup>-</sup> ) electrode 6560S-10C 3200697407 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC	<ul> <li>chloride electrode tip</li> <li>1000mg/L chloride ion standard solution, 50ml</li> <li>100mg/L chloride ion standard solution, 50ml</li> <li>syringe</li> <li>dropper</li> <li>protective pipe</li> <li>water-resistant abrasive sheet</li> <li>manual</li> </ul>	0 - 50	0.35 - 35,000 mg/L Cl <sup>.</sup> (10 <sup>.5</sup> to 1 mol/L Cl <sup>.</sup> )	350 mg/L (10 <sup>-2</sup> mol/L) CI <sup>-</sup> , pH 3 to 11
Fluoride ion (F <sup>-</sup> ) electrode 6561S-10C 3200693774 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC	<ul> <li>fluoride electrode tip</li> <li>1000mg/L fluoride ion standard solution, 50ml</li> <li>100mg/L fluoride ion standard solution, 50ml</li> <li>fluoride electrode filling solution, 50ml</li> <li>fluoride ionic strength adjustor, 50ml</li> <li>dropper</li> <li>protective pipe</li> <li>manual</li> </ul>	0 - 50	0.02 - 19,000 mg/L F <sup>-</sup> (10 <sup>-6</sup> to 1 mol/L F <sup>-</sup> )	0.1 to 1,000 mg/L F <sup>-</sup> , pH 5 to 8
Nitrate ion (NO <sub>3</sub> <sup>-</sup> ) electrode 6581S-10C 3200697408 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC	<ul> <li>nitrate electrode tip, 2pcs</li> <li>1000mg/L nitrate ion standard solution, 50ml</li> <li>100mg/L nitrate ion standard solution, 50ml</li> <li>nitrate electrode filling solution, 50ml</li> <li>nitrate ionic strength adjustor, 50ml</li> <li>syringe</li> <li>dropper</li> <li>protective pipe</li> <li>manual</li> </ul>	0 - 50	0.62 - 62,000 mg/L NO <sub>3</sub> - (10 <sup>-5</sup> to 1 mol/L NO <sub>3</sub> -)	62 mg/L (10 <sup>-3</sup> mol/L) NO <sub>3</sub> -, pH 3 to 7
Potassium ion (K*) electrode 6582S-10C 3200697409 Overall length: 150 mm Diameter of probe: 16 mm Connector: BNC	<ul> <li>potassium electrode tip, 2pcs</li> <li>1000mg/L potassium ion standard solution, 50ml</li> <li>100mg/L potassium ion standard solution, 50ml</li> <li>potassium electrode filling solution, 50ml</li> <li>potassium ionic strength adjustor, 50ml</li> <li>syringe</li> <li>dropper</li> <li>protective pipe</li> <li>manual</li> </ul>	0 - 50	0.39 - 39,000 mg/L K* (10 <sup>-5</sup> to 1 mol/L K*)	3.9 mg/L (10 <sup>.4</sup> mol/L) K <sup>+</sup> , pH 5 to 11



	Selection Coefficient	Replacement Tip	Electrode Filling Solution	100mg/L Standard Solution	1000mg/L Standard Solution	Ionic Strength Adjustor	Applications
		NH <sub>3</sub> electrode membrane caps 3200705774	500-NH3-IFS 3200697173	500-NH4-SL 3200697172	500-NH4-SH 3200697171	500-NH3-ISA 3200697174	Agriculture, Soil, Power Station Water, Fish Tanks, Sea Water, Waste Water, Plating Baths, Air / Stack Gases and Biological Cultures or Samples
1 N	e <sup>3+</sup> = 0.1, Fe <sup>2+</sup> , Zn <sup>2+</sup> = 1, Sr <sup>2+</sup> = 50 Ji <sup>2+</sup> , Cu <sup>2+</sup> = 70, Co <sup>2+</sup> = 350 Mn <sup>2+</sup> = 500, Mg <sup>2+</sup> = 1,000 Ja <sup>+</sup> , K <sup>+</sup> , Ba <sup>2+</sup> , NH <sub>4</sub> <sup>+</sup> = over 1,000	<b>7683S</b> 3200697414	500-CA-IFS 3200697177	500-CA-SL 3200697176	500-CA-SH 3200697175	500-CA-ISA 3200697178	Agriculture / Plant Tissue, Soil, Water Softening Systems, Boiler Feed Water, Drinking / Mineral Water, Biological Cultures, Dental / Clinical Analysis and Dairy / Food / Beverages Applications
e S E	$S_{2}O_{3}^{2-}$ , S <sup>2-</sup> , I <sup>-</sup> , Ag <sup>+</sup> , Hg <sup>2+</sup> = Not icceptable SCN <sup>-</sup> = 0.3, MnO <sub>4</sub> <sup>-</sup> = 0.1 Br = 0.03 $NO_{3}^{-}$ , F <sup>-</sup> , HCO <sub>3</sub> <sup>-</sup> , SO <sub>4</sub> <sup>2-</sup> , PO <sub>4</sub> <sup>2-</sup> = 1,000	<b>7660S</b> 3200697411	500-CL-IFS 3200697169	500-CL-SL 3200697168	500-CL-SH 3200697167	500-CL-ISA 3200697170	Agriculture, River / Tap Water, Plant Tissue, Soils, Boiler Feed Water, Clinical Analysis, Sweat, Urine, Cement, Plating Baths and Dairy / Food / Beverages Samples
r	Possible interference when nultiply-charged ion (ex. Al <sup>3+</sup> , Fe <sup>3+</sup> ) poexisted and foamed the complex.	7661S 3200693606	500-F-IFS 3200697165	500-F-SL 3200697164	500-F-SH 3200697163	500-F-TISAB 3200697166	Dental / Toothpaste / Mouth Wash, Drinking / Seawater, Wastewater, Air / Stack Gases, Acids, Soils, Food, Biological Fluids, Plant Tissue, Coal, Carbonated Beverages and Bone
Λ	$CIO_{4^{-}}$ , I <sup>-</sup> = Not acceptable, Br= 2 $IO_{2^{-}} = 3$ , CI <sup>-</sup> = 300 $ICO_{3^{-}}$ , H <sub>2</sub> PO <sub>4</sub> <sup>-</sup> , SO <sub>4</sub> <sup>2-</sup> =over 1000	<b>7681S</b> 3200697412	500-NO3-IFS 3200697181	500-NO3-SL 3200697180	500-NO3- SH 3200697179	500-NO3-ISA 3200697182	Agriculture / Plant Tissue / Fertilizers, Surface / Seawater / Drinking Water, Sewage Effluent, Soils, Meats, Vegetables, Foods / Beverages
L	Rb* = 0.4, Cs* = 3, NH <sub>4</sub> * = 70 i*, Na*, Mg <sup>2+</sup> , Ca <sup>2+</sup> , Sr <sup>2+</sup> , Ba <sup>2+</sup> = over ,000	7682S 3200697413	500-K-IFS 3200697185	500-K-SL 3200697184	500-K-SH 3200697183	500-K-ISA 3200697186	Agriculture / Plant Tissue, Soils, Wastewater, River / Tap Water, Clinical Analysis, Saliva, Serum, Fertilizers, Soils and Wines, Dairy / Foods / Beverages

Note: Detailed information on standard solutions, ISAs, and filling solutions can be found on page 16

# LAQUA WQ-300 Series Smart Digital Sensors



# HORIBA's Smart digital sensor technology

All sensor heads are offered in 2-m and 5-m cable versions. A 10-m extension cable is also available.

Minimum depth for

temperature sensor

Detachable glass flow cell

2.2cm

#### Maintenance-free, gel-filled pH sensor

No electrolyte refilling required

- KCl gel electrolyte
- Double junction reference
- Porous sintered polyethylene junction
- Built-in temperature sensor
- Rugged polycarbonate body
- Replaceable pH sensor cartridge

#### 2-Cell & 4-Cell conductivity sensors


Wide range of conductivity measurements possible

- 2-Cell conductivity sensor with flow cell is designed for ultra-pure water applications
- From clean water to industrial wastewater, the 4-cell type can measure a variety of samples with different conductivities
- Built-in temperature sensor
- Stainless steel 2-cell cartridge
- Durable epoxy / carbon body 4-cell cartridge
- Replaceable conductivity sensor cartridge

#### Optical dissolved oxygen (DO) sensor

Longer usable life with excellent performance

- Easy to handle not affected with sample flow velocity, not sensitive to hydrogen sulfide, DO sensor cap replacement after 1-2 years\*
- Built-in temperature sensor
- Comes with replaceable DO sensor cap, air calibration bottle and Stainless Steel DO Sensor Protective Guard



#### Maintenance-free, gel-filled ORP sensor

#### No electrolyte refilling required

- Platinum tip attached to glass
- KCl gel electrolyte
- Porous sintered polyethylene junction
- Built-in temperature sensor
- Rugged polycarbonate body
- Replaceable ORP sensor cartridge

#### Ion sensor head

Compatible with conventional ion selective electrodes

- Accepts all combination ion selective electrodes with BNC connector
- Requires sensor head adapter

\*Depending on application, handling, and maintenance

# LAQUA WQ-300 Series Sensor Specifications

pH Sensor Head	d pH / mV / Temp (°C/°F)			
Model	300PH-2	300PH-5		
Part No.	3200812206	3200812207		
pH Range		20.00 рН 20.000 рН		
Resolution	-2.00 to 20. -2.000 to 20.0	00: 0.01 pH 200: 0.001 pH		
Accuracy	-2.00 to 20.00: ±0.01 -2.000 to 20.000: ±0.005			
Calibration Points	Up to 6			
pH Buffer Groups	USA, DIN, NIST, NIST10, Custom			
mV Range	±1000.0 mV			
Resolution	0.1	mV		
Accuracy	±0.1	mV		
Temperature Range		130.0 °C 266.0 °F		
Resolution	0.1 °(	C/°F		
Accuracy	±0.5 °C /	/ ±0.9 °F		
Calibration Option	Y	es		
Body Material	ABS / Polycarbonate			
Length and Diameter	85 x 3	80 mm		
Connector	Push	i-pull		
Cable Length	2 m	5 m		

	ridge pH/mV/Temp (°C/°F)
Model	300-P-C
Part No.	3200786363
pH Range	-2.00 to 20.00 pH -2.000 to 20.000 pH
Temperature Range	0 to 80 °C -32.0 to 176.0 °F
Junction Material	Porous sintered polyethylene
Double Junction	Yes
Temperature Sensor	Built-in
Length and Diameter	110 x 16 mm
Body Material	Polycarbonate, glass bulb

<b>Dissolved Oxygen Senso</b>	<b>r</b> DO (ma/l %) / (	D <sub>2</sub> / Temp (°C/°F)		
Model	300-D-2	300-D-5		
Part No.	3200780940	3200780942		
Dissolved Oxygen (DO) Range		0.00 mg/L 200.0 %		
Resolution	0.01 mg	ı/L, 0.1%		
Accuracy	±0.2 mg	/L, ±2 %		
Salinity Compensation	Auto: by conductivity sensor / Manual: 0.0 to 4 ppt			
Barometric Pressure Compensation		er / Manual: 10.0 to 199.9 Pa		
Calibration Points	Up	to 2		
Oxygen Range	0.0 to 50.0%			
Resolution	0.	1%		
Accuracy	±0	.5%		
Temperature Range		130.0 °C 266.0 °F		
Resolution	0.1 °	C/°F		
Accuracy	±0.5 °C	/ ±0.9 °F		
Calibration Option	Y	es		
Body Material	ABS / Polycarbonate			
Length and Diameter	200 x 16 mm			
Connector	Push-pull			
Cable Length	2 m	5 m		
Sensor cap included	1			

<b>Dissolved Oxygen Senso</b>	r Cap
Model	300-D-M
Part No.	3200781554
DO Range	0.00 to 20.00 mg/L 0.0 to 200.0 %
Temperature Range	0 to 50.0 °C 32.0 to 122.0 °F
Length and Diameter	10 x 16 mm
Body Material	PVC, PMMA

<b>Conductivity Sense</b>	or Head EC/Sal/TC	S / Res / Temp (°C/°F)		
Model	300-C-2	300-C-5		
Part No.	3200784468	3200812202		
Conductivity Range	µS/cm 0.000 to 0.199 0.200 to 1.999 2.00 to 19.99 200 to 199.9 200 to 1999 mS/cm 2.00 to 19.99 20.0 to 19.99 20.0 to 199.9 200 to 2000	μS/m 0.0 to 19.9 20.0 to 199.9 200 to 1999 mS/m 2.00 to 19.99 20.0 to 199.9 200 to 1999 S/m 2.00 to 19.99 20.0 to 19.99 20.0 to 200.0		
Resolution	Auto ranging, up to	4 significant digits		
Accuracy	> 200 mS/cm (20.0 S	le of each range ;/m): ± 1.5% full scale		
Reference Temperature		30°C		
Temperature Coefficient		).00 %/°C		
Calibration Points Units	1 , 7	Up to 5 (Manual)		
		n, S/m 80.00 ppt		
Salinity Range		8.000 %		
Resolution	0.01 ppt, 0.001 %			
Accuracy	•	1 ppt, whichever is greater		
Salinity Curves	,	(UNESCO 1978)		
Calibration Option	Y	es		
Total Dissolved Solids (TDS) Range	Ŭ	200,000 mg/L		
Resolution		significant digits mg/L, whichever is greater		
Accuracy TDS Curves	~	EN27888, 442, NaCl		
Resistivity Range	$\begin{array}{c} \Omega \bullet \text{cm} \\ \Omega \bullet \text{cm} \\ 0.1 \text{ to } 199.9 \\ 200 \text{ to } 199.9 \\ 200 \text{ to } 199.9 \\ 2.00 \text{ to } 19.99 \\ 20.0 \text{ to } 199.9 \\ 200 \text{ to } 199.9 \\ 200 \text{ to } 199.9 \\ 2.00 \text{ to } 199.9 \\ 2.00 \text{ to } 19.99 \\ 2.00 \text{ to } 19.99 \\ 2.00 \text{ to } 19.99 \\ 20.0 \text{ to } 200.0 \end{array}$	$\begin{array}{c} \Omega \bullet m \\ \Omega \bullet m \\ 0.001 \text{ to } 1.999 \\ 2.00 \text{ to } 1.999 \\ 20.0 \text{ to } 199.9 \\ 200 \text{ to } 199.9 \\ 200 \text{ to } 1999 \\ k\Omega \bullet m \\ 2.00 \text{ to } 19.99 \\ 20.0 \text{ to } 19.99 \\ 20.0 \text{ to } 199.9 \\ 200 \text{ to } 2000 \end{array}$		
Resolution	0 0. 1	4 significant digits		
Accuracy	> 200 mS/cm (20.0 S	le of each range ;/m): ± 1.5% full scale		
Temperature Range	-30.0 to 130.0 °C -22.0 to 266.0 °F			
Resolution	0.1 °C / °F			
Accuracy Calibration Option	±0.5 °C / ±0.9 °F			
Calibration Option Body Material	Yes ABS / Polycarbonate			
Length and Diameter		30 mm		
Connector		n-pull		
Cable Length	2 m	5 m		
5				

4-Cell Conductivity Sensor Cartridge				
Model	300-4C-C			
Part No.	3200780928			
Cell Constant	0.172 cm <sup>-1</sup>			
Conductivity Range	10 µS/cm to 2000 mS/cm			
Operating Temperature Range	0 to 100 °C 32.0 to 212.0 °F			
Temperature Sensor	Built-in			
Length and Diameter	110 x 16 mm			
Body Material	Epoxy, carbon			


2-Cell Conductivity Sensor Cartridge				
Model	300-2C-C			
Part No.	3200820579			
Cell Constant	0.1 cm <sup>-1</sup>			
Conductivity Range	0.01 µS/cm to 500 µS/cm			
Temperature Range	0 to 100 °C 32.0 to 212.0 °F			
Temperature Sensor	Built-in			
Length and Diameter	110 x 16 mm			
Body Material	Stainless steel			

Ion Sensor Head Ion / mV / Temp (°C/°F)				
Model	300-1-2	300-1-5		
Part No.	3200812203	3200923560		
Ion Range	(mg/L, mmol/L) 0.000 to 0.999, 1.00 to 9.99, 10.0 to 99.9, 100 to 999, 1000 to 9990, 10000 to 99900			
Resolution	0.001 minimum, 3 significant digits			
Accuracy	±0.3% full scale			
Calibration Points	Up to 5			
mV Range	±1000.0 mV			
Resolution	0.1 mV			
Accuracy	±0.1 mV			
Temperature Range	-30.0 to 130.0 °C -22.0 to 266.0 °F			
Resolution	0.1 °C / °F			
Accuracy	±0.5 °C / ±0.9 °F			
Calibration Option	Yes			
Body Material	ABS / Polycarbonate			
Length and Diameter	85 x 30 mm			
Connector	Push-pull			
Cable Length	2 m	5 m		

ORP Sensor Head ORP / Temp (°C/°F)					
Model	300-0-2	300-0-5			
Part No.	3200812204	3200923561			
ORP Range	-2000 to -	+2000 mV			
Resolution	< ±1000.0 mV: 0.1mV ≥ 1000.0 mV: 1 mV				
Accuracy	< ±1000.0 mV: ±0.1 mV ≥ 1000.0 mV: ±1 mV				
Calibration Option	Yes				
Temperature Range	-30.0 to 130.0 °C -22.0 to 266.0 °F				
Resolution	0.1 °C / °F				
Accuracy	±0.5 °C / ±0.9 °F				
Calibration Option	Yes				
Body Material	ABS / Polycarbonate				
Length and Diameter	85 x 3	30 mm			
Connector	ector Push-pull				
Cable Length	2 m	5 m			

ORP Sensor Cartridge ORP / Temp (°C/°F)				
Model	300-O-C			
Part No.	3200922104			
ORP Range	-2000 to +2000 mV			
Temperature Range	0 to 80 °C -32.0 to 176.0 °F			
Junction Material	Porous sintered polyethylene			
Double Junction	Yes			
Temperature Sensor	Buit-in			
Length and Diameter	110 x 16 mm			
Body Material	Polycarbonate, platinum/glass			

pH / ORP / Conductivity Sensor Head





501-S NIST pH Buffer Solution Kit



502-S USA pH Buffer Solution Kit



503-S Conductivity Standard Solution Kit




500-225 ORP Standard Solution 225 mV



230 Cleaning Solutions

pH Buffer Solution Kits				
Model	Part No.	Description	Volume	
501-S	3999960015	NIST pH Buffer Solution Kit (pH 4.01, 6.86, 9.18 buffers & 3.33M KCI)	250ml each	
502-S	3999960016	USA pH Buffer Solution Kit (pH 4.01, 7.00, 10.01 buffers & 3.33M KCI )	250ml each	
pH Buffer Solu	utions			
Model	Part No.	Description	Volume	
500-2	3999960028	pH 1.68 Buffer Solution at 25°C	500ml	
500-4	3999960029	pH 4.01 Buffer Solution at 25°C	500ml	
500-686	3999960030	pH 6.86 Buffer Solution at 25℃	500ml	
500-7	3999960031	pH 7.00 Buffer Solution at 25°C	500ml	
500-9	3999960032	pH 9.18 Buffer Solution at 25°C	500ml	
500-10	3999960033	pH 10.01 Buffer Solution at 25°C	500ml	
500-12	3999960034	pH 12.46 Buffer Solution at 25°C	500ml	
Conductivity S	Standard Solut	ion Kit		
Model	Part No.	Description	Volume	
503-S	3999960017	Conductivity Standard Solution Kit (84µS/cm, 1413µS/cm, 12.88mS/cm & 111.8mS/cm)	250ml each	
Conductivity S	Standard Solut	ions		
Model	Part No.	Description	Volume	
500-21	3999960035	84 µS/cm Conductivity Standard Solution	500ml	
500-22	3999960036	1413 µS/cm Conductivity Standard Solution	500ml	
500-23	3999960037	12.88 mS/cm Conductivity Standard Solution	500ml	
500-24	3999960038	111.8 mS/cm Conductivity Standard Solution	500ml	
ORP Standard	Solution & Po	wders		
Model	Part No.	Description	Volume	
500-225	4000047848	ORP Standard Solution 225 mV at 25°C	500ml	
160-51	3200043618	ORP Powder 89 mV at 25°C (for 250ml solution)	10 sachets/pack	
160-22	3200043617	ORP Powder 258 mV at 25°C (for 250ml solution)	10 sachets/pack	
pH/ORP Electi	rode Filling Sol	utions		
Model	Part No.	Description	Volume	
525-3	3999960023	3.33M KCI	250ml	
300	3200043640	3.33M KCI	250ml	
pH Electrode	Cleaning Soluti	ions		
Model	Part No.	Description	Volume	
220	3014028653	For removing inorganic residues from glass membrane and liquid junction	2 x 50ml	
230	3200530494	For removing inorganic and organic residues from glass membrane (30ml Solution A & 100ml Solution B)	30ml & 100ml	
250	3200366771	For removing protein residues from glass membrane and liquid junction	400ml	



Ion Standard Solutions					
Model	Part No.			Description	Volume
500-NH4-SH	3200697171 1		1000 mg/L Ammonium Ion Standard Solution		500ml
500-NH4-SL	3200697172		100 m	g/L Ammonium Ion Standard Solution	500ml
500-CA-SH	3200	697175	1000 r	ng/L Calcium Ion Standard Solution	500ml
500-CA-SL	3200	697176	100 m	g/L Calcium Ion Standard Solution	500ml
500-CL-SH	3200	697167	1000 r	ng/L Chloride Ion Standard Solution	500ml
500-CL-SL	3200	697168	100 m	g/L Chloride Ion Standard Solution	500ml
500-F-SH	3200	697163	1000 r	ng/L Fluoride Ion Standard Solution	500ml
500-F-SL	3200	697164	100 m	g/L Fluoride Ion Standard Solution	500ml
500-NO3-SH	3200	697179	1000 r	ng/L Nitrate Ion Standard Solution	500ml
500-NO3-SL	3200	697180	100 m	g/L Nitrate Ion Standard Solution	500ml
500-K-SH	3200	697183	1000 r	ng/L Potassium Ion Standard Solution	500ml
500-K-SL	3200	697184	100 m	g/L Potassium Ion Standard Solution	500ml
Ionic Strength Adjustors					
Model	Part No.			Description	Volume
500-NH3-ISA	3200697174		Ammo	nia Ionic Strength Adjustor 🛛 📀 🚸	500ml
500-CA-ISA	3200697178		Calciu	m Ionic Strength Adjustor	500ml
500-CL-ISA	3200	697170	Chlori	Chloride Ionic Strength Adjustor 500	
500-F-TISAB	3200	697166	Fluorio	Fluoride Ionic Strength Adjustor 500m	
500-NO3-ISA	3200697182		Nitrate Ionic Strength Adjustor		500ml
500-K-ISA	3200	697186	Potassium Ionic Strength Adjustor		500ml
Ion Selective I	Electro	de Fillin	g Solu	tions	
Model	Pai	rt No.		Description	Volume
500-NH3-IFS	3200	697173	Ammo	nia Electrode Filling Solution	500ml
500-CA-IFS	3200	697177	Calciu	m Electrode Filling solution	500ml
500-CL-IFS	3200	697169	Chlori	de Electrode Filling Solution	500ml
500-F-IFS	3200	697165	Fluorio	le Electrode Filling Solution	500ml
500-NO3-IFS	3200	697181	Nitrate	Electrode Filling Solution	500ml
500-K-IFS	3200	697185	Potass	sium Electrode Filling Solution	500ml
Accessories					
Model		Part No.		Description	
3200		320086	1022	Integrated Electrode Stand (Height: 383mm) for 2000 Series	
FA-70A 3200		320064	14455 Integrated Electrode Stand (Height: 338mm) for 1000 Series		1000 Series
FA-70S 32		3200382557		Adjustable, free-standing electrode stand (Height: 384 mm)	
FA-70L		320038	2560	Long, free-standing electrode stand (Height: 450-650mm)	
DP-70S 320052		320052	8474	Electrode stand for 100 Series and D-70, ES-70, OM-70 Series handheld meters (Height: 400mm)	

	1 A-70A	3200044433	Integrated Electiode Stand (Height: SS61111) for 1000 Series
	FA-70S	3200382557	Adjustable, free-standing electrode stand (Height: 384 mm)
	FA-70L	3200382560	Long, free-standing electrode stand (Height: 450-650mm)
	DP-70S	3200528474	Electrode stand for 100 Series and D-70, ES-70, OM-70 Series handheld meters (Height: 400mm)
	- 212	3200373991	Arm for electrode stand FA-70A, FA-70S, & FA-70L
		3200373961	Electrode holders, 2pcs/pack (for mounting electrode with round cap on electrode stand arm)
		3200382482	Electrode protection caps, 5pcs/pack
		3200382482	Electrode protection cap for long electrode (for 9680S-10D, 9480-10C pH Electrode)
		3200044409	Clear pH sensor tip guard (for plastic pH electrodes 9651/9652, 9625, 9630 etc.), 5pcs/pack
		3200828646	Black pH sensor tip guard (for 200 series, 300 series), 3pcs/pack
	Stirrer set	3200897227	Stirrer set for 9521-10D
		3200779640	Electrode adapter
)	300-BNC	3200821465	Sensor head adapter (for WQ-300 pH/ORP/Ion sensor heads)
	Flow Cell	3200844642	Glass flow cell (for 300-2C-C and 9371-10D)
	<b>300-EXT-10</b>	3200921588	10-m sensor head extension cable

### **Technical Tip**

# pH Electrode Care and Maintenance Procedures

Your pH electrode will eventually reach the end of its useful life as its performance naturally degrades over time. To maximize the performance of your pH electrode and extend its life span, proper care and regular maintenance are equally required.




- Part no. 3014028653 Cleaning Solution 220 - contains 10% thiourea and 1% hydrochloric acid (HCI) for removing inorganic residues on glass membrane and junction
- Part no. 3200366771 Cleaning Solution 250 contains less than 0.5% enzyme protease, less than 0.1% sodium azide, and other ingredients (See SDS) for removing protein residues on glass membrane and junction
- Part no. 3999960023 525-3 3.33M KCl pH electrode filling solution (for liquid-filled electrodes)






Refilling

The pH electrode may be filled with either an ionic liquid solution (refillable or liquid-filled pH electrode) or ionic gel solution (sealed or gel-filled pH electrode). Gel-filled pH electrodes do not require routine refilling and typically require less maintenance than liquid-filled electrodes. Liquid-filled pH electrodes are constructed with refilling port, which is securely covered with a slider. The refilling port allows you to fill or empty the reference chamber.



- To top up or re-fill the reference chamber of liquid-filled pH electrode, push the slider upward to uncover the refilling port and insert a dropper containing fresh 3.33M potassium chloride (KCI) solution. The filling solution should reach the bottom of the refilling port.
- The filling solution level must be maintained just below the refilling port and higher than the pH buffer or sample level during calibration and measurement. This creates a positive head pressure forcing the filling solution to leak into pH buffer or sample through the junction and preventing the reverse.
- Bubbles may form and get trapped within the solution of the sensing tip or reference chamber during transportation. This can affect the operation of your pH electrode. To dislodge the bubbles, gently shake the electrode body.
- If the filling solution inside the reference chamber gets contaminated with sample



distilled or deionized water) in a squirt bottle

> or microbial growth or the reading is drifting, change the filling solution. Tilt the pH electrode, uncover the refilling port, and draw out the old solution using a dropper before refilling it with fresh 3.33M KCI solution.

SDSs of HORIBA solutions at www.horiba-

#### Conditioning


Nowadays, combination and 3-in-1 pH electrodes are commonly available. Both types of pH electrodes consist of glass electrode and reference electrode built in one body, but the latter is integrated with temperature sensor for detecting the temperature of the solution being measured.

laqua.com

The glass electrode has a silver-based electrical wire suspended in a neutral solution with KCI contained inside a special glass. The surface of the glass bulb or membrane at the tip of the electrode must be hydrated to function properly. This can be accomplished by immersing the glass membrane in an aqueous solution, where a hydrated layer that is responsible for the pH response of the glass, is developed.

Another component of the pH electrode that must remain hydrated is the junction of the reference electrode. The junction is made of porous material such as ceramic or sintered polyethylene, which allows filling solution of the electrode to leak into the solution being measured. Keeping the reference junction hydrated will prevent precipitation of KCI from the filling solution which may clog it and cause erratic or slow electrode response.

All pH electrodes come with white protective cap. A sponge wet with pure water is positioned at the bottom of the cap to keep the glass membrane and junction moist. If you find KCI salts formed on the junction or refilling port of your pH electrode, simply rinse off using clean water. This KCl creep from the filling solution is normal.



- A dry pH electrode will give inaccurate reading in pH measurement. Condition a dry pH electrode by soaking the glass membrane and junction in pH 7.00, 4.01 buffer, or tap water for at least 1 hour to regenerate the hydrated layer. Note: High salt solutions such as 3.33M KCl and the like are not recommended for conditioning our pH electrodes. After conditioning, rinse the pH electrode with clean water and proceed with calibration.
- Never touch the glass membrane with fingers as oil or dirt may coat the glass and interfere with measurement.

If a liquid-filled pH electrode is in use,

the refilling port must be uncovered and

the filling solution level must be higher

than the pH buffer or sample level. These

conditions will ensure smooth outward flow

of filling solution through the junction during

Before and after measurement, rinse the

pH electrode with clean water and/or with a

portion of the next solution to be measured and blot with soft lint-free tissue to remove

excess water or solution. Rinsing between measurements prevents contamination by

carry-over on the electrodes. Avoid wiping

or rubbing as this can scratch the glass membrane, remove the hydrated layer, and

cause static charge, resulting in inaccurate

Calibrate frequently using at least two fresh pH buffers that bracket the expected


sample pH value. Make sure that the glass

membrane and junction of pH electrode are both immersed in pH buffer or sample.


calibration and measurement.

pH readings.

#### **Calibration and Measurement**








Temperature connector




- To compensate for temperate effect on pH, use either 3-in-1 pH electrode or combination pH electrode and temperature probe. If temperature probe is not available, check the solution temperature using a calibrated thermometer and input the reading into the meter.
- Stir pH buffers and sample at same rate. Stirring provides representative pH value of a solution and faster electrode response. If stirring is not possible due to measurement noise, limited sample volume or other reasons, it may be abandoned in both calibration and measurement.
- There is a wide selection of pH electrodes and each model is designed to suit specific applications. Choose the best pH electrode suitable for your sample.

#### Cleaning

A clean, hydrated glass membrane and free-flowing junction are necessary in performing an accurate measurement of pH. The choice of cleaning solution should effectively remove all contaminants based on sample tested without damaging your pH electrode.





- 1. If the pH electrode is liquid-filled, uncover the refilling port.
  - Clean the tip of your pH electrode using the appropriate cleaning solution. Make sure that the glass membrane and junction are both immersed in cleaning solution.
    - General samples Soak the pH electrode in diluted detergent solution for 5 to 10 minutes, while moderately stirring the solution. A strong cleaning solution is needed for clogged junction, stains, and electrodes exhibiting slow response. Soak the pH electrode in cleaning solution 220 or 0.1M HCl for at least 1 hour.
    - Oily samples Soak the pH electrode in warm, diluted detergent solution for 5 to 10 minutes, while moderately stirring the solution. Alternatively, rinse the pH electrode with methanol or ethanol. Note: Alcohol is only applicable for glass-body electrodes. Never use organic solvents such as alcohol, acetone etc. to clean any plastic-body electrode as they may damage the body and shorten the life span. Use of organic solvents will void the electrode warranty.
    - Protein-containing samples soak the pH electrode in cleaning solution 250 for at least 1 hour.
  - Rinse the pH electrode with clean water.
- If the pH electrode is liquid-filled, draw out the old filling solution from the reference chamber and refill it with fresh 3.33M KCI (See Refilling).
- 5. Condition the pH electrode (See Conditioning).

If calibration with fresh pH buffers failed repeatedly and cleaning failed to restore the performance, replace the pH electrode with a new one.

#### Storage

pH electrodes must be clean before they are stored for any length of time.

- 1. If the pH electrode is liquid-filled, cover the refilling port with the slider to prevent evaporation of filling solution.
  - Wash the protective cap with clean water to wet the sponge and remove KCI salts.
  - Insert the pH electrode into the protective cap with wet sponge. The water will not dissipate easily as the cap fit snugly on the electrode body. This environment is enough to keep the glass membrane and junction moist. It is not necessary to fill the cap with clean water and soak the pH electrode tip.

#### Short-term storage:

Between measurements, the pH electrode can be soaked in pH 7.00 buffer or clean water (e.g., tap, distilled or deionized).


#### Visit HORIBA's website!

#### Water Quality Analyzers

www.horiba-laqua.com

With over 60 years of engineering excellence, HORIBA's diverse range of water quality analyzers and electrodes are ideal for everyday laboratory needs through to the most demanding of applications. Visit our website for a wealth of useful information and water quality measurement tips to help you obtain the best results in your work.







#### **Benchtop Meters**

Developed using extensive feedback from users, our new LAQUA meters deliver the best solution for water quality analysis. Our LAQUA website features an online 'Selection Guide' to enable you to find the perfect LAQUA meter and electrode for your need.



#### **Handheld Meters**

In the lab, in the field or anywhere you need it. LAQUA Handheld meters are designed for use with one hand and with an IP67 waterproof rating and shock-resistant casing. Meters can be used for long periods, even in dark places, making it ideal for field measurements in rivers and lakes.



#### **Pocket Meters**

Analyzing water quality is simplified when using our LAQUAtwin range of meters. Designed to produce accurate and reliable results. Anyone, anywhere, at any time can measure samples easily with a LAQUAtwin meter. See just how good they are at our website.

25



LAQUAtwin pocket meters offer quick and convenient alternative to analyze important parameters with high accuracy. Several application notes are available at (http://goo.gl/znwE6j) detailing the use of LAQUAtwin and the results achieved for the respective applications. Additional application notes will be added when available.



Visit the HORIBA LAQUA Singapore Channel on YouTube and subscribe to see more of our videos.

#### S CE





HORIBA UK Limited

Kyoto Close, Moulton Park, Northampton NN3 6FL

Phone: +44 1604 642500



- · The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.
- It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.
- · Windows is a registered trademark of Microsoft Corporation in the United States and other countries.

#### Asia Pacific

HORIBA Instruments (Singapore) Pte. Ltd. 83 Science Park Drive, #02-02A, The Curie, Singapore 118258 Phone: 65 6908-9660 Fax: 65 6745-8155 e-mail: laqua@horiba.com

Europe, Middle East, & Africa Americas **HORIBA Instruments Incorporated** e-mail: waterguality@horiba.com

9755 Research Drive, Irvine, California 92618 USA Phone: +1 949 250 4811 Fax: +1 949 250 0924, +1 949 468 1890 e-mail: labinfo@horiba.com

Brochure HEA-10-2019E

HORIBA

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific