

HORIBA
Scientific

LAQUA

LAQUA*twin* Ca-11C Waterproof Pocket Bovine Blood Calcium Meter



MADE IN
JAPAN

www.horiba-laqua.com

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

LAQUAtwin

Compact Bovine Blood iCa Checker

Calcium is an important macro mineral for dairy cows. Measurement of ionized calcium (iCa) on site can be of great value to veterinarians and producers. An innovative and lower cost instrument could help with more timely decisions and avoid submission of blood samples to a laboratory. The dairy veterinary field is progressively studying subclinical hypocalcemia. To detect this condition real time using the *LAQUAtwin Ca-11C can make a significant improvement in an animal's health.

*LAQUAtwin Ca-11C is not a medical device.

On-Site Measurements

Collect blood from the animal and place a drop or two on the sensor.

Whole Blood or Serum

The meter can be used to measure not only whole blood but also serum.

Instant Results

The measured ionized calcium value should give a good indication of the animals' health and probability of contracting hypocalcemia or hypercalcemia.

Fast Screening

Short analysis on-site ensures quicker response to the animal's healthcare.

Economical

Monitor animal's health at a fraction of the cost of available methods employing commercial blood analyzers.

Replaceable Sensor

The flat sensor is easy to clean, which is important when handling samples containing proteins. It will last for at least 1,000 measurements and can simply be replaced.

Indicators

Various indicators appear such as a smiley once the value is stable, low battery indicator and measurement mode.

Backlight

For easy viewing in the dark. It also flashes when the batteries run low.

3-Button Operation

The instrument is so easy to use that only 3 buttons are required for operation

IP67 Water- and Dustproof Housing

The instrument is completely water – and dustproof which ensures durability and a long lifetime even in the harshest conditions.

IP67


Accuracy in a single drop

Incorporating the same parts as standard laboratory electrodes, the LAQUAtwin compact meters are built with miniaturised components and unique flat sensor chip, which is less than 1 mm thick, a result of 60 years of HORIBA's sensor engineering technology.

Reference electrode
Liquid junction
Ion selective membrane
Ag/AgCl
Temperature sensor
Thickness Approx. 0.8 mm

Easy Measurement

Calibrate the meter with standard solutions then measure the blood sample.



Scan code to view video online



Complete Set

Perform your measurements wherever you want, whenever is needed.



Package includes:

- LAQUAtwin Ca-11C
- Y052L 1.25mmol/L Standard solution (25ml x 2)
- Y052H 2.50mmol/L Standard solution (25ml x 2)
- 251 Cleaning solution (250ml) x 3
- CR2032 batteries x 2
- Wash bottle
- Instruction & quick manuals
- Carrying case


Meter Specifications

Model	LAQUAtwin Ca-11C
Target Ion	Calcium (Ca ²⁺)
Measurement Principle	Ion Selective Electrode
Sample Volume	0.3 ml (Minimum)
Measurement Range	0.1 to 5.0 mmol/L
Resolution	0.1 or 0.01 mmol/L
Accuracy	± 20% of reading value
Calibration	2 points (1.25 & 2.50 mmol/L)
Recommended Sample Temperature	20 to 25°C
Operating Environment	5 to 40°C; 85% or less relative humidity (no condensation)
IP Rating	IP67 water- and dustproof
Display	Custom (monochrome) digital LCD with backlight
Power & Battery Life	2 x CR2032, 150 hours
Dimensions	164x29x20mm (meter), 274x235x75mm (carrying case)
Weight	50g (meter), 530g (complete set)

Ordering Information

Part No.	Model	Description
3200772940	LAQUAtwin Ca-11C	Includes meter, sensor, batteries, Y052L, Y052H, 251 x 3, wash bottle, manuals, & carrying case
3200772945	S051	Calcium ion sensor
3200774921	Y053	Solution kit – Y052L, Y052H, & 251 x 3
3200774598	Y052L	1.25 mmol/L standard solution (25ml x 2)
3200774599	Y052H	2.50 mmol/L standard solution (25ml x 2)
3200774601	251	Cleaning solution (250ml)

Measurement results shows high correlation with NIST standard.



*An internal investigation



Scan code to view app note



Scan code to view LAQUAtwin brochure




Visit HORIBA's website!

www.horiba-laqua.com

Water Quality Analyzers

With over 60 years of engineering excellence, HORIBA's diverse range of water quality analyzers and electrodes are ideal for everyday laboratory needs through to the most demanding of applications. Visit our website for a wealth of useful information and water quality measurement tips to help you obtain the best results in your work.



Benchtop Meters

Developed using extensive feedback from users, our new LAQUA meters deliver the best solution for water quality analysis. Our LAQUA website features an online 'Selection Guide' to enable you to find the perfect LAQUA meter and electrode for your need.



Handheld Meters

In the lab, in the field or anywhere you need it. LAQUA Handheld meters are designed for use with one hand and with an IP67 waterproof rating and shock-resistant casing. Meters can be used for long periods, even in dark places, making it ideal for field measurements in rivers and lakes.



Electrodes

Various electrodes to match any application. A wide range of products for both benchtop and portable systems are available, including easy and reliable standard models, application-focused models for small samples or large containers, and special electrodes for specific sample characteristics.



Pocket Meters

Analyzing water quality is simplified when using our LAQUAtwin range of meters. Designed to produce accurate and reliable results. Anyone, anywhere, at any time can measure samples easily with a LAQUAtwin meter. See just how good they are at our website.



Application Notes

LAQUAtwin pocket meters offer quick and convenient alternative to analyze important parameters with high accuracy. Several application notes are available at (<http://goo.gl/znwE6j>) detailing the use of LAQUAtwin and the results achieved for the respective applications. Additional application notes will be added when available.



Visit the **HORIBA LAQUA Singapore** Channel on YouTube and subscribe to see more of our videos.



RoHS

- The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.
- It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.
- Windows is a registered trademark of Microsoft Corporation in the United States and other countries.



Asia Pacific
HORIBA Instruments (Singapore) Pte. Ltd.
 83 Science Park Drive, #02-02A,
 The Curie, Singapore 118258
 Phone: 65 6908-9660
 Fax: 65 6745-8155
 e-mail: laqua@horiba.com

Europe, Middle East, & Africa
HORIBA UK Limited
 Kyoto Close, Moulton Park,
 Northampton NN3 6FL
 Phone: +44 1604 642500
 e-mail: waterquality@horiba.com

Americas
HORIBA Instruments Incorporated
 9755 Research Drive, Irvine,
 California 92618 USA
 Phone: +1 949 250 4811
 Fax: +1 949 250 0924, +1 949 468 1890
 e-mail: labinfo@horiba.com

Brochure LPBM-11-2021