

# Instruction Manual COND METER (DS-72G)





This manual describes the operation of the following instrument.

Brand (pet name):LAQUASeries name:Benchtop pH/Water Quality AnalyzerModel:DS-72GModel description:COND METER

Be sure to read this manual before using the product to ensure proper and safe operation of the instrument. Also safely store the manual so it is readily possible whenever necessary.

Product specifications and appearance, as well as the contents of this manual are subject to change without notice.

#### • Warranty and responsibility

HORIBA Advanced Techno Co., Ltd. warrants that the Product shall be free from defects in material and workmanship and agrees to repair or replace free of charge, at option of HORIBA Advanced Techno Co., Ltd., any malfunctioned or damaged Product attributable to responsibility of HORIBA Advanced Techno Co., Ltd. for a period of one (1) year from the delivery unless otherwise agreed with a written agreement. In any one of the following cases, none of the warranties set forth herein shall be extended;

- Any malfunction or damage attributable to improper operation
- Any malfunction attributable to repair or modification by any person not authorized by HORIBA Advanced Techno Co., Ltd.
- Any malfunction or damage attributable to the use in an environment not specified in this manual
- Any malfunction or damage attributable to violation of the instructions in this manual or operations in the manner not specified in this manual
- Any malfunction or damage attributable to any cause or causes beyond the reasonable control of HORIBA Advanced Techno Co., Ltd. such as natural disasters
- Any deterioration in appearance attributable to corrosion, rust, and so on
- · Replacement of consumables

HORIBA Advanced Techno Co., Ltd. SHALL NOT BE LIABLE FOR ANY DAMAGES RESULTING FROM ANY MALFUNCTIONS OF THE PRODUCT, ANY ERASURE OF DATA, OR ANY OTHER USES OF THE PRODUCT.

#### **Trademarks**

Company names and brand names are either registered trademarks or trademarks of the respective companies. (R), (TM) symbols may be omitted in this manual.



EU regulations

Conformable standards

This equipment conforms to the following standards:

C	E

EMC:EN61326-1<br/>Class B, Basic electromagnetic environmentSafety:EN61010-1RoHS:EN50581<br/>9. Monitoring and control instruments

**Warning:** This product is not intended for use in industrial environments. In an industrial environment, electromagnetic environmental effects may cause the incorrect performance of the product in which case the user may be required to take adequate measures.

#### Installation environment

- This product is designed for the following environment.
- Overvoltage category II
- Pollution degree 2

#### Information on disposal of electrical and electronic equipment and disposal of batteries and accumulators


The crossed out wheeled bin symbol with underbar shown on the product or accompanying documents indicates the product requires appropriate treatment, collection and recycle for waste electrical and electronic equipment (WEEE) under the Directive 2012/19/EU, and/or waste batteries and accumulators under the Directive 2006/66/EC in the European Union.

The symbol might be put with one of the chemical symbols below. In this case, it satisfies the requirements of the Directive 2006/66/EC for the object chemical.

This product should not be disposed of as unsorted household waste.

Your correct disposal of WEEE, waste batteries and accumulators will contribute to reducing wasteful consumption of natural resources, and protecting human health and the environment from potential negative effects caused by hazardous substance in products.

Contact your supplier for information on applicable disposal methods.



#### Authorised representative in EU

HORIBA UK Limited

2 Dalston Gardens, Stanmore, Middx HA7 1BQ, UK

### • FCC rules

Any changes or modifications not expressly approved by the party responsible for compliance shall void the user's authority to operate the equipment.

#### Warning

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

# Korea certification

# ● B급 기기 (가정용 방송통신기자재)

이 기기는 가정용(B 급) 전자파적합기기로서 주로 가정에서 사용하는 것을 목적으로 하며, 모든 지역에서 사용할 수 있습니다.

# Taiwan battery recycling mark



#### China regulation

标记的意义 Meaning of Marking マークの意味

> 本标记适用在中华人民共和国销售电器电子产品,标记中央的数 字表示环境保护使用期限的年数。(不是表示产品质量保证期间。) 只要遵守这个产品有关的安全和使用注意事项,从制造日开始算 起在这个年限内,不会给环境污染、人体和财产带来严重的影响。 请不要随意废弃本电器电子产品。



This marking is applied to electric and electronic products sold in the People's Republic of China. The figure at the center of the marking indicates the environmental protection use period in years. (It does not indicate a product guarantee period.) It guarantees that the product will not cause environment pollution nor serious influence on human body and property within the period of the indicated years which is counted from the date of manufacture as far as the safety and usage precautions for the product are observed. Do not throw away this product without any good reason.

本マークは、中華人民共和国で販売される電気電子製品に適用さ れ、マークの中央の数字は環境保護使用期限の年数を意味します (製品の品質保証期間を示すものではありません)。この製品に関 する安全や使用上の注意をお守り頂く限り、製造日から起算するこ の年限内では、環境汚染や人体や財産に深刻な影響を及ぼすこと はありません。本製品をみだりに廃棄しないでください。

# 产品中有害物质的名称及含量

#### Name and amount of hazardous substance used in a product

		有害物质 Hazardous substance						
部件名称 Unit name	铅 Lead (Pb)	汞 Mer− cury (Hg)	福 石 Cad- mium (Cd)	六价铬 Hexa- valent chrom- ium (Cr (VI))	Se 多溴联苯 Poly bromo- biphenyl (PBB)	多溴二苯醚 Poly bromo- diphenyl ether (PBDE)		
本体 Main unit	×	0	0	0	$\bigcirc$	$\bigcirc$		
电池 Battery	×	0	0	0	0	$\bigcirc$		
AC适配器 AC adapter <sup>*1</sup>	×	0	0	0	0	0		
电缆 Cable	×	0	0	0	0	0		
支架 Stand <sup>*2</sup>	0	0	0	0	$\bigcirc$	0		
打印机 Printer <sup>*2</sup>	×	0	0	0	0	0		
电极 Electrode <sup>*2</sup>	×	$\bigcirc$	×	0	$\bigcirc$	0		

本表格依据 SJ/T 11364 的规定编制。

This form is prepared in accordance with SJ/T 11364.

○: 表示该有害物质在该部件所有均质材料中的含量均在 GB/T 26572 规定的限量要 求以下。

Denotes that the amount of the hazardous substance contained in all of the homogeneous materials used in the component is below the limit on the acceptable amount stipulated in the GB/T 26572.

×: 表示该有害物质至少在该部件的某一均质材料中的含量超出 GB/T 26572 规定的 限量要求。

Denotes that the amount of the hazardous substance contained in any of the homogeneous materials used in the component is above the limit on the acceptable amount stipulated in the GB/T 26572.


- \*1: 本部件的环保使用期限为10年。 The environmental protection use period of this product is 10 years.
- \*2: 选配件 Optional products

# For Your Safety

#### Hazard classification and warning symbols

Warning messages are described in the following manner. Read the messages and follow the instructions carefully.

#### Hazard classification



#### • Warning symbols



Description of what should be done, or what should be followed



Description of what should never be done, or what is prohibited

#### • [DEU] Sicherheitsinformation

Lesen Sie vor der Verwendung des Produkts unbedingt diese Anleitung, um den ordnungsgemäßen und sicheren Betrieb des Produkts zu gewährleisten. Bewahren Sie die Anleitung sicher auf, damit sie bei Bedarf jederzeit zur Hand ist.

Die Inhalt dieser Anleitung können ohne Vorankündigung geändert werden.

#### Installationsumgebung

Dieses Produkt ist nicht zum Gebrauch in industriellen Umgebungen, wie in EN61326-1 definiert, vorgesehen.

In einer industriellen Umgebung können die elektromagnetischen Störungen eventuell zu Produktfehlfunktionen führen. Um dieses Produkt unter solchen Umständen verwenden zu können, muss der Benutzer ggf. angemessene Maßnahmen ergreifen.

Das Produkt ist gemäß EN61010-1 für die folgende Umgebung vorgesehen.

- Überspannungskategorie II
- Verschmutzungsgrad 2

#### • [FRA] Informations de sécurité

Veillez à lire le présent manuel avant d'utiliser le produit de manière à garantir son utilisation correcte et sûre. De même, rangez le manuel dans un lieu sûr de manière à pouvoir vous y reporter lorsque cela est nécessaire.

Le contenu du présent manuel peut être modifié sans notification préalable.

#### Environnement d'installation

Ce produit n'est pas destinés à une utilisation dans des environnements industriels, tels que définis dans la norme EN61326-1.

Dans un environnement industriel, les interférences électromagnétiques peuvent entraîner un dysfonctionnement du produit. Pour utiliser le produit dans ce type d'environnements, l'utilisateur peut avoir à prendre des mesures appropriées.

Le produit est conçu pour l'environnement suivant, tel que défini dans la norme EN61010-1.

- Catégorie de surtension II
- Degré de pollution 2

#### • [ITA] Informazioni sulla sicurezza

Leggere attentamente questo manuale prima di utilizzare il prodotto al fine di utilizzarlo in modo sicuro e adeguato. Inoltre, conservare in un luogo sicuro il manuale per poterlo consultare se necessario.

Le contenuti di questo manuale sono soggetti a modifiche senza preavviso.

#### Ambiente di installazione

Questo prodotto non è stati progettati per essere utilizzati in ambienti industriali, secondo la norma EN61326-1.

In un ambiente industriale, le interferenze elettromagnetiche potrebbero causare un malfunzionamento del prodotto. Per utilizzare il prodotto in tali ambienti, all'utente potrebbe essere richiesto di adottare le contromisure necessarie.

Il prodotto è designato per il seguente ambiente, definito nello standard EN61010-1.

- Categoria di sovratensione II
- Livello di inquinamento 2

#### [SWE] Säkerhetsinformation

Se till att du läser denna handbok innan du börjar använda produkten för en korrekt och säker användning av den. Spara sedan handboken på en säker och lättåtkomlig plats så att du kan konsultera den när så behövs.

Innehållet i denna handbok kan komma att ändras utan föregående meddelande därom.

#### Installationsmiljö

Detta produkten är ej avsedda för användning i industriella miljöer enligt riktlinjerna i EN61326-1.

Om den används i industrimiljöer kan de elektromagnetiska störningarna orsaka tekniska fel hos produkten. Om produkten ska användas i sådana miljöer kan användaren behöva vidta lämpliga åtgärder för att lösa dessa problem.

Produkten är utformad för användning i följande miljöer, i enlighet med SS-EN 61010-1.

- Överspänningskategori II
- Föroreningsgrad 2

#### [SPA] Información de seguridad

Asegúrese de leer este manual antes de utilizar el producto para garantizar un uso correcto y seguro del mismo. Asimismo, guarde de forma segura el manual para que esté disponible siempre que sea necesario.

El contenido de este manual están sujetos a cambios sin previo aviso.

#### Entorno de instalación

Este producto está diseñado para su uso en entornos industriales, tal y como se define en EN61326-1.

En un entorno industrial, las interferencias electromagnéticas pueden provocar un funcionamiento incorrecto del producto. Para usar el producto en tales entornos, el usuario debe tomar las medidas adecuadas.

El producto se ha diseñado para el siguiente entorno, definido en EN61010-1.

- Categoría de sobretensión II
- Nivel de contaminación 2

#### [POL] Informacje dotyczące bezpieczeństwa

Przed przystąpieniem do użytkowania tego produktu należy dokładnie zapoznać się z niniejszą instrukcją, aby zapewniona była prawidłowa i bezpieczna eksploatacja produktu. Instrukcję przechowywać w bezpiecznym miejscu, aby w razie potrzeby była zawsze dostępna.

Treść niniejszej instrukcji może ulec zmianie bez wcześniejszego powiadomienia.

#### Środowisko instalacji

Ten produkt nie są przeznaczone do użytkowania w środowisku przemysłowym, zgodnie z definicją określoną w normie EN61326-1.

W środowisku przemysłowym zakłócenia elektromagnetyczne mogą powodować nieprawidłowe działanie produktów. Możliwe, że aby użytkować produkt w takich środowiskach, użytkownik będzie musiał podjąć stosowne środki zaradcze.

Produkt jest przeznaczony do użycia w poniższym środowisku zdefiniowanym w normie EN61010-1.

- Kategoria przepięciowa II
- Stopień zanieczyszczenia 2

#### [NLD] Veiligheidsinformatie

Lees deze handleiding voordat u dit product gebruikt zodat u het op de juiste manier en veilig kunt gebruiken. Bewaar de handleiding goed zodat u hem wanneer nodig kunt raadplegen.

De inhoud van deze handleiding kunnen zonder voorafgaande kennisgeving worden gewijzigd.

#### Installatieomgeving

Dit product is niet bedoeld voor gebruik in een industriële omgeving zoals gedefinieerd in EN 61326-1.

In een industriële omgeving kan de elektromagnetische interferentie de werking van dit product storen. Voor gebruik van het product in een dergelijke omgeving moet de gebruiker mogelijk maatregelen treffen om de storing te verhelpen.

Het product is ontworpen voor de volgende omgeving, gedefinieerd in EN 61010-1.

- Overspanningscategorie II
- Vervuilingsgraad 2

### ● [JPN] 安全情報

ご使用になる前に、本書を必ずお読みください。お読みになった後は 必要なときにすぐに取り出せるように大切に保管してください。

本書に記載されている内容は予告なく変更される場合があります。あ らかじめご了承ください。

#### ● 設置環境

本製品は、EN61326-1 で定義される工業環境で使用することを想定した製品ではありません。

工業環境においては、電磁妨害の影響を受ける可能性があり、その場 合には使用者が適切な対策を講ずることが必要となることがありま す。


本製品は、EN61010-1 で定義される以下の環境用に設計されています。

過電圧カテゴリー ||

• 汚染度 2

# Safety Precautions

This section provides precautions to enable you to use the product safely and correctly and to prevent injury and damage. The terms of DANGER, WARNING, and CAUTION indicate the degree of imminency and hazardous situation. Read the precautions carefully as it contains important safety messages.



0	Harmful chemicals Some ion electrodes are used with hazardous standard solutions. Handle them with care. If the internal solution comes in contact with the skin, wash it off immediately. If it gets into eyes, flush with plenty of water and then consult a doctor.
0	Harmful chemicals The internal solution of an electrode is highly concentrated potassium chloride (3.33 mol/L KCI). If the internal solution comes in contact with the skin, wash it off immediately. If it gets into eyes, flush with plenty of water and then consult a doctor.
0	<b>Broken glass</b> Broken glass may cause injury. The outer tube and tip of an electrode are made of glass. Handle them with care.

# **CAUTION**

O not use the cable of serial communication, USB, or AC adapter under wet or humid conditions.

Otherwise, it may cause a fire, electric shock, or breakage.

# Product Handling Information

### Operational precautions

- Only use the product including accessories for their intended purpose.
- Do not drop, crash, or give any physical impact on the instrument.
- Do not immerse the instrument into alcohol, organic solvent, strong acid, strong alkaline, or the like. The instrument body contains ABS resin, acrylic resin, and some rubber parts.
- If the instrument is dropped into water or gets wet, wipe it using soft cloth. Do not heat to dry it with a hair-dryer (or the like).
- Use fingers to press the operation keys or the touch panel. Do not use a hard object like a metal stick or rod.
- Be careful not to let water into the instruction inside. The instrument is not water-proof.
- To disconnect an electrode or interface cable, hold the connector and pull it off. If you pull at the cable, it may cause a breakage.
- The touch panel is capacitance-type. Make sure to turn OFF the power before cleaning the panel. If you wipe it with the power ON, it may cause instrument malfunction.
- RS-232C or USB communication between the instrument and a personal computer may fail because of environmental conditions, such as (radio/electromagnetic) noise.
- Make sure to use the provided power supply cable to power this product.

### Environmental conditions for use and storage

- Temperature: 0°C to 45°C
- Humidity: under 80% in relative humidity and free from condensation

#### Avoid the following conditions:

- Dusty environment
- Strong vibration
- Direct sunlight
- Corrosive gas environment
- Close to an air-conditioner
- Direct wind

#### • Transportation

When transporting the instrument, repackage it in the original package box. Otherwise, it may cause instrument breakage.

### Disposal

Standard solution used for the calibration must be under neutralized before the disposal. As for the disposal of the meter, treat it as an industrial waste.

# Manual Information

### Description in this manual

NOTE

This interprets the necessary points for correct operation and notifies the important points for handling the product.



This indicates the part where to refer for information.

— HINT!—

This indicates reference information.

### Original language

This is the English translation of an original Japanese document.

# Contents

■ Preface	
Regulations	
■ For Your Safety	
<ul> <li>Product Handling Information</li> </ul>	
<ul> <li>Manual Information</li> </ul>	
	XIV
Chapter 1 Overview	1
1.1 Description of Each Part	
1.1.1 Rear	
1.1.2 Display	
1.1.3 Left Side	
1.1.4 Right Side	
1.1.5 Accessories	
1.1.6 Identification of Manufacturing Date	
1.1.7 Operation Keys	
1.1.8 Icons (Icon Bar)	
1.1.9 Status Icons	
1.1.10 Meas Screen	
1.2 Basic Operation of Touch-Panel and Tou	-
1.3 Function and Operation of the Meas Scr	
1.4 Assembling the Electrode Stand	
1.5 Connecting the Electrode	
1.5.1 Electrode Connector	
1.5.2 Temperature Connector	
1.6 Connecting the Power Source	
1.7 Connecting the Printer	
1.8 Connecting the Personal Computer	
1.9 Turn ON the Power	15
	10
Chapter 2 Before Measurement (Meter SET)	
2.1 Meter SET Screen	
2.2 Auto Hold Setting	
2.3 Custom Setting	
2.4 Sample Name Setting	
2.5 Interval Memory Setting	
2.6 USB Memory Setting	
2.7 Printer Setting	
2.8 Screen Settings	
2.9 Sound Setting	
2.10 Language Setting	
2.11 Security Setting	
2.12 User Entry/Info Change/Delete	
2.13 Date Setting	
2.14 Analog Output Adjustment	
2.15 Temperature Sensor Calibration	

# Contents

2.16 Resetting to Factory Defaults	34
Chapter 3 COND (Conductivity) Measurement	35
3.1 COND Calibration	35
3.1.1 Automatic Calibration Setting	
3.1.2 Calibration of Standard Solution	
3.2 COND Measurement Setting	39
3.2.1 Cell Constant Setting	
3.2.2 COND Measurement Unit Setting	
3.2.3 Temperature Setting	
3.2.4 Temperature Conversion Function Setting	
3.2.5 Alarm Setting	
3.2.6 Electrode Model Setting	45
3.2.7 Electrode Lot No. Setting	46
3.3 COND Measurement	47
Chapter 4 SAL (Salinity) Measurement	48
4.1 Measurement Target Selection	48
4.2 SAL Calibration Setting	
4.3 SAL Measurement Setting	
4.4 SAL Measurement Unit Setting	
4.4 SAL measurement onit Setting	
4.6 Alarm Setting	
4.6.1 Input Upper or Lower Limit Values	
4.7 Electrode Model Setting	
4.8 SAL Measurement	52
Chapter 5 Resist (Resistivity) Measurement	53
5.1 Resist Measurement Setting	53
5.2 Resist Measurement Unit Setting	53
5.3 Temperature Setting	53
5.4 Alarm Setting	54
5.4.1 Input Upper or Lower Limit Values	54
5.5 Electrode Model Setting	55
5.6 Resist Measurement	
Chapter 6 TDS (Total Dissolved Solids) Measurement	56
6.1 TDS Measurement Setting	
6.2 TDS Measurement Mode Setting	
6.2.1 Input TDS Linear Value when Selecting LINEAR	
6.3 Temperature Setting	
6.4 Alarm Setting	
6.4.1 Input Upper or Lower Limit Values	58
6.5 Electrode Model Setting	58

# Contents


	6.6 TDS Measurement	59
Chapte	r 7 Application Mode	60
	7.1 Pharmacopeia Mode	60
	7.1.1 Shift to Pharmacopeia Mode	
	7.1.2 Measured by USP (Stage 1)	
	7.1.3 Measured by USP (Stage 2)	
	7.1.4 Measured by EP	63
	7.1.5 Measured by JP (OFF-LINE)	
	7.1.6 Measured by JP (0mL-10mL (in container))	
	7.1.7 Measured by JP (10mL- (in container))	
	7.1.8 Measured by PPRC (CP) (Stage 1)	
	7.1.9 Measured by PPRC (CP) (Stage 2)	
	7.1.10 Temperature and Conductivity Requirements	69
Chapte	r 8 Periodic Inspection Mode	70
	8.1 COND Periodic Inspection Mode Setting	70
	8.1.1 Pharmacopoeia Mode	70
	8.1.2 COND Checker (X-52) Mode	
	8.2 Comment Input	
Chante	r 9 Data	74
Onapic		
	9.1 Measured data_All	
	9.2 Deleting Saved Data	
	9.3 Measured data_latest50	
	9.4 Measured data_search	
	9.5 Copy all meas. Data	
	9.6 Delete all meas. Data	76
Chapte	r 10 Specifications	77
	10.1 Model Information	77
	10.2 Measuring Object	
	10.3 Default Settings	
	10.3.1 Meter Default Settings	
	10.3.2 Measurement Condition Default Settings	
	(Can Be Set per Operator)	80
	10.4 Options	82

# **Chapter 1 Overview**


This chapter describes functions and basic operations of the instrument.

# 1.1 Description of Each Part


### 1.1.1 Rear



# 1.1.3 Left Side



# 1.1.4 Right Side



# 1.1.5 Accessories

Name	Function
AC adapter*	Used to power the instrument.
Electrode stand	Used to move and set electrodes during measurement.
Rubber cover	Protects the instrument side surfaces.
Instruction manual	Instructs the usage of the instrument.
Quick manual	Instructs the operations of calibration and measurement.

\*: The AC adapter includes 6 plug adapters. Referring to the following table, attach the appropriate plug adapter to the AC adapter depending on the country to be used.

① Australia	2 China	③ Europe	④ Korea	⑤ U.K., Singapore	⑥ USA, Canada, Taiwan
		(mm)	(mm)		

# **1.1.6 Identification of Manufacturing Date**

Manufacturing date can be identified from MFG No. described in the ID label on the backside of the instrument.

Third number from the left in the MFG No. indicates manufacturing year.


Forth alphabet from the left in the MFG No. indicates manufacturing month.

The alphabet is assigned to month according to the table below.

Ex.: ID: AA6A0000 means the device manufactured in 2016 January.

JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
А	В	С	D	E	F	G	Н	J	K	L	М

### 1.1.7 Operation Keys



Opera	ation key	Function
Ċ	POWER	Turns ON or OFF the power. (Press and hold for 2 seconds or more.)
	CAL	Displays the calibration screen (CAL screen).
MEAS	MEAS	Displays the measurement screen (MEAS screen).
	DATA	Displays the data screen (DATA screen).

NOTE

The POWER key does not work for 10 seconds after the AC adapter is connected. Wait for a while after connecting AC adapter.

# 1.1.8 Icons (Icon Bar)


The icons displayed on the bottom of the touch panel allow you to set up the instrument, check calibration information, and print out and save data.



	Icon	Function
₿	Menu	Used to perform measurement, display the Meter SET screen, and switch to the inspection and application modes.
Ģ	Information	Used to check calibration information on the MEAS or CAL screen, and application information on the Meter SET screen.
	User's guide	Used to check operation instructions and information about measurement and maintenance.
	Printer	Used to print out measurement or calibration values or saved data when a printer is connected.
€ <del>€</del>	Save in USB	Used to save measured data into a USB memory.
Ŷ	Save data	Used to save measurement values displayed on the screen into the instrument.
Ū	Trash box	Used to delete calibration data or the data saved in the instrument.
START	Operation	Used to start and stop the operations of measurement and calibration, and to change to the instantaneous value display. The icon label depends on the corresponding operation.

### 1.1.9 Status Icons


The icons displayed on the top of the touch panel show information on the instrument.



S	tatus icon	Function
A	Auto hold	Shows that the automatic hold function is ON, and that the end point is determined automatically according to input signals from the electrode based on the pre-selected stability criterion of measurement values. Refer to "2.2 Auto Hold Setting" (P.16).
M	Manual hold	Shows that the manual hold function is ON, and that the end point is determined manually. Refer to "2.2 Auto Hold Setting" (P.16).
•	USB1 <sup>*1</sup>	Shows that the instrument is connected with a personal computer via a USB cable.
•	USB2 <sup>*1</sup>	Shows that the instrument is connected with a USB data storage media.
	Printer	Shows that the instrument is connected with a printer with a dedicated printer cable.


\*1: These icons appear when a USB cable is connected, but it does not always mean that the communication is conducted.


# 1.1.10 Meas Screen



Indicator	Name	Description	
ļ	Temperature compensation electrode connection indicator	Displayed: Not displayed:	A temperature compensation electrode is connected. The displayed temperature is the electrode temperature (ATC). The displayed temperature is preset value (MTC).
Hold	HOLD indicator	Not displayed: Blinking: Lighting up:	An instantaneous value is displayed. In-process for HOLD HOLD completed.

The instrument has touch panel and keys and you can operate it by touching the screen. The 3 basic operations of Tap, Flick, and Drag allow you to operate the instrument intuitively. This section describes the basic operations.






The MEAS screen has three display methods to check variation and tendency of measurement values.


You can shift the screen to the digital, graph or analog screen by flicking it.

Digital display




If arrows, like and and appear when you touch the screen, you can flick in the directions to switch the screen display.

#### Graph display




Analog display



#### • Graph display

You can change the scale of the vertical axis in the graph display. It allows you to check a small change in measurement values.




Tap on the screen after the above operations, and the latest data will be displayed in optimized range.

#### Analog display


You can change the scale of the vertical axis in the analog display. It allows you to check a small change in measurement values.

The displayed range narrows. It allows you to check the detailed variation of measurement values.




wide-ranging variation of measurement values.

Tap on the screen after the above operations, and the latest data will be displayed in optimized range.




- **1.** Attach the stand shaft to the stand base.
- **2.** Attach the stopper and the stand arm to the stand shaft.

# 1.5.1 Electrode Connector



# 1.5.2 Temperature Connector




**1.** Insert the temperature connector into the jack socket on the instrument.

#### NOTE

If the temperature connector is unconnected or the connection is wrong, the MTC set temperature is displayed as the liquid temperature.

# **1.6 Connecting the Power Source**




1. Insert the AC adapter cable by fitting with the connector socket of in the instrument.

NOTE

Do not insert the cable with force when the connectors do not match.

AC adapter connector

# **1.7 Connecting the Printer**



Printer connector

The following printer is available.

#### Printer

CITIZEN CBM-910-24RJ120 V: plain paper type (Parts No.: 3014030146) CITIZEN CBM-910-24RJ230 V: plain paper type (Parts No.: 3014030147) Optional printer cable (Parts No.: 3014030148) is required.

NOTE


- Make sure to use the appropriate cable for the printer.
- Make sure to power OFF the instrument before connecting a printer.
- When you do not connect a printer with the instrument, disconnect the printer cable and put the rubber cap firmly on the connector socket on the instrument.

#### • Setting the Printer

Refer to the instruction manual of the printer for settings and operations of the printer.

- 1. Set the DIP switch No. 6 to ON and No. 7 to OFF, and then set printer paper and ink ribbon. Keep the LF key held down.
- Keep the SEL key held down.
   The printer prints output when the SEL key is being pressed.

# **1.8 Connecting the Personal Computer**




USB connector for personal computer communication

- Use designated cables to connect with a personal computer. Designated cable Parts name: USB cable (1 m) Parts No.: 3200373941
- Make sure that the transfer formats of the measuring instrument and personal computer are same. Otherwise, communication may fail due to a communication error or the online mode start failure. If you change the transfer formats, power OFF both of the instrument and the personal computer once, and then turn ON them again.
- For the details, register with our website and see the free download page of manuals.
- The communication software is subject to change without notice. Use the latest version of communication software uploaded on our website.

The latest version of software is 1.47 at this manual is issued.

Press and hold the POWER key for 2 seconds or longer.

Following the startup screen, the MEAS screen will be displayed.



#### NOTE


- The POWER key does not work for 10 seconds after the AC adapter is connected. Wait for a while after connecting AC adapter.
- If the following message appears on the screen during operation, disconnect the AC adapter and then connect it again and power ON the instrument.

==DS-7X series memory manager== Exception failure occurred. Please detach AC adapter and restart.

# Chapter 2 Before Measurement (Meter SET)

This chapter explains the procedures of the instrument condition setting, which should be performed before measurement.

# 2.1 Meter SET Screen



- Tap and tap Meter SET.
   Meter SET items are displayed.
   You will see the remaining items by dragging.
- **2.** Select items and set the conditions.

The setting procedures for each item are explained below.

# 2.2 Auto Hold Setting

set 🗛	↔2	2011/01/01 08:24
AUTO HOLD		NORMAL >
Sample name		>
Interval memory	OFF 💌	
USB Memory		>
Printer		>
Screen settings	1	>
<b>■ 0 </b> ■		

In the AUTO HOLD mode, the instrument judges potential stability automatically to the measurement values. This instrument allows you to select one among the 6 type criteria of potential stability.

- **1.** Change the auto hold settings, tap > on the right of the AUTO HOLD item.
- **2.** Tap  $\bigtriangledown$  on the right of the HOLD TYPE item.
- **3.** Select the measurement stability condition of the 6 types (EXACT, NORMAL, BRIEF, TIME, CUSTOMIZE, Manual) in the AUTO HOLD selection screen.

To cancel the operation, tap  $\mathbf{X}$  to return to the previous screen.

Each HOLD condition is described below.

Stability condition		Function						
Auto hold In the AUTO HOLD mode, the instrument judges potential state automatically to set the measurement values.								
	Measuring	Content						
Mode	target	Time (s)	Temperature (°C)	Criteria	【Default】			
	COND, Resis	st		Minimum display digit: 1 digit				
EXACT	SAL	10	2.0	0.30 ppt (0.03%)				
	TDS			10 mg/L				
	COND, Resis	st	2.0	Minimum display digit: 3 digits	Default setting of auto hold			
NORMAL	SAL	10		1.00 ppt (0.10%)				
	TDS			30 mg/L				
	COND, Resis	st	2.0	Minimum display digit: 5 digits				
BRIEF	SAL	10		3.00 ppt (0.30%)				
	TDS			100 mg/L				
TIME	Common	-	-	Arbitrarily set at 2 s to 999 s.	【10 s】			
	COND	Arbitrary	2.0	Arbitrarily set at 0.001 mS/cm to 0.100 mS/cm (0.1 mS/m to 10.00 mS/m).	【0.001 mS/cm (0.1 mS/m)】			
CUSTOMI ZE	SAL	setting 2 s to 60 s		Arbitrarily set at 0.10 PPT to 10.00 PPT (0.01% to 1.00%).	【0.3 PPT (0.03%)】			
	Resist	【10 s】		Setting value of COND is reflected.				
	TDS			Arbitrarily set at 0.1 mg/L to 100 mg/L.	【0.1 mg/L】			
Mai	nual hold	Determine a (Tap START		•				

#### 2.3 Custom Setting

set ∽	A 4	2011/10/20 09:58
<	AUTO HOLD	
Hold type		CUSTOMIZE 💌
Time		10 sec. 💌
рH		0.005 pH 💌
ION		0.3 mV 💌
mV		0.3 mV 💌
<b>B Ø N</b>		


We will explain the procedures of CUSTOMIZE setting taking the AUTO HOLD item as an example.

- **1.** Select the CUSTOMIZE of the Hold type to set the stability condition time and the stability condition value.
- 2. Use the numeric-key screen to enter measurement variations as HOLD criteria for each measurement item.

Tap < to return to the previous screen.

AUTO HOLD	2011/01/01 08:24 NORMAL >	Y	ou can set sample names for CH2.
Sample name	$\overline{\mathbf{b}}$		
Interval memory USB Memory Printer	OFF	1. 2.	Tap $>$ on the right of the Sample name item. Tap $\checkmark$ on the right of the item in the CH2 to enter the sample name.
Screen settings		3.	Tap [A1] to switch the keyboard entry screen of Alphabet> Numerical/Symbol. Tap [SHIFT] to input in lower-case alphabets. Up to 10 characters can be input.
		4.	Tap ENTER . The setting applies.
			To cancel the settings, tap $\left[ {f X}  ight]$ .
			Tap $\leq$ to return to the previous screen.
IINT ! To delete a register ng, and tap ENTER	•	tap	■ ■ on the right of the sample name, enter noth

\_\_\_\_



The measured data can be stored at set time intervals.

**1.** Tap **▼** on the right of the Interval memory item and select ON.

#### **Enter Interval Time**

set 🚓	2011/01/01 08:27
AUTO HOLD	NORMAL >
Sample name	$\triangleright$
Interval memory	ON 💌
Time	30 sec. 💌
USB Memory	$\geq$
Printer	$\triangleright$
<b>B 0 m</b>	

- Display the Time item when select ON.
   Tap ▼ on the right of the Time item.
- **2.** Enter the interval time in the numerical key screen.

(Setting range: 1 second to 999 seconds)

**3.** Tap ENTER . The setting applies.

To cancel the settings, tap  $[\mathbf{X}]$ .

#### 2.6 USB Memory Setting

՝Դ SET 🖪 ↔₂	2011/01/01 08:24
AUTO HOLD	NORMAL >
Sample name	$\triangleright$
Interval memory	OFF 💌
USB Memory	$\triangleright$
Printer	$\triangleright$
Screen settings	$\triangleright$
₿ 🖗 🛤	

Memory data can be written into a USB memory.

- **1.** Tap > on the right of the USB Memory item.
  - The USB memory setting screen is displayed.
    - Tap < to return to the previous screen.


#### **Simultaneously Memory**

set 🖪	↔₂ 20	11/01/01	08:30
< USE	8 Memory		
Simultaneously me	mory	OFF	
Eject			
Format			
Update			
<b>B Ø M</b>			

When a USB memory is inserted into the instrument, the data can be written into both the instrument and USB memories at the same time.

**1.** Tap ▼ on the right of the Simultaneously memory item and select ON.

#### Eject



Use this item to eject the USB memory from the instrument.

- Tap ▼ on the right of the Eject item and tap OK in the execution confirmation screen.
   To cancel the operation, tap CANCEL .
- **2.** When the ejection is completed, a notice message will appears. Tap OK.


#### NOTE

If you remove a USB memory from the instrument in a way other than mentioned above, data may not be saved correctly or data may be corrupted.

#### Format

★ SET         ▲         ←₂         2011/01/01 08:27           ✓         USB Memory         USB Memory           Simultaneously memory         OFF ▼	Use this item to format a USB memory in FAT16. Note that formatting deletes all stored data.
Eject  Format Update	<ol> <li>Tap ▼ on the right of the Format item and tap</li> <li>OK in the execution confirmation screen.</li> <li>To cancel the operation, tap CANCEL.</li> </ol>
<b>B () II</b>	A message that formatting is in progress appears during formatting. Do not remove the USB memory and do not disconnect the instrument power while this message appears. The instrument and USB memory are being accessed.
	<b>2.</b> When the format is completed, a notice message will appears. Tap OK.

#### 2.7 Printer Setting



The Printer item allows you to set printing contents, etc. effective only when a printer is connected with the instrument.


- Tap > on the right of the Printer item. The printer setting screen is displayed.
  - Tap < to return to the previous screen.

### When a printer is connected with the instrument, this item allows you to perform a printer test.

Tap 
 on the right of the Printer test item.
 The printing test is executed.
 Printout contents

!"#\$%&'()*+,/0123 456789:;<=>?@ABCDEFG HIJKLMNOPQRSTUVWXYZ[ ¥]^_`abcdefghijklmno
456789:;<=>?@ABCDEFG
HIJKLMNOPQRSTUVWXYZ[
¥]^_`abcdefghijklmno
pqrstuv#xyz{ }


#### Auto Printout



When a printer is connected with the instrument, this item allows you to perform an automatic printer test after measurement or calibration completion.

**1.** Tap **▼** on the right of the Auto printout item and select ON.


#### **Printout Layout**



This item allows you to change printing contents.

- **1.** Tap **▼** on the right of the Printout Layout item. The printing format screen is displayed.
  - Tap < to return to the previous screen.

#### When selecting CUSTOMIZE



CUSTOMIZE allows you to select items you want to print out among Measurement date/time, User name, Settings, Sensor info•Cal. User.

- 1. Select CUSTOMIZE from Printout Layout, and
  - tap  $\checkmark$  on the right of each printing item.
  - $\checkmark$  is ON: The item is selected.
  - is OFF: The item is not selected.

#### • Printout example

The following are the examples of BRIEF, NORMAL and GLP printouts.

Contents of results or conditions follow colon mark (:) of each item name.

If they exceed 10 characters, the exceeded part is displayed on the next line with right alignment.

When selecting CUSTOMIZE, you can select items that you want to print out among the GLP printing contents. But measurement values are always printed.

	RIEF easurement	-	RMAL easurement	•	ISTOMIZE) easurement	
Date Time COND HOLD Temperature	:2010/12/18 :11:36 : 1.001mS/cm :AUTO :25.0°C ATC	Date Time COND HOLD Temperature Operater Sample Inst. model Inst. SN Elect. model Elect. lot Temp. Coef Cal. Operator	∶1234567 ∶2. 00%/°C	GLPFormat Date Time COND HOLD Temperature Operater Sample Inst. model Inst. SN Elect. model Elect. lot	:2011/01/01 :10:10 :1.001mS/cm :AUTO :25.0°C ATC :«GUEST* :COND Sol :DS-72 :1234567 :3551-10D :1234567	Measurement date Measurement –value (Not be omitted) Measurement operator –Settings
				Cell Temp. Coef Cal. Operator Calibration o Cal.data Cal.Time		Electrode —Calibration operator —Calibration data

Signature

Signature

#### 2.8 Screen Settings

🍾 SET 🖪 ↔2	2011/01/01 08:24
AUTO HOLD	NORMAL >
Sample name	>
Interval memory	OFF
USB Memory	$\mathbf{>}$
Printer	$\mathbf{>}$
Screen settings	$\triangleright$
₿ Ø ₩	

You can change screen settings.

- **1.** Tap ≥ on the right of the Screen settings item. The screen settings screen is displayed.
  - Tap < to return to the previous screen.

#### **Screen Theme**

set 🕆	$\mathbf{A}$	2011/01/01 08:30
$\mathbf{<}$	Screen set	tings
Screen tl	neme	STANDARD
Display I	prightness	5 💌
Power sa	ving mode	OFF 💌

You can select one among 4 type screen themes; STANDARD, COOL, MONOTONE and KYOTO.

- **1.** Tap  $\bigtriangledown$  on the right of the Screen theme item.
- **2.** Select screen theme.

To cancel the operation, tap  $\mathbf{X}$  to return to the previous screen.

### ★ SET ▲ ←₂ 2011/01/01 08:30 ✓ Screen settings Screen theme STANDARD Display brightness 5 Power saving mode OFF


**Display Brightness** 

You can adjust the display brightness by tapping
[+] or $[-]$ , or by dragging on the scale.

- **1.** Tap **▼** on the right of the Display brightness item.
- 2. When the screen becomes the desired brightness, tap ENTER.

To cancel the operation, tap  $\mathbf{X}$  to return to the previous screen.

#### **Power Saving Mode**



You can set the time for power saving mode.

1. Tap ▼ on the right of the Power saving mode item and select ON.

When selecting ON for Power saving mode, the Back light off time item is displayed.

- **1.** Tap **▼** on the right of the Back light off time item.
- **2.** Enter the desired time on the numerical key screen. (Setting range: 1 minute to 999 minutes)
- Tap ENTER .
  The set time applies.
  To cancel the settings, tap X.

#### - HINT!

Display brightness

Power saving mode

Back light off time

During the power saving mode, the LED lamp above the POWER key lights up. Press the POWER key to exit the power saving mode.

5 💌

ON 💌

60 min 💌

#### 2.9 Sound Setting

set A 🔩	2011/01/01 08:31
Screen settings	
Sound settings	$\triangleright$
Language	English 💌
Security	$\triangleright$
Date/Time setting	
Analog output adj.	$\triangleright$
<b>B 0 M</b>	

You can change sound settings.

- **1.** Tap > on the right of the Sound settings item. The sound settings screen is displayed.
  - Tap < to return to the previous screen.

#### Sound Theme

set ∽	А	↔2	2011/01/01	08:32
<	Soun	ıd setti	ngs	
Sound them	ie		STANDARD1 💌	)
Sound volu	ime set	tings	5 💌	)
	-			

You can select one among 4 type sound themes; STANDARD, COOL, MONOTONE and KYOTO.

- **1.** Tap  $\bigtriangledown$  on the right of the Sound theme item.
- **2.** Select sound theme.

To cancel the operation, tap  $\mathbf{X}$  to return to the previous screen.

#### **Volume Setting**



You can adjust the sound volume by tapping + or -, or by dragging on the scale. When the sound volume is set to 0, the instrument is in the mute mode.

- **1.** Tap **▼** on the right of the Sound volume settings item.
- 2. When the screen becomes the desired volume, tap ENTER.

To cancel the operation, tap  $\mathbf{X}$  to return to the previous screen.

#### 2.10 Language Setting

set 🖪 🔶	2011/01/01 08:32
Language	English 💌
Security	$\mathbf{>}$
Date/Time setting	$\mathbf{>}$
Analog output adj.	$\mathbf{>}$
Temp. calibration	$\mathbf{>}$
Meter initialization	
80 ₩	

You can change language settings.

- **1.** Tap  $\bigtriangledown$  on the right of the Language item.
- 2. Select the language.To cancel the operation, tap X.

#### 2.11 Security Setting



Security setting allows you to set a password for an administrator of the instrument. After the setting is ON, the instrument requires you to select an operator name at the startup. Security setting, Date/Time setting, Analog output adj., Temp. calibration and Meter initialization are restricted to the administrator. To change the administrator or operator when the Security setting is ON, power OFF the instrument. At the next startup, the user selection screen appears to allow you change it. 25 administrators or operators in total can be registered.

**1.** Tap ≥ on the right of the Security item. The User management screen is displayed.

To cancel the operation, tap  $|\langle \rangle|$  to return to the previous screen.

**2.** Tap > on the right of the User management item and select ON.

#### When using the Security setting, administrator registration is required.

Administrator entry Require to resister administrator to activate user management function User name Password CANCEL ENTER	set 🕆	Α	↔2	2011/01/01	08:33
to activate user management function User name Password		Admini	strator	entry	
	to activate user management function User name				

- **1.** Tap the blank area at the right of "User name" to display the letter entry screen.
- **2.** Enter the operator name, and tap **ENTER**.

Tap the A1 to switch the keyboard entry screen of Alphabet and Numerical/Symbol. Tap the SHIFT to input in lower-case alphabets. Up to 12 characters can be input.

- **3.** Tap the blank area at the right of "Password" to display the numerical screen.
- **4.** Enter the password, and tap ENTER . The password can be set between 2 and 10 characters.
- **5.** Tap ENTER to set.

#### NOTE

When the Security setting is ON, at least 1 administrator is required for the instrument. Administrators have to keep the password. We recommend registering 2 or more administrators.

Administrator's names are marked with a star on the user selection screen.



When an operator is registered, the operator name can be put in measurement/calibration information, data printouts, data memory.

 Tap > on the right of the User entry/info change/delete item, when user registration, change password and user deletion. The User entry/info change/delete screen is displayed.

To cancel the operation, tap  $\leq$  to return to the previous screen.

#### **User Registration**

SET ▲ 😚 20	011/01/01 10:30
✓ User entry/info change,	/delete
User entry	>
User info change	>
User delete	>
Administrator entry	>
8 O H	
- <i>r</i>	

You can register operators.

- **1.** Tap > on the right of the User entry item.
- **2.** Tap the blank area at the right of "User name" to display the letter entry screen.
- **3.** Enter the operator name, and tap ENTER .

Tap the A1 to switch the keyboard entry screen of Alphabet and Numerical/Symbol. Tap the SHIFT to input in lower-case alphabets. Up to 12 characters can be input.

- **4.** Tap the blank area at the right of "Password" to display the numerical screen.
- **5.** Enter the password, and tap ENTER . The password can be set between 2 and 10 characters.
- 6. Tap ENTER to set.

#### NOTE

When the Security setting is ON, at least 1 administrator is required for the instrument. Administrators have to keep the password. We recommend registering 2 or more administrators.

Administrator's names are marked with a star on the user selection screen.

#### **User Information Changing**

Operators can change the password.

- **1.** Tap > on the right of the User info change item.
- 2. Enter the password, and tap ENTER .
- **3.** Tap the current password at the right of "Password" to display the numerical-key screen.
- **4.** Enter the password, and tap ENTER . The password can be set between 2 and 10 characters.

Only administrators can deregister an operator.

- **1.** Tap > on the right of the User delete item.
- **2.** Tap > on the right of the operator item.
- **3.** Tap OK.
  - Tap CANCEL, when do not deleting.

#### **Administrator Registration**

set 🕆	Α	↔2	2011/01/01 10:31
<	Admini	strator	entry
- LAQUA $\star$			~

Only administrators can assign/remove an operator as an administrator.

- **1.** Tap > on the right of the Administrator entry item.
- Tap to add a new administrator at the Administrator entry screen. Then, the lights up to show it is in the state of being selected. Tap to change the current administrator to operator.

At this time, the  $\checkmark$  lights out to show it is in the state of being unselected.

#### NOTE

When the Security setting is ON, at least 1 administrator is required for the instrument. Administrators have to keep the password. We recommend registering 2 or more administrators.

Administrator's names are marked with a star on the user selection screen.

**User Deleting** 

#### 2.13 Date Setting

Language English 💌	You can set the date and time.
Security Date/Time setting Analog output adj. Temp. calibration Meter initialization	<ul> <li>Tap &gt; on the right of the Date/Time setting item.</li> <li>The Date/Time setting screen is displayed.</li> <li>Tap &lt; to return to the previous screen.</li> </ul>
Date	
★ SET ▲ ←2 2011/01/01 08:34 ✓ Date/Time setting	You can set the date.

#### Time

set ≤	► ←₂ Date/Time set	2011/01/01 08:34 tting	You	can set
Year, mont	th, day		<b>1</b> та	ap (▼) c
hour/min			<b>2.</b> Ta	ap 🕂 C ap 🕂 C ap Entei
80.				o cancel evious :

et the time.

- on the right of the hour/min item.
- or 🗕 to set the time.
- R.

I the operation, tap  $\mathbf{X}$  to return to the screen.



Voltage output can be acquired from the analog output connector located at the instrument side

 Tap > on the right of the Analog output adj. item.

The Analog output adj. screen is displayed.

Tap < to return to the previous screen.

#### How to Analog Output Adj.

set 🕆	Α	$\Leftrightarrow_2$	2011/01/0	01 04:25
<	Analog	output	adj.	
CH2 Outpu	t 2V			
CH2 Outpu	t OV			•
CH2 Temp	2V			•
CH2 Temp	0V			•
,				

Connect the instrument with a digital multimeter, digital recorder, pen recorder or the like using a designated cable (analog output cable: Parts No.3014030152), and check and adjust the analog output value of the instrument.

- Tap ▼ on the right of the analog output item. The Output value adjustment screen is displayed.
- **2.** Tap ▲ or ▼ to adjust the analog output voltage.
- **3.** Tap ENTER .

To cancel the operation, tap  $\mathbf{X}$  to return to the previous screen.

#### 2.15 Temperature Sensor Calibration

in Set 🖪 🛩	2011/01/01 08:32
Language	English 💌
Security	>
Date/Time setting	$\mathbf{>}$
Analog output adj.	$\triangleright$
Temp. calibration	$\mathbf{>}$
Meter initializatio	n 💌
<b>B Ø M</b>	

You can perform calibration of the temperature sensor.

**1.** Tap > on the right of the Temp. calibration item.

The Temp. calibration setting screen is displayed.

To cancel the operation, tap < to return to the previous screen.

- 2. Display the measured temperature by the temperature sensor connected to the instrument. Display "-----", when not connecting the temperature sensor.
- **3.** Tap > on the right of the temperature sensor's channel item.
- **4.** Enter the temperature with the numerical screen and tap **ENTER**.
  - Tap  $|\mathbf{X}|$  when do not reflect the setting.

#### 2.16 Resetting to Factory Defaults



You can reset the instrument to the factory default conditions.

- **1.** Tap **▼** on the right of the Meter initialization item.
- **2.** Tap OK in the execution confirmation screen. Tap CANCEL, when do not resetting.
- **3.** Display the confirmation screen again, and tap OK.

Tap CANCEL, when do not resetting.

- **4.** Restart after the Meter initialization was finished. Press the POWER key to turn OFF.
- **5.** Press and hold the POWER key for 2 seconds to turn ON.

NOTE

If you disconnect the AC adapter after powering OFF, the POWER key does not work for 10 seconds after the AC adapter is reconnected. Wait for a while after reconnecting AC adapter.

#### Chapter 3 COND (Conductivity) Measurement


#### 3.1 COND Calibration

This section describes the procedures to set the conditions of COND calibration. Set the condition of temperature compensation before COND calibration according to "3.2 COND Measurement Setting" (P.39).

The cell constants of COND electrodes are different.

When using the conductivity electrode for the first time, set the cell constant written on the electrode into the instrument before use according to "3.2.1 Cell Constant Setting" (P.39).

#### 3.1.1 Automatic Calibration Setting



#### 3.1.2 Calibration of Standard Solution

A verified cell constant is written on a COND electrode label.

However, the actual cell constant may fluctuate depending on the usage circumstances and it is desirable to calibrate the cell constant in that case.

The procedures of cell constant calibration are mentioned below.

NOTE

- Perform "3.1.1 Automatic Calibration Setting" (P.35), before the following operations.
- •Make sure that the temperatures of the standard solution and the electrode are stable before the following operations. If you perform the operations with unstable temperatures, the calibration result may be incorrect.
- Immerse the electrode into the standard solution at the proper depth and stir it slowly with a stirrer. Do not return the used standard solutions into the original container. Dispose of them.
- Tapping **(**) on the COND CAL screen allows you to check the current calibration data. To clear the calibration data, tap **(**).
- •When auto calibration is turned OFF in "3.1.1 Automatic Calibration Setting" (P.35), it is necessary to enter the concentration value of the standard solution in accordance with the settings in " Manual calibration" (P.37).

#### Automatic calibration



**1.** Tap the channel and measurement item in the measurement screen to set "CH2" and "COND".

- 2. Tap START to start the calibration. When the calibration is completed, the HOLD indicator is lit up and the calibration state is displayed.
- **3.** Tap CLOSE after checking the calibration result to return to the CAL screen.

To start the COND measurement, press the MEAS key.

#### **Manual calibration**

CAL	Α	↔2	2015/10/30 15:37
CH2		K=1	24.98°C
			Set: 19.00 mS/cm
COND			9. 00 mS/cm
🔺 Press	START t	o start c	alibration
80			START

- **1.** Tap the value next to "Set:" to display the numerical-key screen.
- **2.** Enter the conductivity value of the standard solution, and then tap **ENTER**.
- Select the auxiliary unit of the standard solution for calibration using mS/m ▼, and then enter the concentration of the standard solution. After that, tap ENTER.
- **4.** Tap <u>START</u> to start the calibration. When the calibration is completed, the HOLD indicator is lit up and the calibration state is displayed.
- Tap CLOSE after checking the calibration result to return to the CAL screen.
   To start the COND measurement, press the MEAS key.

NOTE

The conductivity value of the standard solution used in the calibration process is the compensated value into the calibrating temperature by the temperature coefficient 2%/°C from the 25°C value. For more precise measurement, it is recommended to operate the calibration process at 25°C.

Temp.	Conductivity value at 25°C				
(°C)	84.00 (μS/cm)	1413 (μS/cm)	12.88 (mS/cm)	111.8 (mS/cm)	
0	64.01	776	7.15	65.4	
5	65.00	896	8.22	74.1	
10	67.00	1020	9.33	83.2	
15	68.00	1147	10.48	92.5	
16	70.00	1173	10.72	94.4	
17	71.00	1199	10.95	96.3	
18	73.00	1225	11.19	98.2	
19	74.00	1251	11.43	100.2	
20	76.00	1278	11.67	102.1	
21	78.00	1305	11.91	104.0	
22	79.00	1332	12.15	105.9	
23	81.00	1359	12.39	107.9	
24	82.00	1386	12.64	109.8	
25	84.00	1413	12.88	111.8	
26	86.00	1440	13.13	113.8	
27	87.00	1467	13.37	115.7	
28	89.00	1494	13.62	117.7	
29	90.00	1521	13.87	119.7	
30	92.00	1548	14.12	121.8	
31	94.00	1575	14.37	123.9	

#### • Conductivity standard values at various temperature

This section describes the procedures to set the conditions of COND measurement.

#### CH2 ► 2015/10/30 15:37 CH2 ► 25.0°C COND ► 0000 µS/cm Press START to start measurement. E • ■ ■ 말 START

- 1. Tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "COND".
- 2. Tap = and tap "CH2 MEAS SET".
- **3.** COND measurement setting items are displayed.

You will see the remaining items by dragging.

**4.** Select items and set the conditions.

The setting procedures for each item are explained below.


#### 3.2.1 Cell Constant Setting



- 1. Tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "COND".
- 2. Press the CAL key.
- 3. Tap 📰 and tap "CH2 CAL SET".

A screen to select Cell Constant and Auto Calibration is displayed.

**4.** Tap **▼** on the right of the Cell Constant item. The cell constant setting screen is displayed.



- **5.** Tap the left side numerical value to display the numerical screen.
- **6.** Enter the numerical value written on the COND electrode.
- 7. Tap ▲ or ▼ to enter the digit written on the COND electrode.
- 8. Tap ENTER.
  Reflect the setting.
  To cancel the settings, tap CANCEL.

#### NOTE

The unit indication of the cell constant depends on the electrode. Convert the unit to the one for the meter before input.  $100 \text{ m}^{-1} = 1 \text{ cm}^{-1}$  $1000 \text{ m}^{-1} = 10 \text{ cm}^{-1}$  $10 \text{ m}^{-1} = 0.1 \text{ m}^{-1}$ 


#### 3.2.2 COND Measurement Unit Setting

SET CH2 A 😽	2011/01/01 15:06
Unit	S/cm 💌
TEMP setting	S/m
TEMP conversion	Man S/cm FIX
Coefficient	2.00 %/°C 🔽
Reference	25 °C 💌
Alarm, upper limit	OFF

You can select S/m, S/cm or FIX (Unit is fixed at mS/cm as the COND measurement unit.

- **1.** Tap  $\bigcirc$  on the right of the Unit item.
- **2.** Select S/m, S/cm or FIX. The selected unit applies.

#### 3.2.3 Temperature Setting



There are two types of temperature setting for COND measurement; Automatic Temperature Compensation (ATC) and Manual Temperature Compensation (MTC). In ATC, the instrument detects the solution temperature with the connected temperature sensor, and performs temperature compensation for the COND values of the standard solutions used for calibration. In MTC, measure the solution temperature and enter the value in advance. The instrument performs temperature compensation using the entered temperature.

#### NOTE


If the temperature terminals of the instruction and electrode are not connected, temperature setting is performed in MTC even when ATC is set.

#### Solution Temperature Entry in MTC (Manual Temperature Compensation)

🖧 SET CH2 🖪 😽	2011/01/01 14:36
Unit	S/cm 💌
TEMP setting	MTC 💌
Temperature	25.0 °C 💌
TEMP conversion	OFF 💌
Alarm, upper limit	OFF 💌
Alarm, lower limit	OFF 💌

- Display the Temperature item when select MTC. Tap ▼ on the right of the Temperature item.
- **2.** Enter the solution temperature on the numericalkey screen.
- Tap ENTER . The setting applies. To cancel the settings, tap X.

#### 3.2.4 Temperature Conversion Function Setting



The measured COND value of a sample varies with the temperature. In addition, the change degree with temperature depends on the sample property.

If the change degree (temperature coefficient) of the sample is known, set this item to ON to display COND values converted at 25°C. If the temperature coefficient is unknown, set this item to OFF.

- **1.** Tap **▼** on the right of the TEMP conversion item.
- **2.** Select the temperature conversion method.

#### - HINT!

When select the pure water mode or the natural water mode, the temperature conversion conforms to the following standards. Pure water: ASTM D 1125-91 Table3

Natural water: ISO7888:1985 (JIS K0400-13-10:1999)

#### **Input Temperature Conversion Factor**

🔊 SET CH2 🖪 🔸	2015/10/30 15:39
Unit	S/cm 💌
TEMP setting	MTC 💌
Temperature	25.0 °C 💌
TEMP conversion	Manual 💌
Coefficient	2.00 %/°C 💌
Reference	25 °C 💌

- **1.** Tap **▼** on the right of the TEMP conversion item.
- 2. Select "Manual" on the TEMP conversion screen.
- **3.** Tap  $\overline{\bullet}$  on the right of the Coefficient item.
- **4.** Enter the temperature conversion factor on the numerical-key screen.
- 5. Tap ENTER . The setting applies.
  To cancel the settings, tap X.

#### **Reference Temperature Setting of Temperature Conversion**

🐼 SET CH2 🖪 🔸	2015/10/30 15:39
Unit	S/cm 💌
TEMP setting	MTC 💌
Temperature	25.0 °C 💌
TEMP conversion	Manua l 💌
Coefficient	2.00 %/°C 💌
Reference	25 °C 💌

- 1. Tap ▼ on the right of the TEMP conversion item.
- **2.** Select "Manual" on the TEMP conversion screen.
- **3.** Tap  $\bullet$  on the right of the Reference item.
- **4.** Enter any temperature value between 15°C and 30°C on the numerical-key screen.
- **5.** Tap ENTER. The setting applies.

To cancel the settings, tap  $[\mathbf{X}]$ .


#### 3.2.5 Alarm Setting

When the measurement values exceed the set upper or lower limit, the instrument detects it to display the notice on the screen or to output the signal from the external output terminal.

If the measurement values exceed the alarm range, the color of the pertinent channel "CH" is changes on the MEAS screen.

Set the upper limit alarm to ON for the upper limit control of measurement value. Set the lower limit alarm to ON for the lower limit control of measurement value.


Upper limit value



Lower limit value

SET CH2 🖪 🔶	2011/0	01/01 15:0
Reference	25	°C 💌
Alarm, upper limit	(	DFF 💌
Alarm, lower limit	(	DFF 💌
Electrode model	Custom	ON
Customize		OFF
Electrode lot		

#### Input Upper or Lower Limit Values Upper limit value entry



- When selecting ON the Alarm, upper limit item, the Upper limit value, tap ▼ on the right of the Upper limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.

To change the unit (mS/m,  $\mu$ S/m, etc.), tap on the unit change key on the right of the numerical-key screen.

**3.** Tap ENTER . The setting applies.

To cancel the settings, tap  $[\mathbf{X}]$ .

#### Input Upper or Lower Limit Values Lower limit value entry

- 🔊 SET CH2 🗛 🔸 2011/01/01 12:00 Alarm. upper limit ON 💌 Upper limit value 1.999 S/cm 🔻 ON 💌 Alarm, lower limit Lower limit value 0.003 µS/cm 🔽 Electrode model Customize 💌 Customize -田
- When selecting ON the Alarm, lower limit item, the Lower limit value, tap ▼ on the right of the Lower limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.

To change the unit (mS/m,  $\mu$ S/m, etc.), tap on the unit change key on the right of the numerical-key screen.

**3.** Tap ENTER . The setting applies.

To cancel the settings, tap  $\mathbf{X}$ .

NOTE

Even if changing units (S/m, S/cm, FIX), the alarm set value is not changed.

#### 3.2.6 Electrode Model Setting

When an electrode model is set, the model name can be displayed on data printouts or recorded in saved data.

Select the electrode model to be used for measurement.

You can set a desired name with up to 10 characters by selecting the Customize item.

## SET CH2 ▲ 2011/01/01 12:00 Alarm, upper limit ON Upper limit value 1.999 S/cm Alarm, lower limit ON Lower limit value 0.003 µS/cm Electrode model Customize

- **Electrode Model Selection**
- Tap ▼ on the right of the Electrode model item. The electrode model selection screen appears.

To cancel the settings, tap  $|\mathbf{X}|$ .

 Select the electrode model to be use. Tap an electrode model name, and the selected model applies.

#### **Electrode Model Entry**

🔊 SET CH2 🖪 🔸	2011/01/01 15:48
Electrode model	<b>X</b>
3573-10C	
3574-10C	
3582-10D	
6861-10D	
9382-10D	
Customize	

You can set a desired name with up to 10 characters.

- **1.** Tap "Customize" in the electrode model selection screen.
- When selecting Customize for the Electrode model item, the Customize item is displayed.
   Tap ▼ on the right of the Customize item.
- **3.** Enter an electrode model name using the keyboard screen.

Тар	<b>A</b> 1	to s	witch	the k	eyboa	rd e	entry	screen	of
Alpha	abet	>	Nume	erical/	′Symb	ol.	Тар	SHIFT	to
input	in lo	wer-	case	alpha	bets.				
		-							

Up to 10 characters can be input.


**4.** Tap ENTER . The setting applies.

To cancel the settings, tap  $|\mathbf{X}|$ .

```
— HINT! -
```

To delete a registered electrode model name, tap ▼ on the right of the electrode model name, enter nothing, and tap ENTER.

#### 3.2.7 Electrode Lot No. Setting



#### — HINT! -

To delete a registered electrode model name, tap ▼ on the right of the electrode model name, enter nothing, and tap ENTER.

This section describes the procedures of COND measurement.

# MEAS ▲ 2015/10/30 15:40 CH2 25.0°C COND 199.04 µS/cm Press START to start measurement. Image: Constant measurement. Image: Constant measurement. Image: Constant measurement.

- 1. Press the MEAS key, and tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "COND".
- 2. Tap START to start measurement. The measurement value is displayed, and the HOLD indicator blinks until the reading stabilizes.

To stop measurement tap **STOP** while the HOLD indicator blinks.

When the reading stabilizes, the value is held and HOLD indicator lights up.

During instantaneous value measurement, or when a measurement value is held, you can

store the measurement values by tapping **L** on the bottom of the screen.

**3.** After the measurement is completed, tap **STOP** to proceed to the next measurement.

#### Chapter 4 SAL (Salinity) Measurement

#### 4.1 Measurement Target Selection

Select measurement target for SAL measurement from sea water or other liquid.

CH2 SAL	2011/01/01 07:22 25.0℃	1.	Tap in t "S/
Press START to start measur	<b>).</b> 96 PPT		
▤◑◾ ⊎			
SET CH2 A +2 Unit SAL type Alarm, upper limit Alarm, lower limit	2015/10/30 15:42 PPT  NaCl OFF OFF OFF	2. 3.	Tar The dis Tar
8 11			
SAL type       SAL type       V       NaCL       SEA water	2015/10/30 15:42	4. 5.	tap

p the channel setting and measurement item the measurement screen to set "CH2" and AL".

- p Ħ and tap "CH2 MEAS SET". he SAL measurement setting items are splayed.
- $\mathbf{\nabla}$  on the right of the SAL type item.

SET (	CH2 🗛 🔶	2015/10/30 15:
SAL ty	/pe	<u> </u>
<b>~</b>	NaCl	
	SEA water	
	-	

- hen the target for measurement is sea water, p "SEA water."
- ıр 🗙 . ne setting applies.

This section describes the procedures to set the conditions of SAL calibration.

A SAL (salinity) value is obtained by conversion of a COND (conductivity) value. However, you can perform calibration using standard solutions.

Make sure to perform the calibration at the temperature specified on the standard solution label. The procedures are mentioned below.

NOTE
------

- •Before SAL calibration, do the unit settings of "3.1.1 Automatic Calibration Setting" (P.35) and "4.3 SAL Measurement Setting" (P.50).
- Tapping ① on the SAL CAL screen allows you to check the current calibration data. To clear the calibration data, tap 前.
- For the sea water measurement, perform the procedure in "4.1 Measurement Target Selection" (P.48) before SAL calibration.



CAL	l	λ	$\Leftrightarrow_2$	2011	/01/01 07:23
CH2					25.00°C
				Set:	0.96 PPT
SAL					00
				U	96
				•	PPT
					FFI
Press	START	to	start	calibrat	tion

- 1. Tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "SAL".
- **2.** Press the CAL key to display the SAL CAL screen.
- **3.** Wash the COND electrode with pure water (ion exchange water), and wipe it off by filter paper or tissue paper.
- **4.** Open the internal solution filler port of the COND electrode.
- **5.** Immerse the COND electrode into a beaker of the standard solution.
- **6.** Tap the numerical value at the right of "Set:" to display the numerical-key screen.
- **7.** Enter the salinity value of standard solution, and tap ENTER.

The conductivity value of standard solution used for calibration applies.

- **8.** Tap START to start the calibration. When the calibration is completed, the HOLD indicator is lit up and the calibration result is displayed.
- 9. Tap CLOSE after checking the calibration result to return to the CAL screen.
   To start SAL measurement, press the MEAS key.

This section describes the procedures to set the conditions of SAL measurement.

Salinity concentration is calculated (Practical Salinity Scale (UNESCO 1978)) from the measured value of conductivity.

Therefore, when the cell constant is set in conductivity measurement, there is no need to input the cell constant. If no cell constant is set, refer to "3.1.1 Automatic Calibration Setting" (P.35).




- 1. Tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "SAL".
- 2. Tap 📰 and tap "CH2 MEAS SET".

The SAL measurement setting items are displayed.

**3.** Select items and set the conditions.

The setting procedures for each item are explained below.

#### 4.4 SAL Measurement Unit Setting



- The measurement unit of SAL measurement values select either PPT or %.
- **1.** Tap **v** on the right of the Unit item.
- **2.** Select PPT or %. The selected unit applies.

#### 4.5 Temperature Setting

The settings of temperature compensation and temperature conversion in COND measurement apply for SAL measurement (refer to "3.2.3 Temperature Setting" (P.41) and "3.2.4 Temperature Conversion Function Setting" (P.42)).

#### 4.6 Alarm Setting


When the measurement values exceed the set upper or lower limit, the instrument detects it to display the notice on the screen or to output the signal from the external output terminal.

If the measurement values exceed the alarm range, the color of the pertinent channel "CH" is changes on the MEAS screen.

Set the upper limit alarm to ON for the upper limit control of measurement value.

Set the lower limit alarm to ON for the lower limit control of measurement value.

Upper limit value



Lower limit value


🖧 SET CH2 🖪 🔶	2015/10/30 15:42
Unit	PPT 💌
SAL type	NaCl 💌
Alarm, upper limit	ON 💌
Upper limit value	80.00 PPT 💌
Alarm, lower limit	OFF
	ON
•	OFF

#### 4.6.1 Input Upper or Lower Limit Values

#### Upper limit value entry

🔊 SET CH2 🖪 🔶	2015/10/30 15:43
Unit	PPT 💌
SAL type	NaCl 💌
Alarm, upper limit	ON 💌
Upper limit value	80.00 PPT 💌
Alarm, lower limit	OFF 💌

#### Lower limit value entry



- When selecting ON the Alarm, upper limit item, the Upper limit value, tap ▼ on the right of the Upper limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.
- Tap ENTER . The setting applies. To cancel the settings, tap X.
- When selecting ON the Alarm, lower limit item, the Lower limit value, tap ▼ on the right of the Lower limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.
- **3.** Tap ENTER . The setting applies.

To cancel the settings, tap  $\mathbf{X}$ .

The electrode model setting in COND measurement applies for SAL measurement (refer to "3.2.6 Electrode Model Setting" (P.45)).

#### 4.8 SAL Measurement

This section describes the procedures of SAL measurement.



- 1. Press the MEAS key, and tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "SAL".
- 2. Tap <u>START</u> to start measurement. The measurement value is displayed, and the HOLD indicator blinks until the reading stabilizes.

To stop calibration tap **STOP** while the HOLD indicator blinks.

When the reading stabilizes, the value is held and HOLD indicator lights up.

During instantaneous value measurement, or when a measurement value is held, you can store the measurement values by tapping de on the bottom of the screen.

**3.** After the measurement is completed, tap **STOP** to proceed to the next measurement.

# Chapter 5 Resist (Resistivity) Measurement

This section describes the procedures to set the conditions of Resist measurement.

# ΜΕΑS Δ 42 2011/01/01 14:38 CH2 25.0°C Resist 1.000 Ω・cm Press START to start measurement. Ξ I START

# 5.1 Resist Measurement Setting

- 1. Tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "Resist".
- 2. Tap I and tap "CH2 MEAS SET".

The Resist measurement setting items are displayed.

**3.** Select items and set the conditions.

The setting procedures for each item are explained below.

## 5.2 Resist Measurement Unit Setting

The measurement units ( $\Omega$ •m or  $\Omega$ •cm) of the Resist measurement values are reflecting the setting units (S/m or S/cm) of COND measurement setting ("3.2.2 COND Measurement Unit Setting" (P.40)).

# 5.3 Temperature Setting

The settings of temperature compensation and temperature conversion in COND measurement apply for Resist measurement (refer to "3.2.3 Temperature Setting" (P.41) and "3.2.4 Temperature Conversion Function Setting" (P.42)).

# 5.4 Alarm Setting


When the measurement values exceed the set upper or lower limit, the instrument detects it to display the notice on the screen or to output the signal from the external output terminal.


If the measurement values exceed the alarm range, the color of the pertinent channel "CH" is changes on the MEAS screen.

Set the upper limit alarm to ON for the upper limit control of measurement value.

Set the lower limit alarm to ON for the lower limit control of measurement value.

Upper limit value





### 5.4.1 Input Upper or Lower Limit Values

#### Upper limit value entry

🔊 SET CH2 🖪 😽	2011/01/01 14:39
Alarm, upper limit	ON 💌
Upper limit value	19.99 kΩ·cm 💌
Alarm, lower limit	ON 💌
Lower limit value	0.000 Ω·cm 💌

- When selecting ON the Alarm, upper limit item, the Upper limit value, tap ▼ on the right of the Upper limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.

To change the unit (M $\Omega$ •m, k $\Omega$ •m etc.), tap on the unit change key on the right of the numerical-key screen.

**3.** Tap ENTER . The setting applies.

To cancel the settings, tap  $\mathbf{X}$ .

### Lower limit value entry

🖧 SET CH2 🖪 😽	2011/01/01 14:39
Alarm, upper limit	ON 💌
Upper limit value	19.99 kΩ·cm 💌
Alarm, lower limit	ON 💌
Lower limit value	0.000 Ω·cm 💌

- When selecting ON the Alarm, lower limit item, the Lower limit value, tap ▼ on the right of the Lower limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.

To change the unit (M $\Omega$ •m, k $\Omega$ •m etc.), tap on the unit change key on the right of the numerical-key screen.

**3.** Tap ENTER . The setting applies. To cancel the settings, tap X.

# 5.5 Electrode Model Setting

The electrode model setting in COND measurement applies for Resist measurement (refer to "3.2.6 Electrode Model Setting" (P.45)).

## 5.6 Resist Measurement

A MEAS	↔2	2011/01/01	14:38
CH2		25	5.0°C
Resist ►			
	- 4	0	0
	_ 1	- U	U
	- 1	Ω	• cm
Press START to start	measure	ment.	
	Ŷ	ST	ART

- This section describes the procedures of Resist measurement.
  - 1. Press the MEAS key, and tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "Resist".
  - 2. Tap START to start measurement. The measurement value is displayed, and the HOLD indicator blinks until the reading stabilizes.

To stop calibration tap **STOP** while the HOLD indicator blinks.

When the reading stabilizes, the value is held and HOLD indicator lights up.


During instantaneous value measurement, or when a measurement value is held, you can

store the measurement values by tapping do n the bottom of the screen.

# Chapter 6 TDS (Total Dissolved Solids) Measurement


This section describes the procedures to set the conditions of TDS measurement.

## 6.1 TDS Measurement Setting



The setting procedures for each item are explained below.

## 6.2 TDS Measurement Mode Setting



ላጉ SET CH2 🖪 😚	2015/10/30 15:44
TDS type	LINEAR
TDS linear value	0.50 💌
Alarm, upper limit	OFF 💌
Alarm, lower limit	OFF 💌

### 6.2.1 Input TDS Linear Value when Selecting LINEAR

- **1.** Tap  $\checkmark$  on the right of the TDS linear value item.
- Enter the TDS linear value on the numerical-key screen and tap ENTER.
   The setting applies.

To cancel the settings, tap  $[\mathbf{X}]$ .

# 6.3 Temperature Setting

The settings of temperature compensation and temperature conversion in COND measurement apply for TDS measurement (refer to "3.2.3 Temperature Setting" (P.41) and "3.2.4 Temperature Conversion Function Setting" (P.42)).

# 6.4 Alarm Setting

When the measurement values exceed the set upper or lower limit, the instrument detects it to display the notice on the screen or to output the signal from the external output terminal.

If the measurement values exceed the alarm range, the color of the pertinent channel "CH" is changes on the MEAS screen.

Set the upper limit alarm to ON for the upper limit control of measurement value. Set the lower limit alarm to ON for the lower limit control of measurement value.

#### Upper limit value

🖧 SET CH2 🖪 😽	2015/11/11 22:33
TDS type	EN27888 💌
Alarm, upper limit	OFF 💌
Alarm, lower limit	ON
	OFF

Lower limit value

🔊 SET CH2 🖪 😽	2015/10/30 15:45
TDS type	EN27888
Alarm, upper limit	OFF
Alarm, lower limit	OFF
	ON
	OFF

### 6.4.1 Input Upper or Lower Limit Values

## Upper limit value entry


🐼 SET CH2 🖪 🔸	2015/10/30 15:46
TDS type	EN27888 💌
Alarm, upper limit	ON 💌
Upper limit value	1.00 g/L 💌
Alarm, lower limit	OFF 💌

- When selecting ON the Alarm, upper limit item, the Upper limit value, tap ▼ on the right of the Upper limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.

To change the unit, tap on the unit change key on the right of the numerical-key screen.

Tap ENTER .
The setting applies.
To cancel the settings, tap X.

## Lower limit value entry



- When selecting ON the Alarm, lower limit item, the Lower limit value, tap ▼ on the right of the Lower limit value item.
- **2.** Enter an upper limit value on the numerical-key screen.

To change the unit, tap on the unit change key on the right of the numerical-key screen.

Tap ENTER .
The setting applies.
To cancel the settings, tap X.

# 6.5 Electrode Model Setting

The electrode model setting in COND measurement applies for TDS measurement (refer to "3.2.6 Electrode Model Setting" (P.45)).

This section describes the procedures of TDS measurement.

# MEAS ▲ 2011/01/01 07:34 CH2 25.0°C TDS One Oogong/L Press START to start measurement. Image: Comparison of the start measurement. Image: Comparison of the start measurement. Image: Comparison of the start measurement.

- 1. Press the MEAS key, and tap the channel setting and the measurement item in the MEAS screen to set "CH2" and "TDS".
- 2. Tap START to start the measurement. The measurement value is displayed, and the HOLD indicator blinks until the reading stabilizes.

To stop calibration tap **STOP** while the HOLD indicator blinks.

When the reading stabilizes, the value is held and HOLD indicator lights up.

During instantaneous value measurement, or when a measurement value is held, you can

store the measurement values by tapping **L** on the bottom of the screen.

**3.** After the measurement is completed, tap **STOP** to proceed to the next measurement.

# **Chapter 7 Application Mode**

The application mode enables the measurement for the pharmaceutical water inspection methods under various Pharmacopeias by conductivity measurement in conformance to specific measurement methods. By simply submerging the electrode to a sample, the instrument will walk you through the process and will determine the result. This chapter explains about the settings and procedures of measurement using the pharmaceutical water inspection methods under various Pharmacopeias by conductivity measurement.

### 7.1 Pharmacopeia Mode

In this mode, evaluation of pharmaceutical water (purified water and injection syringe water) in conformity with US Pharmacopeia (USP), European Pharmacopeia (EP), Japanese Pharmacopeia (JP), and Pharmacopoeia of the People's Republic of China (PPRC) can be evaluated. This mode enables evaluation of pharmaceutical water that is measured based on the standards in accordance with the Pharmacopeia regulations in each country.


This mode has the function to indicate "Exceeded Limit" which shows that the sample does not conform to the specifications when a measurement value does not satisfy the Pharmacopeia regulations during measurement. When a measurement value is out of the specification after the measurement, the non-conformity is indicated in the measurement results. This applies for printouts.

This mode, you can save the measurement results only into a USB memory and print out them. If you need to save or print out the data, turn ON the "Simultaneously Memory" of "2.6 USB Memory Setting" (P.20) or "Auto Printout" of "2.7 Printer Setting" (P.22) in advance.

AS MEAS	Α	€2	2011/01	/01	07:48
COND EP				25	.0°C
Exceeded Limit		1	8		6
				mS	/cm
Press START	to start	t measure	ment.		
		€÷	(	ST	ART


AS MEAS		2011/01/	/01 07:48
E	P Measuremen	t result	
Meas. Value:	18.6mS/c	m	25.0 °C
Spec :	1.3 μS/	cm or less	25.0 °C
	EP Exceeded	limit	
Settings			
Cell const.:	1.000×	1 cm-1	
Temp. conv.:	0FF		
	•	- F	CLOSE
_		C	

#### 7.1.1 Shift to Pharmacopeia Mode



- **1.** Tap **H** and tap "Application".
- 2. Tap > on the right of the COND pharmacopeia mode item and select a desired Pharmacopoeia from USP, EP, JP, and CP (PPRC).

### 7.1.2 Measured by USP (Stage 1)



Evaluation is conducted based on the "7.1.10 Temperature and Conductivity Requirements" (P.69).

If the measured temperature is between the indicated temperatures, the value at temperature lower than the measured temperature is applied as the permissible conductivity.

- **1.** Select the USP in the COND pharmacopeia mode screen.
- **2.** Tap > on the right of the Stage 1 item.
- **3.** Before measurement, set the temperature conversion to OFF in accordance with the regulation prescribed by USP, and the setting of the unit is automatically changed to S/cm. The changed settings are applied.
- **4.** Tap OK to proceed to next the procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

#### 7.1.3 Measured by USP (Stage 2)

AS MEAS	А	€2	2011/01	/01 03:24
$\leq$	COND pha	irmacopoe	eia mode	
Stage 1				>
Stage 2				>
		•4		
_				

In this mode the value when the measured temperature is at 25°C  $\pm$ 1°C and the conductivity change for 5 minutes is 0.1  $\mu$ S/cm or less is judged whether it exceeds the evaluation standard, 2.1  $\mu$ S/cm or not.

- 1. Select the USP in the COND pharmacopeia mode screen.
- **2.** Tap > on the right of the Stage 2 item.
- **3.** Before measurement, set the temperature conversion to OFF in accordance with the regulation prescribed by USP, and the setting of the unit is automatically changed to S/cm. The changed settings are applied.
- **4.** Tap OK to proceed to the next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

#### 7.1.4 Measured by EP

AS MEAS	А	€2	2011/01/01 03:24
$\leq$	COND pha	irmacopoe	eia mode
USP			$\triangleright$
EP			$\triangleright$
JP			>
CP			>
		•¢	

Evaluation is conducted based on the "7.1.10 Temperature and Conductivity Requirements" (P.69).

If the measured temperature is between the indicated temperatures, the value at temperature lower than the measured temperature is applied as the permissible conductivity.

- **1.** Select the EP in the COND pharmacopeia mode screen.
- Before measurement, set the temperature conversion to OFF in accordance with the regulation prescribed by EP, and the setting of the unit is automatically changed to S/cm. The changed settings are applied.
- **3.** Tap OK to proceed to the next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.


#### 7.1.5 Measured by JP (OFF-LINE)

A MEAS	€2	2011/01/01	03:24
COND p	oharmacopoe	eia mode	
OFF-LINE		$\geq$	>
OmL - 10mL		$\mathbf{>}$	•
10mL-			•
	₩¢		

In this mode the value when the measured temperature is at  $25^{\circ}C \pm 1^{\circ}C$  and the conductivity change for 5 minutes is 0.1  $\mu$ S/cm or less is judged whether it exceeds the evaluation standard, 2.1  $\mu$ S/cm or not.

- **1.** Select the JP in the COND pharmacopeia mode screen.
- **2.** Tap > on the right of the OFF-LINE item.
- **3.** Before measurement, set the temperature conversion to OFF in accordance with the regulation prescribed by JP, and the setting of the unit is automatically changed to S/cm. The changed settings are applied.
- **4.** Tap OK to proceed to the next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

#### 7.1.6 Measured by JP (0mL-10mL (in container))



This is the test procedure for purified water, sterile purified water or water for injection contained in a container of 10 mL or less. The value when the measured temperature is at  $25^{\circ}C \pm 1^{\circ}C$  and the conductivity change for 5 minutes is 0.1  $\mu$ S/cm or less is judged whether it exceeds the evaluation standard, 25  $\mu$ S/cm or not.

- **1.** Select the JP in the COND pharmacopeia mode screen.
- **2.** Tap > on the right of the 0mL-10mL item. The changed settings are applied.
- **3.** Tap OK to proceed to the next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

#### 7.1.7 Measured by JP (10mL- (in container))

A MEAS	$\Leftrightarrow_2$	2011/01/01	03:24
COND pha	armacopoe	ia mode	
OFF-LINE			>
OmL-10mL			>
10mL-			>
	÷		

This is the test procedure for purified water, sterile purified water or water for injection contained in a container of 10 mL or more. The value when the measured temperature is at  $25^{\circ}C \pm 1^{\circ}C$  and the conductivity change for 5 minutes is 0.1  $\mu$ S/cm or less is judged whether it exceeds the evaluation standard, 5.0  $\mu$ S/cm or not.

- **1.** Select the JP in the COND pharmacopeia mode screen.
- **2.** Tap ≥ on the right of the 10mL- item. The changed settings are applied.
- **3.** Tap OK to proceed to the next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

## 7.1.8 Measured by PPRC (CP) (Stage 1)

A MEAS	↔2	2011/01/01 03:24
COND p	harmacopoe	ia mode
Stage 1		$\triangleright$
Stage 2		>
	· <del>(</del> +	

Evaluation is conducted based on the "7.1.10 Temperature and Conductivity Requirements" (P.69).

If the measured temperature is between the indicated temperatures, the value at temperature lower than the measured temperature is applied as the permissible conductivity.

- 1. Select the CP in the COND pharmacopeia mode screen.
- **2.** Tap > on the right of the Stage 1 item.
- **3.** Before measurement, set the temperature conversion to OFF in accordance with the regulation prescribed by PPRC, and the setting of the unit is automatically changed to S/cm. The changed settings are applied.
- **4.** Tap  $\bigcirc \mathsf{K} \upharpoonright$  to proceed to next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

### 7.1.9 Measured by PPRC (CP) (Stage 2)

AS MEAS	Α	€2	2011/01/01 03:24
<	COND pha	rmacopoe	eia mode
Stage 1			>
Stage 2			>
		*	

In this mode the value when the measured temperature is at  $25^{\circ}C \pm 1^{\circ}C$  and the conductivity change for 5 minutes is 0.1  $\mu$ S/cm or less is judged whether it exceeds the evaluation standard, 2.1  $\mu$ S/cm or not.

- **1.** Select the CP in the COND pharmacopeia mode screen.
- **2.** Tap > on the right of the Stage 2 item.
- **3.** Before measurement, set the temperature conversion to OFF in accordance with the regulation prescribed by PPRC, and the setting of the unit is automatically changed to S/cm. The changed settings are applied.
- **4.** Tap OK to proceed to the next procedure.
- Immerse the COND electrode in sample solution and tap START to start the measurement. When measurement is completed, the conductivity of the sample solution and the measurement condition are displayed as a measurement result.

## 7.1.10 Temperature and Conductivity Requirements

(for non-temperature compensated conductivity measurement)

Temperature (°C)	Required maximum (µS/cm)
0	0.6
5	0.8
10	0.9
15	1.0
20	1.1
25	1.3
30	1.4
35	1.5
40	1.7
45	1.8
50	1.9
55	2.1
60	2.2
65	2.4
70	2.5
75	2.7
80	2.7
85	2.7
90	2.7
95	2.9
100	3.1

Corresponding to USP (Stage1), EP, PPRC (CP) (Stage 1).

# **Chapter 8 Periodic Inspection Mode**

This chapter explains about the function to periodically check performance of the instrument and the electrode in COND measurements using.

We recommend that you perform the check once every 3 months. Setting conditions are described individually in each COND measurement item.

# 8.1 COND Periodic Inspection Mode Setting

There are two modes for the COND periodical check: Pharmacopoeia mode, or Checker (X-52) mode.

#### Pharmacopoeia mode

This mode conforms to the Japanese Pharmacopoeia 16th edition. You can check the cell constant and assess the conformity of the instrument.

#### Checker (X52) mode

Only the instrument check can be performed using the optional COND checker (X-52).

#### NOTE

The Pharmacopoeia mode is based on the corresponding regulations, but not fully compliant with the regulations. Note that the modes may not follow the regulations if the regulations are revised or amended.

#### 8.1.1 Pharmacopoeia Mode

You can perform the inspection compliant with the 16th edition of the Japanese Pharmacopoeia; checking the cell constant (within 5% difference between the actual cell constant and the value written on the COND electrode), measuring standard solutions 3 times to 5 times to check the error (within 5%) from the standard values and relative standard deviation (within 2%).

Before the operation, set the cell constant written on the COND electrode referring to "3.1.1 Automatic Calibration Setting" (P.35).

In this mode, the settings are changed as follows automatically.

Unit: S/cm

Temperature conversion: OFF Temperature setting: MTC, 20.0°C



TEST A	↔2	2011/01/01 03:26		
COND test	(Pharmaco	poeia) set		
Std. sol.	133.0	μS/cm		
Meas. number	3	time(s)		
Cell const.	1.065 * 0	0.1 cm-1		
Change a cell const. if this cell const. is different from it of COND EL.				
RESET	)	ENTER		

- **1.** Select COND periodic inspection in the check mode screen.
- **2.** Tap > on the right of the COND periodic inspection.
- **3.** Tap the Std. sol. value to display the numericalkey screen, and enter value of the standard solution used for the inspection.
- 4. Tap the Meas. number value, and use ▲ and
   ▼ to select measurement times (3 times to 5 times) for checking relative standard deviation.
- 5. After the setting is completed, tap ENTER.
   To return the set value to default, tap RESET.

According to the operation guide, perform the check. When the measurement and check is completed, the result data is displayed.

#### **Result data output**

- Measurement values
- Cell constant (calculated from the measured standard solution values)
- Error (difference between the cell constant written on the electrode and the calculated cell constant (regulated value: within 5%))
- Repeated measurement average
- Error (difference between the setting standard solution value and the repeated measurement average (regulated value: within 5%))
- Relative standard deviation (relative standard deviation at the repeated measurement (regulated value: within 2%))

NOTE


An accurate thermometer is required for the measurement. Prepare an accurate thermometer and perform the measurement at  $20^{\circ}C \pm 0.1^{\circ}C$ .

The cell constant calculated in this check does not apply for cell constant calibration.

#### 8.1.2 COND Checker (X-52) Mode

In this mode, the instrument operations are checked by using the optional checker X-52. Refer also to the instruction manual of the checker X-52 before the operation. When the COND periodical check mode (X-52) starts, the following items are set automatically as follows.

- Unit: S/m
- Cell constant:  $1.000 \times 100 \text{ m}^{-1}$
- Temperature setting: ATC



Follow the guidance to check.

#### Span check

19.00 S/m 1.900 S/m 190.0 mS/m 19.00 mS/m 1.900 mS/m Linearity check 10.00 S/m 1.000 S/m 100.0 mS/m 10.00 mS/m 1.000 mS/m 0.000 mS/m **Temperature check** 0.0°C 30.0°C 60.0°C 100.0°C

#### NOTE

The conductivity measurement values displayed during in the above operations are not concerned with measurement.

When all the check is completed, the result is displayed automatically.

#### Span check result

Criteria:  $\pm 0.5\% \pm 1$  digit of the full scale  $\pm 1.5\% \pm 1$  digit of the full scale only when 19.00 S/m is entered.

#### Linearity check result

Criteria:  $\pm 0.5\% \pm 1$  digit of the full scale  $\pm 1.5\% \pm 1$  digit of the full scale only when 10.00 S/m is entered.

#### Temperature check result

Indication error for each entry (regulated value ±0.4°C).

# 8.2 Comment Input

A comment can be entered up to 100 characters. Use this function to record periodical checks, etc.

Tap [INPUT] to use the function.

To delete the content input previously, tap [ALL DEL].



# Chapter 9 Data


The DATA screen allows you to check and delete saved measurement data, check the calibration data, save data into a USB memory, and delete all measurement and calibration data.

🖪 DATA 🖪 🔩 2	2011/01/01 04:25
Archive recorded of	data
Measured data_All	>
Measured data_latest50	>
Measured data_search	>
Copy all meas. Data	>
Delete all meas. Data	$\mathbf{>}$

You can search saved data by measurement item, operator, or sample name.

1. Press the DATA key to display the DATA screen.

# 9.1 Measured data\_All




**1.** Tap > on the right of the Measured data\_All item.

One item of measured data is displayed in one line. Data can be checked by dragging the item in order. 100 items of data can be viewed on 1 page.


- **2.** Tap **◀ ▶** to check other pages. The next 100 items are displayed.
- **3.** Tap > of the each data to check details of the data.
- **4.** Flick on a detailed data screen, and the previous/next detailed data screen is displayed.

# 9.2 Deleting Saved Data



You can check just the latest 50 data. The data are sorted in descending order of measurement data.

## 9.4 Measured data search



You can search saved data by one of measurement item, operator, or sample name. (You cannot use multiple search conditions at a time.)


- **1.** Tap |>| on the right of the Measured data\_search item.
- **2.** Search by measurement date Enter measurement date in the measurement date search screen, and tap ENTER . Search by measurement item Tap > on the right of each measured item. Search by operator name Enter operator name in input screen, and tap ENTER . Search by sample ID Enter sample name in input screen, and tap ENTER . **3.** When you select Measured on, enter the measurement date and tap **ENTER** on the next
- screen.

When you select Measuring mode, tap > on a measurement item on the next screen.

When you select User name, enter operator name and tap **ENTER** on the next screen.

When you select Sample ID, enter sample name and tap **ENTER** item on the next screen.

Search is performed and the result is displayed.



You can save the copy of the measurement data saved in the instrument into a USB memory. To execute the copy, connect a USB memory to the instrument.

- **1.** Tap > on the right of the Copy all meas. Data.
- **2.** Tap YES to copy the all measurement data. To cancel the operation, tap NO.
- **3.** Tap OK in the Copy all meas. Data completion screen.

NOTE

Before copying data, make sure that sufficient capacity is available in the USB memory. If the copy stops in the middle, turn OFF the power and reboot the instrument, and then execute the copy again.

#### 9.6 Delete all meas. Data



You can delete all measurement data saved in the instrument.

- **1.** Tap > on the right of the Delete all meas. Data
- **2.** Tap YES to delete the all measurement data. To cancel the operation, tap NO.
- **3.** Tap OK in the Delete all meas. DATA screen.

# **10.1 Model Information**

Item	Description
Brand (pet name)	LAQUA
Series name	Benchtop pH/Water Quality Analyzer
Model	DS-72G
Model description	COND METER

# **10.2 Measuring Object**

Measuring object	Item	Description
	Measuring principle	Thermistor method
	Display range	-30.0°C to 130.0°C
Temperature	Measuring range	0.0°C to 100.0°C
	Resolution	0.1°C
	Repeatability	±0.1°C ±1 digit
	Measuring principle	2 AC bipolar method
		Cell constant 1 cm <sup>-1</sup> : 0.00 µS/cm to 199.9 mS/cm
Conductivity	Measuring range	Cell constant 0.1 cm <sup>-1</sup> : 0.000 $\mu$ S/cm to 19.99 mS/cm
(COND)		Cell constant 10 cm <sup>-1</sup> : 0.0 µS/cm to 1.999 mS/cm
	Resolution	0.05% of full scale
	Repeatability	±0.5% ±1 digit of full scale
	Measuring principle	Conversion from conductivity value
		Cell constant 1 cm <sup>-1</sup> : 0.000 k $\Omega$ •cm to 19.99 M $\Omega$ •cm
Resistivity	Measuring range	Cell constant 0.1 cm <sup>-1</sup> : 0.00 k $\Omega$ •cm to 199.9 M $\Omega$ •cm
(Resist)		Cell constant 10 cm <sup>-1</sup> : 0.0 $\Omega$ •cm to 1.999 M $\Omega$ •cm
	Resolution	0.05% of full scale
	Repeatability	±0.5% ±1 digit of full scale
Colligit	Measuring principle	Conversion from conductivity value
Salinity (SAL)	Measuring range	0.00 PPT to 80.00 PPT (0.000% to 8.000%)
	Resolution	0.01 PPT (0.001%)

Measuring object	Item	Description
	Measuring principle	Conversion from conductivity value
TDS	Measuring range	0.01 mg/L to 1000 g/L
	Resolution	0.01 mg/L

# 10.3.1 Meter Default Settings

Item		Selection item/Setting range	Default values
Security Security management function		Enable, Disable	Disable
Hold condition	Hold setting mode	EXACT, NORMAL, BRIEF, TIME, CUSTOM, OFF (Manual)	NORMAL
In selecting "TIME"	Time setting value	2 seconds to 999 seconds	10 seconds
	Time setting value	2 seconds to 60 seconds	10 seconds
In selecting	Conductivity variation width	1 digit to 100 digits	1 digit
"CUSTOM"	Salinity variation width	0.10 PPT to 10.00 PPT	0.30 PPT
	Resistivity variation width	1 digit to 100 digits	1 digit
	TDS variation width	0.1 mg/L to 100.0 mg/L	100.0 mg/L
Interval	Interval memory function	Enable, Disable	Disable
memory	Time setting value	1 second to 999 seconds	30 seconds

## 10.3.2 Measurement Condition Default Settings (Can Be Set per Operator)

Item			Selection item/ Setting range	Default values
	Alarm	Upper limit value setting	Enable, Disable	Disable
		Lower limit value setting	Enable, Disable	Disable
	condition	Upper limit value	0.003 μS/cm to 1.999 S/cm	1.999 S/cm
		Lower limit value	0.003 μS/cm to 1.999 S/cm	0.003 μS/cm
	Measurement	value unit	S/m, S/cm, FIX	S/cm
Conductivity	Temperature	Temperature setting	ATC (Automatic temperature compensation), MTC (Manual temperature compensation)	ATC
measurement condition	Setting	Temperature input value in selecting "MTC"	0.0°C to 100.0°C	25.0°C
	Temperature conversion	Temperature conversion function	Pure water, Natural water, Manual, Disable	Manual
		Temperature conversion coefficient	0.00%/°C to 10.00%/°C	2.00%/°C
		Reference temperature	15.0°C to 30.0°C	25.0°C
	Electrode data	Model		None
		lot No.		None
	Alarm condition	Upper limit value setting	Enable, Disable	Disable
Salinity measurement condition		Lower limit value setting	Enable, Disable	Disable
		Upper limit value	0.00 PPT to 80.00 PPT (0.000% to 8.000%)	80.00 PPT
		Lower limit value	0.00 PPT to 80.00 PPT (0.000% to 8.000%)	0.00 PPT
	Measurement value unit		PPT, %	PPT
	SAL type		NaCl, SEA water	NaCl

	Item		Selection item/ Setting range	Default values
Resistivity measurement	Alarm condition	Upper limit value setting	Enable, Disable	Disable
		Lower limit value setting	Enable, Disable	Disable
condition		Upper limit value	0.000 Ω•cm to 19.99 kΩ•cm	19.99 kΩ•cm
		Lower limit value	0.000 Ω•cm to 19.99 kΩ•cm	0.000 Ω•cm
		Upper limit value setting	Enable, Disable	Disable
	Alarm	Lower limit value setting	Enable, Disable	Disable
TDS measurement condition	condition	Upper limit value	0.00 mg/L to 1000.0 g/L	1.00 g/L
Condition		Lower limit value	0.00 mg/L to 1000.0 g/L	0.00 mg/L
	TDS measurement mode		LINEAR, 422, EN27888, NaCl	LINEAR
	TDS linear value		0.40 to 1.00	0.50
Sample ID				None
	Language		Japanese, English, Chinese, Korean, Vietnamese	English
	Screen setting	Screen theme	STANDARD, COOL, MONOTONE, KYOTO	STANDARD
		Brightness	1 to 10	5
		Power saving mode	Enable, Disable	Disable
Interface		Back light off time	1 minute to 999 minutes	60 minutes
condition		Volume	0 to 9	5
	Sound setting	Sound theme	STANDARD1, STANDARD2, AQUA, KYOTO	STANDARD1
	Printer setting	Automatic printing	Enable, Disable	Disable
		Printing format	BRIEF, NORMAL, GLP, CUSTOMIZE	NORMAL
	USB memory	Simultaneous memory	Enable, Disable	Disable

# 10.4 Options

This section lists spare and optional parts for LAQUA series. These parts are possible through our representatives in your region. Place an order specifying their name, model, and part number.

Part name		Part number	Remarks
AC adapter		3200647413	With 6 plug adapters
Plain paper printer	Printer (USA, 120 V)	3014030146	Printer cable sold separately
	Printer (EU, 230 V)	3014030147	
	Printer cable	3014030148	1.5 m
	Roll paper	3014030149	20 rolls/set
	Ink ribbon	3014030150	5 pcs/set
USB cable		3200373941	1 m
Serial cable		3014030151	
Analog (alarm) output cable		3014030152	
Electrode stand (Standard type)		3200382557	
Electrode stand (Long type)		3200382560	
Stand arm		3200373991	
Sensor holder		3200373961	
X-51 Digital Simulator		_	For pH, mV, ION, and DO
X-52 Digital Simulator		-	For COND



http://www.horiba-adt.jp

For any questions regarding this product, please contact your local agency, or inquire from the following website. http://global.horiba.com/contact\_e/index.htm

