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CODE: I042661000

Metal (ORP) Electrode Operation Manual

Thank you for purchasing a HORIBA ORP electrode.
Read this Operation Manual carefully before using the electrode.
© Copyright HORIBA, Ltd. 1997

参 考
Reference Only

1 Cautions during handling

- The D-20 Series of pH meters can be used in combination with this electrode to provide water-resistant construction (conforming to IP-67). When measuring, do not immerse the cap, lead or connector in the sample.
- Do not allow the electrode to come in contact with any hard surface.
- Always use #300 internal solution with the reference electrode.
- If the liquid junction of the electrode is dry, immerse the tip of the electrode in pure (ion exchange) water prior to use.
- The electrode connector requires a high degree of insulation. Do not allow the connector to come in contact with water or dirty hands.

2 Contents

- Electrode 1 pc.
- Operation Manual 1 copy

3 Specifications and names of parts

• Model: 9300-10D

Specifications	Measurement range	0 - ± 1,999 mV
	Usable temp. range	0 - 60 °C
	Storage temp. range	0 - 50 °C
	Metal electrode	Pt
	Internal electrode	Silver/silver chloride
	Reference electrode internal solution	3.33 mol/L - KCl solution (#300)
	Liquid junction material	Porous ceramic
	Wetted material	Glass, ceramic
	Dimensions:	
	Electrode length	150 mm (incl. cap)
	Height of wetted part	11 mm (position from electrode tip)
	Ext. diam. of wetted part	12 mm
	Lead length	1 m
Names of parts		

4 Preparations

4.1 Connecting to pH meter

1 Insert the electrode connector into the connector port sleeve on the meter, after aligning with the pin. Do not insert the connector unless it is aligned properly with the connector port.

2 Press the electrode connector into the connector port on the meter, while turning the connector to the right.

3 Slide the connector cover over the connector. Then, screw the cover to the right until it comes in light contact with the meter case.

4 Insert the temperature connector into the jack on the meter. Insert the connector firmly, until the O-ring on the connector can no longer be seen.

4.2 Preparing electrode

Caution

Chemical solution

Highly concentrated potassium chloride (3.33 mol/L KCl) is used in the internal solution of the electrode. If the internal solution in the electrode comes in contact with your hands or skin, wash immediately with water. If the internal solution comes in contact with your eyes, flush immediately with large amounts of water and seek treatment by a physician.

Caution

Glass fragments

Glass fragments cause injury. The outer tube of the electrode is made from glass. Use care not to break it.

1 Hold the cap and remove the protective cap.

2 Fill the electrode with new internal solution until it nears the rubber stopper.

3 Close the opening with the rubber stopper, to apply pressure to the internal solution and confirm that the internal solution is seeping from the liquid junction.

4 Wash the electrode tip well with pure (ion exchange) water, then wipe it dry with filter or tissue paper.

Note ▶▶▶ Replace the electrode internal solution after using or storing the electrode for a long period of time (approx. two or three months).

- Use a plunger to remove the electrode internal solution.

- Fill the electrode with new internal solution until it nears the rubber stopper in the filler port.

5 Points to remember when measuring (calibrating)

1 Wash the electrode in pure (ion exchange) water and dry it with filter or tissue paper, each time prior to immersing it in the standard solution or sample.

2 Leave open the rubber stopper to the internal solution filler port of the electrode.

3 Immerse the electrode in the sample, to at least three centimeters from the electrode tip. Make sure that the surface of the internal solution in the electrode is higher than the surface of the sample.

4 After immersing the electrode in the sample, lightly stir the electrode around to both get it used to the sample and remove any air bubbles.

Ref. ▶▶▶ Electrode calibration

- The measurement potential with ORP standard solution is used as a guide for knowing the dirtiness, etc. of the metal electrode being used. The measurement potential cannot be used like pH standard solution to calibrate the meter.
- The operation check performed using ORP standard solution is used to diagnose the electrode or meter, according to whether or not readings are within ± 15 mV of the values shown in the below table.

ORP standard solution readings (mV) at various temperatures

Standard solution (°C)	160-22 Phthalate + quinhydrone	160-51 Neutral phosphate + quinhydrone
0	+277.5	+116.5
5	+274.2	+111.9
10	+270.9	+106.9
15	+266.8	+101.0
20	+262.5	+95.0
25	+257.6	+89.0
30	+258.5	+82.7
35	+248.6	+76.2
40	+243.6	+69.0

Ref. ▶▶▶ Problems with sample being measured

- Note that general responsiveness and reproducibility may be poor when measuring the ORP of samples with extremely low concentrations of oxides or reductants, such as tap water, well water, and pure water.
- When alkali ion water is allowed to sit, it absorbs Carbon Dioxide in the air and its ORP undergoes large changes. For this reason, the ORP of alkali ion water should be measured quickly.

6 Maintenance

- If the electrode is very dirty and cannot be washed clean using pure (ion exchange) water, wash it using the appropriate method below. Then, rinse the electrode well using pure (ion exchange) water.
- To maintain electrode performance, replace the internal solution approximately every two or three months. Refer to the section entitled, "Preparing electrode."

- **General/oily dirt**
Wipe the electrode using a cotton gauze containing neutral cleansing agent.

- **Inorganic or other dirt**
Immerse the electrode in diluted nitric acid (a solution of nitric acid diluted by a factor of two).

Note ▶▶▶ If brown colored silver chloride is detected in the internal solution of the reference electrode, replace the internal solution (#300) to prevent clogging of the liquid junction.

7 Storage

1 Wash the electrode well using pure (ion exchange) water, to remove any sample still clinging to the electrode.

2 Put the rubber stopper into the internal solution filler port.

3 Wash the inside of the protective cap with pure (ion exchange) water, then add enough pure (ion exchange) water to soak the sponge.

4 Attach the protective cap.

Note ▶▶▶ Storing the electrode for an extended period of time while the inside of the protective cap is dry may lead to a decline in electrode responsiveness and sensitivity.

Note ▶▶▶ Avoid storing the electrode in hot place or places with high humidity. Store the electrode indoors, out of direct sunlight.

Technical questions regarding this product should be directed to the following HORIBA Customer Support Center.

Call toll free (in Japan): 0120 - 37 - 6045
Hours: 9:00 a.m. - 5:00 p.m., Mon. - Fri.

HORIBA, Ltd.