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Electrode Cleaning Solution Instruction Manual

230 (for electrode for low conductivity solution)

CAUTION	
	<p>Chemical solution</p> <p>When you handle the cleaning solution, make sure to use latex gloves, protective glasses, and a protective mask. Cleaning solution must not come in contact with the skin, or consumed. This solution stimulates eyes, the skins and the lungs. If it comes in contact with your skin or eyes, flush immediately with running water. If you drink this solution by accident, your mouth, esophagus, and gastrointestinal will be damaged. If swallowed, seek medical advice.</p>

● Specification

Model	230	Usable temperature range	5°C to 30°C
Name	Electrode cleaning solution	Expiration date	Described on the label.
Cleaning solution A composition	0.1 mol/L ammonium hydrogen fluoride Deionized water	Storage	In a cool and dark place (5°C to 30°C), Storage is impossible after opening.
Cleaning solution B composition	0.01 mol/L hydrochloric acid Deionized water	Content	Cleaning solution A: 30 mL (1 bottle) Cleaning solution B: 100 mL (1 bottle)
Packaged contents	Cleaning solution A, cleaning solution B, measuring cup		


● Application

Use this cleaning solution for measurements of the low conductivity or the low buffer capacity solution at the following cases.

- For conditioning the electrode before the measurement
- When the response of the electrode becomes slow

● Usage

1. Put 10 mL or more of the cleaning solution A into the measuring cup.
2. Open the internal solution filler port of the electrode.
3. Immerse the electrode in the cleaning solution A, and keep it immersed for three minutes (refer to the figure on the right).
4. Wash the electrode with pure water (or deionized water).
Dispose the waste solution using a bottle made of plastic.
5. Immerse the electrode in the cleaning solution B, and keep it immersed for 12 hours to 24 hours.
6. Wash the electrode with pure water (or deionized water).
7. Remove the protective guard from the tip of the electrode, and then wipe the responsive glass membrane with filter paper moistened with pure water (or deionized water).
8. Be sure to perform calibration prior to measurement.

● **Storage**

- Store the cleaning solution in a cool and dark place.
- Use the cleaning solution immediately after opening the bottle.

● **Caution**

- Be sure to use the accompanying measuring cup when using the cleaning solution A. Do not use the glass container.
- Do not immerse the electrode in the cleaning solution A over 90 minutes.
- If the performance of electrode does not recover even after cleaning, it may be the end of electrode life time.
- When disposing of the cleaning solution, follow the related laws and/or regulations of your country for disposal of the product.