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## Electrode Cleaner Instruction Manual 250 (for cleaning protein dirt)

| <b>CAUTION</b>   |
|--|
| <p><b>!</b> Chemical solution</p> <p>Cleaner must not come in contact with the skin, or consumed. This solution stimulates eyes and the skins. If it comes in contact with your skin or eyes, flush immediately with running water. If you drink this solution by accident, your mouth, esophagus, and gastrointestinal will be damaged. If swallowed, seek medical advice. When you handle the cleaner, make sure to use latex gloves, protective glasses, and a protective mask.</p> |

### ● Specification

|             |  |                          |  |
|-------------|--|--------------------------|--|
| Model       | 250  | Usable temperature range | 5°C to 30°C  |
| Name        | Electrode cleaner  | Expiration date          | Described on the container   |
| Composition | Sodium chloride: 0.5% to 1.0%、<br>Tris (hydroxymethyl) aminomethane:<br>0.1% to 0.5%<br>Surface active agent: 0.1% to 0.5%<br>Enzyme: 0.1% to 0.5%<br>Sodium azide: < 0.1% | Storage                  | In a dark place (5°C to 30°C),<br>Storage is impossible after opening. |
| Quality     | pH 8 to pH 10  | Content                  | 400 mL (1 bottle)  |

### ● Application

Use this cleaner when response and sensitivity declines due to the following conditions.

- Clog by the solution containing inorganic component on the liquid junction
- Dirt on the responsive membrane

### ● Usage

1. **Open the internal solution filler port of the electrode.**
2. **Immerse the electrode in this cleaner up to 3 cm or more from the electrode tip, and keep it immersed for an hour or more.**  
If clog or dirt is not removed after an hour, continue immersing it in the cleaner half to all day.
3. **Rinse the electrode with deionized water, and replace the internal solution.**
4. **Immerse the electrode in deionized water, and keep it immersed for an hour or more.**
5. **Perform calibration again before measurement.**

### ● Storage

- Store the cleaner in a dark place.
- Use the cleaner immediately after opening the bottle.

### ● Caution

- If the performance of electrode does not recover even after cleaning, it may be the end of electrode life.
- This cleaner is ineffective on inorganic component dirt.
- This cleaner contains sodium azide, which is acutely toxic. But the cleaner is not designated as a poisonous substance, because the content is 0.1% or less.
- When disposing of the cleaner, follow the related laws and/or regulations of your country for disposal of the product.