

Aquaculture		Companion Animals			
Shrimp <ul style="list-style-type: none"> AHPND / EMS Plasmid AHPND / EMS Toxin 1 Baculovirus penaei CMNV EHP IHHNV IMNV MBV MNV NHPB PvNV SHIV TSV Vibrio fluvialis Vibrio harveyi WSSV YHV 	Fish <ul style="list-style-type: none"> Aeromonas salmonicida Carp Edema Virus (CEV) Flavobacterium psychrophilum IRIDO-Megalocytivirus IRIDO-Ranaviruses ISAV KHV Streptococcus agalactiae Streptococcus iniae SVCV Tilapia Lake Virus (TiLV) VHSV Vibrio alginolyticus VNN 	Canine <ul style="list-style-type: none"> Anaplasma platys Babesia gibsoni Bordetella bronchiseptica Canine Adeno Virus 2 Canine Babesia Canine Influenza Virus Canine Leishmaniasis CDV CHV CPIV CRCoV D. immitis 	<ul style="list-style-type: none"> E. granulosus / E. multilocularis Echinococcus canadensis Ehrlichia canis Influenza H3N8 Leptospirosis lipL32 Lyme disease Mec A gene for MRSA Neospora caninum Parvovirus Parvovirus (wild type) 	Feline <ul style="list-style-type: none"> Bartonella henselae Bordetella bronchiseptica Candidatus Mycoplasma haemominutum Candidatus Mycoplasma turicensis Chlamydia felis Cytauxzoon felis D. immitis FCoV FCV Feline Parvovirus FeLV FHV 	<ul style="list-style-type: none"> FIV Leptospirosis lipL32 Lyme Disease Mec A gene for MRSA Mycoplasma felis Mycoplasma haemofelis Parvovirus (wild type) Toxoplasma gondii
Equine					
<ul style="list-style-type: none"> Anaplasma phagocytophilum EAV EHV-1 EHV-3 EHV-4 EIAV Influenza H3N8 Leptospirosis lipL32 Pseudomonas aeruginosa Rotavirus A Salmonella spp. Streptococcus equi Taylorella equigenitalis (CEM) Theileria equi 					
Livestock					
Bovine <ul style="list-style-type: none"> Bovine Leukemia Virus Bovine Papular Stomatitis Virus Bovine Tuberculosis Brucella abortus Brucella spp. BVDV-1 BVDV-2 Campylobacter fetus 	<ul style="list-style-type: none"> Infectious Bovine Rhinotracheitis Lumpy Skin Disease Virus Q Fever Rotavirus Tritrichomonas foetus 	Poultry <ul style="list-style-type: none"> ALV-J Avian Influenza H9 Avian Metapneumovirus Avian Reovirus (Asia & America regions only) Chicken Infectious Anemia Virus Duck Virus / Enteritis Virus Fowl Adenovirus (Pan Adv, type 4, type 8b) 	<ul style="list-style-type: none"> Infectious Bronchitis Virus Infectious Bursal Disease Virus Infectious Laryngotracheitis Influenza A Influenza H5 (Asian Lineage) Influenza H5 (American Lineage) Influenza H7 	<ul style="list-style-type: none"> Influenza HP-H7 with insertion Influenza N9 Marek's Disease Marek's Disease serotype 3 (HVT) Mycoplasma gallisepticum Mycoplasma synoviae NDV-class 1 NDV-class 2 	<ul style="list-style-type: none"> NDV-LaSota REV Salmonella spp.
Swine <ul style="list-style-type: none"> Actinobacillus Pleuropneumonia Africa Swine Fever Virus Brachyspira hyodysenteriae Chlamydia psittaci CSFV FMDV JEV Lawsonia intracellularis 	<ul style="list-style-type: none"> Mycoplasma hyopneumoniae Mycoplasma suis PCV2 PDCoV PEDV Porcine Parvovirus PRRSV-NA PRRSV NADC30-like PRRSV-CN 	<ul style="list-style-type: none"> PRRSV-EU Pseudorabies Virus -gB gene Pseudorabies Virus -gE gene SVA Swine Influenza Virus (Type A) Haemophilus parasuis Streptococcus suis 	Ruminant <ul style="list-style-type: none"> Bluetongue Virus Brucella abortus Brucella melitensis Brucella spp. Crimean-Congo Hemorrhagic Fever FMDV Infectious Bovine Rhinotracheitis PPRV Q Fever 	Emerging Threats <ul style="list-style-type: none"> Bacillus anthracis PL3 Bacillus anthracis pXO1 Bacillus anthracis pXO2 Brucella abortus Brucella melitensis Brucella spp. Chagas Disease 	<ul style="list-style-type: none"> FV3 Rana Virus MERS-CoV Mink Enteritis Virus Paenibacillus larvae R. communis (Ricin) Rift Valley Fever Virus Salmonella spp.

PACKIT Central





In-house Automated PCR Testing Enables Fast Diagnostics in Veterinary Care

Over 190 Assays available!
with more under development coming soon....

www.horiba.com/medical

HORIBA UK Limited | Kyoto Close, Moulton Park, Northampton, NN3 6FL United Kingdom
Tel 01604 542 650 | Fax 01604 542 699 | E-Mail enquiries.uk@horiba.com

POCKIT Central In-house PCR Improving Veterinary Care



The **POCKIT Central** is a bench-top iiPCR Analyser with a small footprint - suitable for fast in-house testing at your Veterinary Practice.

The instrument is very easy to set up and operate, with fully automatic walk-away molecular testing using proven iiPCR

technology for up to 8 samples simultaneously in the same batch, for single or multiple pathogens.

The in-built touch screen allows for user-friendly input and operation, as well as the direct display of test results post analysis.

Once analysis has been initiated, sample handling is fully automated, allowing combined nucleic acid extraction, amplification and detection.

Following each sample analysis run, the **POCKIT Central** automatically initiates a self-cleaning cycle, using ultraviolet to ensure it is ready to process the next batch of samples in quick succession.

The **POCKIT Central** in-house PCR analyser...


Sample in - Result Out

...within 85 minutes!

POCKIT Central



■ Preloaded ready-to-use cartridges



■ Built-in UV lamp for decontamination



■ Barcode scanning available



■ Result printing available