



Experience the power of  
**FluoroMax Plus**  
in your lab



Expand  
your  
**FluoroMax Plus**  
with more  
capability!



...just add a  
2nd detector  
for NIR  
capability



Expand your lab's capabilities with **Near Infrared fluorescence detection** and enjoy extended performance that allows for the analysis of NIR dyes, sensors, photoluminescence samples and so much more.

#### Here's how:

The latest generation of the original, high performance tabletop FluoroMax, can be expanded with NIR fluorescence detection simply by adding a 2nd detector to expand detection from 900-1700nm. Multiple options are available. HORIBA's FluoroMax Plus expansion package\* includes all the NIR accessories, packing, shipping, installation and warranties so that you can acquire high-end research performance with benchtop convenience!

**Contact us for information on expanding your FluoroMax Plus or upgrading an older fluorometer.**

**fluoromax.com**  
email: [info.sci@horiba.com](mailto:info.sci@horiba.com)

\*This upgrade is only for FluoroMax Plus instruments. The older Fluoromax-4 is not compatible, as it is not equipped to add the additional detector port.

