


LabRAM Odyssey® 50th Anniversary Special Edition



Prebook your
LabRAM Odyssey® limited
to 50 units

Pioneers from 1968

50 Years
of **Raman**
Spectroscopy



LabRAM Odyssey®

Commemorating 50 years of Raman Spectroscopy Innovation

Features :

- Ultrafast Raman imaging with Duoscan™ and SWIFT™ XS with built-in EMCCD and high QE proprietary gratings
- Diffraction Limited Spatial Resolution with motorized True Confocal Pinhole continuously adjustable
- Unequalled Spectral Resolution combining high groove density holographic gratings with 800mm focal length spectrograph
- Unique Achromatic Design for an unrivaled Spectral Coverage from 200nm to 2100nm without changing optics
- Fully automated System controlled by Simply Powerful Labspec6™ Spectroscopy Software suite, with built-in EasyNav™, ParticleFinder™, Multivariate Analysis, KIA Horiba Edition™
- Flexible Correlative Microscopy platform ready for upgrade for MicroPL, NanoRaman™, NanoPL, SEM-Raman, Photocurrent Mapping, TCSPC...

For prebooking your **LabRAM Odyssey®** limited to 50 units:

Bring your completed form directly to our HORIBA booth,
or send it to us by email: **info-sci.fr@horiba.com**

- I wish to receive more information about the **LabRAM Odyssey®** 50th anniversary special edition.

Date: _____

First name / Last name:

Organism:

Adresse:

Tel :

Email :

Web :

www. _____

- I agree to be contacted by HORIBA using my contact information above.*
- I agree to receive communications on the services and solutions offered by HORIBA*