



HORIBA Instruments Inc.

Scientific

Kiwan Park

Pondering Piezo Particles

5/10/2016



What we'll talk about today

- Fun Facts
- The Fundamentals
- Applications
- The How's
 - How is Lead Zirconate Titanate (PZT) manufactured?
 - How is particle sizing crucial?
 - 1. Lead (IV) Oxide
 - 2. Lead (II) Oxide
 - 3. Lead Zirconate
 - 4. Lead Zirconate Titanate (PZT): soft and hard
- Q&A

_




Fun Facts

- "Piezo", derived from the Greek piezein, which means to squeeze or press
- It's pronounced as pee-eh-zo, it has nothing to do with a pie
- Piezoelectricity was discovered in 1880 by French physicists Jacques and Pierre Curie.
- Radioactivity was discovered by Marie Curie, Pierre Curie's wife
- 5 Nobel Prizes




The Fundamentals

- Direct piezoelectric current
- Inverse piezoelectric current
- Cyclic





Application – DPI







Image source: https://upload.wikimedia.org/wikipedia/commons/a/ae/Astma-medication.png




Application – Ultrasound




© 2016 HORIBA, Ltd. All rights reserved. Explore the future Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Application - Igniter



https://upload.wikimedia.org/wikipedia/commons/0/02/Lighter_diagram.PNG




Tire pressure monitoring systems (TPMS) embedded inside



3 A

2-4 Weeks

How is PZT manufactured?



Particle sizing required | *Particle size and shape required

HORIBA




Measurement Method

- Laser diffraction
- Mie Scattering theory
- Dispersion method: water with surfactant (e.g. 0.1wt% of Darven C or 0.1wt% of Sodium Pyrophosphate solution)
- Refractive Index @ 650nm: 2.32-3.22
- RI is important for particles <50um (ISO13320)




Lead (II) oxide






Lead (II) oxide





Lead (IV) oxide







© 2016 HORIBA, Ltd. All rights reserved.



HORIBA

Lead Zirconate







Soft PZT / Hard PZT

"Hard" PZT

- "hard" because its dipole or domain is restricted by composition.
- Can with stand high electrical, high temperature and mechanical stresses.
- ✓ Used in applications when high power is required. E.g. welders, ultrasound cleaners.


"Soft" PZT

- "Soft" because dipole or domain isn't restricted by composition
- Has high permittivity, high coupling, high Curie temperature.
- ✓ Used in applications when high sensitivity is required. E.g. medical ultrasound, gas igniters.

Explore the future








© 2016 HORIBA, Ltd. All rights reserved.



RI Assessment





17



Acknowledgment

Special thanks to Wayne Lee of Harris Corp for providing general knowledge of the PZT manufacturing process

www.exelispzt.com

2016 HORIBA, Ltd. All rights reserved.





Email us your questions:

Labinfo@horiba.com

Download this webinar:

http://www.horiba.com/us/en/scientific/products/particlecharacterization/download-center/webinars/

View our full range of particle characterization analyzers:

www.horiba.com/us/particle



Thank you very much for your attention.

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Thank you

감사합니다

Cảm ơn

ありがとうございました

Dziękuję

धन्यवाद

Grazie

Merci

谢谢

நன்ற

ขอบคุณครับ

Obrigado

Σας ευχαριστούμε

Tack ska ni ha

Большое спасибо

Danke

Gracias



Omoshiro-okashiku

Joy and Fun



© 2016 HORIBA, Ltd. All rights reserved.